

+ Wound Management Catalogue

Smith+Nephew

Shaping what's possible
in wound care




Table of Contents

Holistic Wound Management	5
--	---

T.I.M.E. Clinical Decision Support Tool	6
--	---

Products

ACTICOAT [◇] & ACTICOAT [◇] 7	7
---	---

ACTICOAT [◇] Flex.....	8
---------------------------------	---

ALGISITE [◇] M.....	8
------------------------------	---

ALLEVYN [◇] Adhesive	9
-------------------------------------	---

ALLEVYN [◇] Thin.....	10
--------------------------------	----

ALLEVYN [◇] Gentle Border	10
--	----

ALLEVYN [◇] Skin Tear Kit	11
--	----

ALLEVYN [◇] LIFE	12
---------------------------------	----

ALLEVYN [◇] Non-Adhesive	13
---	----

ALLEVYN [◇] Gentle	13
-----------------------------------	----

ALLEVYN [◇] Ag Adhesive & Non-Adhesive.....	14
--	----

ALLEVYN [◇] Ag Gentle Border.....	15
--	----

BACTIGRAS [◇]	16
------------------------------	----

BIOBRANE [◇]	17
-----------------------------	----

CICA-CARE [◇]	17
------------------------------	----

CUTICERIN [◇]	18
------------------------------	----

CUTIFILM [◇] Plus.....	18
---------------------------------	----

CUTILIN [◇]	19
----------------------------	----

CUTIPLAST [◇] / CUTIPLAST [◇] STERIL	19
--	----

DURAFIBER [◇]	20
------------------------------	----

DURAFIBER [◇] Ag.....	21
--------------------------------	----

EXU-DRY [◇] Dressing	22
-------------------------------------	----

FLAMAZINE [◇] Cream	23
------------------------------------	----

INTRASITE [◇] Conformable	24
--	----

INTRASITE [◇] Gel.....	24
---------------------------------	----

IODOSORB [◇] Ointment.....	25
-------------------------------------	----

IODOSORB [◇] Powder.....	25
-----------------------------------	----

IODOSORB [◊] Dressing.....	25
IV3000 [◊]	26
JELONET [◊]	27
LEUKOSTRIP [◊]	27
MELOLIN [◊]	28
MELOLITE [◊]	28
OPSITE [◊] FLEXIFIX [◊]	29
OPSITE [◊] FLEXIFIX [◊] Gentle.....	29
OPSITE [◊] FLEXIGRID	30
OPSITE [◊] Incise.....	30
OPSITE [◊] Post-Op	31
OPSITE [◊] Post-Op Visible.....	32
OPSITE [◊] Spray.....	32
PICO [◊] 7.....	33
PICO [◊] 14.....	34
PRIMABAND [◊]	35
PRIMACREPE [◊]	35
PRIMAFIX [◊] Plus.....	36
PRIMAGAUZE [◊]	37
PRIMAPORE [◊]	37
PRIMASWAB [◊]	38
PRIMATAPE [◊]	38
PROFORE [◊]	40
PROFORE [◊] Lite.....	40
PROSHIELD [◊]	41
REMOVE [◊]	41
RENASYS [◊] TOUCH	42
RENASYS [◊] Foam	43
RENASYS [◊] Gauze	44

RENASYS [◇] Kit Accessories	45
SECURA [◇]	47
S+N Skin Care Lotion.....	48
SKIN-PREP [◇]	48
SOLOSITE [◇] Gel.....	49
UNI-SOLVE [◇]	49
VERSAJET [◇] II.....	50
ZOFF [◇]	51

Training and education from Smith+Nephew	52
---	-----------

Wound Management Resources

Skin Tear Management with ISTAP* Classification Tool	55
IAD Differentiation, Classification and Treatment.....	56
Biofilm management protocol	57
Infection Management Pathway	58
PICO sNPWT hard-to-heal wounds pathway.....	59

Holistic Wound Management



The **T.I.M.E. CDST (clinical decision support tool)** includes holistic patient assessment and supports a multidisciplinary team (MDT) approach.

- The T.I.M.E. principles of wound assessment and management were established by a group of experts and published in 2003¹
- T.I.M.E. has been widely adopted into clinical practice; however, its main focus was on wound management rather than holistic patient assessment, accurate diagnosis and ongoing evaluation of treatment outcomes
- A recent survey of delegates attending the European Wound Management Association (EWMA) annual meeting in 2018 highlighted the need for an easy to use, accurate and practical assessment framework for all healthcare professionals based on T.I.M.E. principles²
- A draft tool was developed based on the results of the EWMA survey, which was assessed and refined by an international group of expert clinicians to form the T.I.M.E. CDST³

Key Points

- TIME is T- Tissue non-viable, I - Infection/inflammation, M- Moisture imbalance, E- Edge of wound
- The experts considered the essential elements of the T.I.M.E. CDST to be holistic patient assessment and involvement of the MDT, followed by wound assessment and control of systemic causes, treating the elements of T.I.M.E. that are impeding wound healing and wound management
- T.I.M.E. CDST uses an evidence-based, 'A, B, C, D, E' approach that incorporates T.I.M.E. principles into a framework that nurses and other professionals can follow within their procedural guidelines
 - Assess patient, wellbeing and wound
 - Bring in MDT and informal carers to promote holistic care
 - Control or treat underlying causes and barriers to wound healing
 - Decide appropriate treatment
 - Evaluate and reassess the treatment and wound management outcomes
- The flow diagram on the following page explains T.I.M.E. in more detail and provides a brief practical guide to holistic wound management. Following each product, you will see a circle symbol. The highlighted portions of the circle correlate to where that product fits within the T.I.M.E. framework and can guide your clinical decision making.

For more information contact your Smith+Nephew representative or National Customer Service on 13 13 60 (Australia), 0800 807 663 (NZ).

1. Schultz GS, Sibbald RG, Falanga V, et al. Wound bed preparation: a systematic approach to wound management. *Wound Repair Regen.* 2003;11 (Suppl 1):S1-28. 2. Ousey K, Gilchrist B, Jaimes H. Understanding clinical practice challenges: a survey performed with wound care clinicians to explore wound assessment frameworks. *Wounds International.* 2018;9(4):58-62. 3. Moore Z, Dowsett C, Smith G, et al. TIME CDST: an updated tool to address the current challenges in wound care. *J Wound Care.* 2019;28(3):154-161.

T.I.M.E. clinical decision support tool

Assess patient, wellbeing and wound

Establish diagnosis and baseline characteristics for appropriate support and comorbidities that may impact healing. Record wound type, location, size, wound bed condition, signs of infection / inflammation, pain location and intensity, comorbidities, adherence / concordance to treatment


Bring in multi-disciplinary team and informal carers to promote holistic patient care

Record referral to others such as surgical team, wound specialist nurse, dietician, pain team, vascular and diabetes team, podiatrist, physiotherapist, family carers and trained counsellor

Control or treat underlying causes and barriers to wound healing

Record management plan for: systemic infection, diabetes, nutritional problems, oedema, continence, mobility, vascular issues, pain, stress, anxiety, non-adherence / concordance with offloading and compression, lifestyle choices

Decide appropriate treatment



*Use appropriate secondary dressing as per your local protocol

Evaluate and reassess the treatment and wound management outcomes

Evaluate: Record wound progression within given timelines. **Flag** if no change, go back to A, B, C and change treatment where indicated

Developed with the support of Glenn Smith³ and Moore *et al.* 2019⁴

*Always ensure adequate blood supply before debriding necrotic tissue. †NPWT: Negative Pressure Wound Therapy. ‡Level of exudate for wounds suitable for NPWT. SECURA Range includes SECURA moisturising creams and lotions, SECURA Barrier Cream D, SECURA No Rinse Cleanser, SECURA No Sting Barrier Film, SECURA No Sting Skin Prep; ALLEVYN GENTLE BORDER Range includes ALLEVYN GENTLE BORDER and ALLEVYN GENTLE BORDER LITE. INTRASITE Range includes INTRASITE Gel and INTRASITE CONFORMABLE. ACTICOAT Range includes ACTICOAT and ACTICOAT FLEX. Reference: 1. Schultz GS, Sibbald RG, Falanga V, *et al.* Wound bed preparation: a systematic approach to wound management. *Wound Rep Reg* (2003);11:1-28. 2. Leaper DJ, Schultz G, Carville K, Fletcher J, Swanson T, Drake R. Extending the TIME concept: what have we learned in the past 10 years? *Int Wound J* 2012; 9 (Suppl. 2):1-19. 3. Smith G, Greenwood M, Searle R. Ward nurse's use of wound dressings before and after a bespoke educational programme. *Journal of Wound Care* 2010;19(9):396-402. 4. Moore Z, Dowsett C, Smith G, *et al.* TIME CDST: an updated tool to address the current challenges in wound care. *Journal of Wound Care* 2019;28(3):154-161.

The products used in the T.I.M.E. clinical decision support tool may vary in different markets. Not all products referred to may be approved for use or available in all markets. Please consult your local Smith & Nephew representative for further details on products available in your market. Intended for healthcare professionals outside of the US only.

Smith & Nephew does not provide medical advice. The information presented is not, and is not intended to serve as, medical advice. For detailed device information, including indications for use, contraindications, precautions and warnings, please consult the product's Instructions for Use (IFU) prior to use. It is the responsibility of healthcare professionals to determine and utilise the appropriate products and techniques according to their own clinical judgment for each of their patients.

Smith+Nephew

ACTICOAT[◊]

Antimicrobial barrier dressing



ACTICOAT dressings contain nanocrystalline silver which kills a broad range of bacteria in as little as 30 minutes.¹⁻⁴ ACTICOAT protects the wound and assists in rapid healing by:

- Killing over 150 wound pathogens, including antibiotic-resistant bacteria⁵
- Shielding against bacterial penetration
- Reducing the risk of bacterial colonisation which may lead to infection
- Helping maintain a moist wound environment in the presence of exudate

ACTICOAT has benefits for both patient and carer:

- Long-acting efficacy means fewer dressing changes

Indications

Full & partial thickness wounds such as:

- Leg ulcers
- Burns
- Skin grafts (Excludes ACTICOAT FLEX)
- Pressure injuries
- Diabetic ulcers
- Infected wounds


ACTICOAT & ACTICOAT 7

ACTICOAT consists of two outer layers of nanocrystalline silver-coated polyethylene net and an inner absorbent layer. ACTICOAT 7 consists of three nanocrystalline silver-coated layers and two inner layers, for seven day efficacy.

Moisten dressing with sterile water or INTRASITE[◊] Gel (do not use saline) and cut to shape.



Code	Product description	Items per unit
ACTICOAT		
66000808	5cm x 5cm	Box/5
66000789	10cm x 10cm	Box/5
66000791	10cm x 10cm	Box/12
66000792	10cm x 20cm	Box/12
66000793	20cm x 40cm	Box/6
66000794	40cm x 40cm	Box/6
66000795	10cm x 120cm	Box/6
ACTICOAT 7		
66000809	5cm x 5cm	Box/5
66000796	10cm x 12.5cm	Box/5
66000797	15cm x 15cm	Box/5



1.Wright JB, et al. Wound management in an era of increasing bacterial antibiotic resistance: a role for topical silver treatment. Am J Infect Cont 1998; 26: 572-577. 2.Wright JB, et al. The comparative efficacy of two antimicrobial barrier dressings: in vitro examination of two controlled release silver dressings. Wounds 1998; 10: 179-188. 3. Wright JB, et al. Efficacy of topical silver against fungal burn wound pathogens. Am J Infect Cont 1999; 27: 344-350. 4.Yin HQ, et al. Comparative evaluation of the antimicrobial activity of ACTICOAT antimicrobial barrier dressing. J Burn Care Rehab 1999; 20:195-200. 5. Data on file. Broad-spectrum antibacterial activity of nanocrystalline silver. 2003; report 0109003.


ACTICOAT[◇] Flex

Antimicrobial barrier dressing

ACTICOAT Flex is a knitted polyester weave, coated in nanocrystalline silver. The dressing is highly flexible and conformable and stretches to allow patient movement. The open weave design allows easy fluid and wound exudate migration. ACTICOAT Flex can be used as a wound contact layer with Negative Pressure Wound Therapy including PICO and RENASYS.

If wound has only low levels of exudate, moisten the dressing with water or INTRASITE Gel (DO NOT use saline) before application.

Code	Product description	Items per unit
ACTICOAT Flex 3		
66800396	5cm x 5cm	Box/5
66800398	10cm x 10cm	Box/5
66800399	10cm x 10cm	Box/12
66800409	10cm x 20cm	Box/12
66800419	20cm x 40cm	Box/6
66800432	40cm x 40cm	Box/6
66800435	10cm x 120cm	Box/6
66801290	4cm x 15cm	Box/ 5
66801291	4cm x 25cm	Box/ 5
66801292	4cm x 35cm	Box/ 5
66801293	20cm x 20cm	Box/ 5
ACTICOAT Flex 7		
66800395	5cm x 5cm	Box/5
66800397	10cm x 12.5cm	Box/5
66800420	15cm x 15cm	Box/5
66800400	20cm x 40cm	Box/6
66800401	40cm x 40cm	Box/6
66800545	2.5cm x 60cm (Ribbon)	Box/5



ALGISITE[◇] M

Calcium alginate dressing

ALGISITE M forms a soft gel when it comes in contact with wound exudate, providing a moist wound environment to promote rapid healing. ALGISITE M absorbs excess exudate, is easy to remove and retains its shape when wet, protecting delicate healing tissues even during dressing changes.

In vitro studies¹ have shown that ALGISITE M can absorb up to ten times its own weight in exudate, which means it can be left in place for up to 4 days.^{2,3} The fibres in ALGISITE M have been specially bonded to reduce fibre shed. ALGISITE M is soft and highly conformable to wound contours and is easily removed in one piece.

Indications

Wounds with moderate to heavy exudate such as:

- Pressure injuries
- Venous leg ulcers
- To facilitate the control of minor bleeding

Code	Product description	Items per unit
66000519	5cm x 5cm	Box/10
66000520	10cm x 10cm	Box/10
66000521	15cm x 20cm	Box/10
66000522	2cm x 30cm (rope/strip)	Box/5



1. McCann A. An in vitro comparison of ALGISITE M versus other alginate dressings, and a hydro-fibre dressing, Aquacel. Data on file; 2000: report PS/WR/00/05003a. 2. Smith+Nephew 2021.PMCF Survey Data for ALGISITE M Calcium Alginate Dressing. Internal Report. EO.AWM.PCS77.002.v1. 3. Smith+Nephew 2020.PMS 303: ALGISITE M Calcium Alginate Dressing: PMCF Activity Summary Report. Internal Report. EO.AWM.PCS77.001.v1.

ALLEVYN[◊] Range



Foam dressings

ALLEVYN is a comfortable, highly absorbent dressing that keeps wounds moist and clean. The dressing can remain in place for up to seven days. Rapid healing combined with infrequent dressing changes makes ALLEVYN the right choice for patients and budgets.

ALLEVYN 3 layer construction provides a moist environment ideal for healing:

- Breathable moisture responsive film allows extra moisture to evaporate away from the skin, while providing a waterproof barrier that prevents bacterial contamination.
- Hydrocellular pad absorbs and retains exudate for up to seven days, reducing the risk of maceration and maintaining a moist wound environment ideal for healing.
- Perforated polyurethane wound contact layer allows even viscous exudate to pass into the dressing, while ensuring the dressing holds its shape and minimising disturbance of the wound on removal.

ALLEVYN[◊] Adhesive

Foam dressing

Adhesive, but non-adherent, dressing for wounds with moderate to high levels of exudate. Provides soft, cushioned protection for up to seven days.

Indications

Wound management by secondary intention on shallow, granulating wounds, chronic and acute exudative wounds, full and partial thickness wounds such as:

- Pressure injuries
- Leg ulcers
- Diabetic foot ulcers
- Surgical wounds
- Donor sites
- Fungating wounds
- Can be used in conjunction with INTRASITE[◊] Gel for necrotic or sloughy wounds



Code	Product description	Items per unit
66000043	7.5cm x 7.5cm; pad: 5cm x 5cm	Box/10
66000599	10cm x 10cm; pad: 7.5cm x 7.5cm	Box/10
66000044	12.5cm x 12.5cm; pad: 10cm x 10cm	Box/10
66000045	17.5cm x 17.5cm; pad: 15cm x 15cm	Box/10
66000046	22.5cm x 22.5cm; pad: 20cm x 20cm	Box/10
66000451	22.5cm x 22.1cm; sacrum – large	Box/10
66000700	17cm x 17cm; sacrum – small	Box/10

ALLEVYN[◊] Thin

Foam dressing

ALLEVYN Thin dressings offer all the proven benefits of moist wound healing, without any breakdown of the dressing caused by contact with exudate. This promotes early, clean and trouble-free healing of wounds with light to moderate levels of exudate.

Indications

Indicated for the management of surface wounds with low to moderate exudate such as minor burns and abrasions. May also be used under the direction of a healthcare professional for wounds such as:

- Post-op wounds
- Donor sites

Code	Product description	Items per unit
66047576	5cm x 6cm	Box/10
66047578	10cm x 10cm	Box/5
66047579	15cm x 20cm	Box/3

ALLEVYN[◊] Gentle Border

Foam dressing

ALLEVYN Gentle Border is a hydrocellular foam dressing with a silicone gel adhesive contact layer. It is indicated for wounds on patients with fragile or sensitive skin.

With a wide range of shapes and sizes, ALLEVYN Gentle Border dressings are designed to conform to the body,¹⁻³ helping to enable patient comfort¹ during wear. The dressing can be cut too.⁴

The silicone gel adhesive allows reapplication as necessary (if using in pressure injury prevention) and reduces trauma on removal. The ALLEVYN foam prevents leakage and may help reduce the risk of maceration of the wound and surrounding skin meaning fewer painful dressing changes and less stress for patients.


The breathable top film acts as a protective barrier, reducing the risk of infection and contamination. The film is moisture responsive, reacting to changing skin conditions – breathability increases as exudate levels increase to prevent unnecessary moisture retention.

Indications

Wound management by secondary intention on shallow, granulating wounds, chronic and acute exudative wounds, full and partial thickness wounds such as:

- Pressure injuries. Can be used for pressure injury prevention on intact skin, including pressure injuries caused by medical devices, as part of a pressure injury prevention protocol
- Leg ulcers
- Diabetic foot ulcers
- Infected wounds*^
- Malignant wounds^
- Surgical wounds
- Superficial epidermal to mid dermal burns
- Donor sites^
- Skin tears
- Fungating wounds^
- Can be used in conjunction with INTRASITE[◊] Gel for necrotic or sloughy wounds

Please turn the page for order codes →





Code	Product description	Items per unit
ALLEVYN Gentle Border		
66800269	7.5cm x 7.5cm	Box/10
66800270	10cm x 10cm	Box/10
66800900	10cm x 20cm	Box/10
66800264	10cm x 25cm	Box/10
66800265	10cm x 30cm	Box/10
66800272	12.5cm x 12.5cm	Box/10
66800273	17.5cm x 17.5cm	Box/10
ALLEVYN Gentle Border Heel		
66800506	23cm x 23.2cm	Box/5
ALLEVYN Gentle Border Multisite		
66800959	17.1cm x 17.9cm	Box/10
ALLEVYN Gentle Border Sacrum		
66800898	16.8cm x 17.1cm	Box/10
ALLEVYN Gentle Border Lite		
66800833	5cm x 5cm	Box/10
66800836	5.5cm x 12cm	Box/10
66801013	8.6cm x 7.7cm Oval	Box/10
66801041	8cm x 8.4cm Multisite	Box/10

*Where the product is used on infected wounds the infection should be inspected and treated as per local clinical protocol. ^Not indicated for ALLEVYN Gentle Border Lite.

1. Hurd T, et al. A multi-centre in-market evaluation of Allevyn GB. Wounds UK 2009; 5(3): 32–44. 2. Palmer S and Mistry C. Smith & Nephew Wound Management, 2011 Data on File . Report Reference OR-DOF/017. An open, prospective, randomised, comparative study to compare the performance of ALLEVYN, Gentle Border Multisite with an alternative silicone adhesive dressing. 3. Callaghan R et al. Case series evaluation: the use of Allevyn Gentle Border Multisite on chronic wounds. Wounds UK 2012; 8(4):110-118. 4. Smith & Nephew Data on File DS/14/318/R - Cutting of Allevyn Variants.

ALLEVYN^o Skin Tear Kit

Skin tear management

ALLEVYN Skin Tear Kit is a simple but effective tool to help with treatment standardisation and to help implementation of evidence-based principles. It encourages clinicians to assess and manage skin tears effectively.

Convenient and portable, ALLEVYN Skin Tear Kit contains a skin tear assessment tool – the iSTAP guideline together with sterile tools including an ALLEVYN silicone dressing that can be left in place for up to 7 days¹, while maintaining the optimal balance in fluid for moist wound healing.²



Code	Product description	Items per unit
36361382	ALLEVYN Gentle Border Lite Oval Skin Tear Kit	Each

1. Volunteer Trial Report Reference VTR/HVT 050. Results from an open, prospective, randomised, within volunteer comparison of Dakota and Mepilex Border. Palmer S and Smith G. December 2008. 2. Smith & Nephew Wound Management Data on File Report Reference DS/10/100/R1. ALLEVYN Gentle Border Lite Physical Properties. Tompkins L. August 2010.



Watch the video:
How to use the Skin Tear Kit with ALLEVYN GENTLE BORDER Lite



Foam dressing

ALLEVYN LIFE is a multi-layered dressing incorporating a hydrocellular foam, a hyper-absorber lock away core and a masking layer. It is indicated for wounds on patients who have shallow, granulating wounds, chronic and acute exudative wounds, full and partial thickness wounds including pressure injuries, skin tears, diabetes-related foot ulcers and venous leg ulcers. ALLEVYN LIFE can be left in place for up to seven days and its design minimises wastage and dressing changes:

- Quadrilobe shape* and wider border designed to fit the contours of the human body so it conforms securely, and allows patients to shower.¹
- Lock away layer provides peace of mind. It has excellent fluid locking under pressure², leakage prevention³ and wound odour absorption when tested in-vitro on representative compounds.⁴
- Masking layer which minimises the visual impact of exudate. The masking layer avoids unnecessary dressing changes and wastage (see diagram on the right).
- ALLEVYN LIFE's multilayered design (hydrocellular foam, hyper-absorbent lock away core and masking layer) provides cushioning and redistributes pressure.⁵

Indication

Wound management by secondary intention on shallow, granulating wounds, chronic and acute exudative wounds, full and partial thickness wounds including:


- Pressure injuries
- Leg ulcers
- Diabetic foot ulcers
- Surgical wounds
- Superficial epidermal to mid dermal burns
- Donor sites
- Skin tears
- Fungating wounds
- Can be used in conjunction with INTRASITE^o GEL for necrotic or sloughy wounds

ALLEVYN LIFE can be used for pressure injury prevention on intact skin as part of a pressure injury prevention protocol.⁶

Code	Product description	Items per unit
66801067	Small – 10.3cm x 10.3cm	Box/10
66801068	Medium – 12.9cm x 12.9cm	Box/10
66801069	Large – 15.4cm x 15.4cm	Box/10
66801070	Extra Large – 21cm x 21cm	Box/10
66801304	Heel 25cm x 25.2cm	Box/5
66801306	Sacrum 17.2cm x 17.5cm	Box/10
66801307	Sacrum 21.6cm x 23cm	Box/10

1. Palmer S, Dharma H. An open, prospective, randomised, comparative volunteer trial to compare the performance of silicone adhesive dressings. Data on file; 2012: report OR-DOF/020. 2. Roberts S. Fluid handling properties of ALLEVYN Life. Data on file; 2012; report DS/12/124/DOF. 3. Roberts S. Physical properties of ALLEVYN Life. Data on file; 2012; report DS/12/123/DOF. 4. Roberts S. Odour reducing properties of ALLEVYN Life. Data on file; 2012: report DS/12/127/DOF. 5. Stephen-Haynes J, et al. The clinical performance of a silicone foam in an NHS community trust. J Comm Nurs 2013; 27(5): 50-59. 6. Forni C, D'Alessandro F, Gallerani P, et al. Effectiveness of using a new polyurethane foam multi-layer dressing in the sacral area to prevent the onset of pressure ulcer in the elderly with hip fracture: A pragmatic randomized controlled trial. Int Wound J. 2018;1-8.

*ALLEVYN LIFE is available in different shaped dressing shapes such as sacrum and heel.



No border coverage dressing can remain in place with exudate masked



50% border coverage consider changing dressing



75% border coverage change dressing



ALLEVYN^o Non-Adhesive



Foam dressing


For general use on wounds with moderate to high levels of exudate. Ideal for use where surrounding skin is friable. Provides soft, cushioned protection for up to seven days.

Indications

Wound management by secondary intention on shallow, granulating wounds. Chronic and acute exudative wounds, full and partial thickness wounds such as:

- Pressure injuries
- Leg ulcers
- Diabetic foot ulcers
- Superficial epidermal to mid dermal burns
- Donor sites
- Fungating wounds
- Can be used in conjunction with INTRASITE^o Gel for necrotic or sloughy wounds

ALLEVYN Heel contributes to a pressure relieving protocol when used in conjunction with pressure relieving devices.



Code	Product description	Items per unit
66007643	5cm x 5cm	Box/10
66007637	10cm x 10cm	Box/10
66007335	10cm x 20cm	Box/10
66000093	15cm x 15cm	Box/10
66007638	20cm x 20cm	Box/10
66007630	10.5cm x 13.5cm (heel)	Box/5
66007640	9cm x 9cm (tracheostomy)	Box/10*

*Not available in New Zealand

ALLEVYN^o Gentle

Foam dressing

ALLEVYN Gentle is a silicone gel adhesive dressing that provides the optimal moist wound healing environment to promote the healing of moderately and highly exuding wounds. The absorbent pad comprises of hydrocellular foam. The silicone adhesive layer is gentle, even to fragile skin and makes it easy to apply, reposition and remove. A breathable outer film prevents bacterial contamination and is showerproof when used with appropriate secondary retention.

Indications

Chronic and acute full thickness, partial thickness or shallow granulating, exuding wounds such as:

- Leg ulcers
- Diabetic foot ulcers
- First and second degree burns
- Skin graft donor sites
- Skin tears



Code	Product description	Items per unit
66802128	5cm x 5cm	Box/10
66802129	10cm x 10cm	Box/10
66802131	15cm x 15cm	Box/10
66802133	20cm x 50cm	Box/10

ALLEVYN[◇] Ag

Antimicrobial foam dressings



The ALLEVYN Ag range combines effective exudate management and antimicrobial action in one dressing. ALLEVYN Ag is a unique mix of the existing ALLEVYN triple layer foam technology with silver sulfadiazine (SSD) providing long wear time, up to seven days antimicrobial properties.¹

The SSD is released into the wound bed for up to seven days:

- Effective against a broad range of pathogens, including bacteria, yeast, fungi and antibiotic-resistant strains.¹
- Easy to use – silver is released when in contact with wound fluid.

Indications

Chronic and acute full thickness, partial thickness or shallow granulating, exuding wounds such as:

- Pressure injuries
- Venous ulcers
- Diabetic ulcers
- Burns
- Fungating/malignant wounds
- Surgically dehisced wounds
- May be used on infected wounds

ALLEVYN[◇] Ag Adhesive & Non-Adhesive


ALLEVYN Ag is a highly absorbent antimicrobial foam dressing designed for use on exuding wounds.

It has the same fluid handling technology of ALLEVYN, which kills bacteria by releasing levels of silver (44-77ppm) for up to seven days.²⁻⁴

This easy-to-use dressing is available in a non-adhesive and adhesive format.

1. Woodmansey E. Antimicrobial activity of Allevyn Ag Non-Adhesive against a broad spectrum of microorganisms. Data on file; 2007: report 0703006. 2. Daubney L. ALLEVYN Ag dressings silver release testing. Data on file; 2007: report 0706038. 3. Driffield K. Antimicrobial activity of ALLEVYN Ag dressings in comparison to other silver-based antimicrobial dressings against MRSA using a dynamic shake flask method. Data on file; 2007: report 0708055. 4. Woodmansey E. Antimicrobial activity of ALLEVYN Ag Non-Adhesive and ALLEVYN Ag Gentle against a broad spectrum of wound pathogens. Data on file; 2009: report 0903004.

Please turn the page for order codes →



ALLEVYN[®] Ag Dressings

Antimicrobial foam dressing

ALLEVYN Ag Dressings are antimicrobial hydrocellular foam dressings with and without a silicone gel adhesive. They are indicated for wounds on patients with fragile or sensitive skin where the risk of infection is suspected or needs to be managed.

ALLEVYN Ag Gentle Border protects sensitive patients in three ways:

- The silicone gel adhesive reduces trauma on removal for patients with fragile skin
- The silver sulfadiazine (SSD) provides sustained antimicrobial protection
- The ALLEVYN foam helps to prevent leakage and minimise the risk of maceration of the wound and surrounding skin – fewer painful dressing changes means less stress for your patients



Code	Product description	Items per unit
ALLEVYN Ag Non-Adhesive		
66800083	5cm x 5cm	Box/10
66800086	10cm x 10cm	Box/10
66800089	15cm x 15cm	Box/10
66800091	20cm x 20cm	Box/10
ALLEVYN Ag Adhesive		
66800073	7.5cm x 7.5cm (Pad 5cm x 5cm)	Box/10
66800075	10cm x 10cm (Pad 7.5cm x 7.5cm)	Box/10
66800078	12.5cm x 12.5cm (Pad 10cm x 10cm)	Box/10
66800081	17.5cm x 17.5cm (Pad 15cm x 15cm)	Box/10
ALLEVYN Ag Gentle Border		
66800460	7.5cm x 7.5cm	Box/10
66800461	10cm x 10cm	Box/10
66800462	12.5cm x 12.5cm	Box/10
66800463	17.5cm x 17.5cm	Box/10

BACTIGRAS[◇]

Medicated Tulle Gras

BACTIGRAS is a gauze of leno weave impregnated with White Soft Paraffin BP containing 0.5% Chlorhexidine Acetate BP in a sterile presentation.

BACTIGRAS is effective against a wide range of gram-positive and gram-negative wound pathogens, including MRSA.^{1,2} BACTIGRAS is for use as a topical treatment only. It is indicated for a wide range of wounds where there is a risk of infection. For already infected wounds, BACTIGRAS may be used in conjunction with systemic antibiotics.

The BACTIGRAS tulle gras dressing, made from an open-weave gauze, has low adherence and allows the wound to drain freely into an absorbent secondary dressing.

Indications


It is indicated for a wide range of wounds such as:

- Minor burns and scalds
- Lacerations, abrasions and other skin loss wounds
- Donor and recipient graft sites

Code	Product description	Items per unit
7456	5cm x 5cm	Box/50
66003650	10cm x 10cm	Box/50*
7457	10cm x 10cm	Box/10
7461	15cm x 20cm	Box/10*
66003661	10cm x 40cm	Box/10*
66007505	15cm x 1m	Roll/1

*Not available in New Zealand

1. Holland KT, Davis W. A note on an in vitro test system to compare the bacterial properties of wound dressings. *J App Bact* 1985; 59: 61-63. 2. Kobayashi H, et al. Bactericidal effects antiseptics and disinfectants against Methicillin- resistant *Staphylococcus aureus*. *Infect Cont Hosp Epidem* 1989; 10: 562-564.



BIOBRANE[◇]

Temporary Biosynthetic Skin Dressing

BIOBRANE is a biosynthetic dressing made from an ultrathin, semipermeable, perforated silicone membrane that is mechanically bonded to a flexible knitted trifilament nylon fabric.

Denatured porcine dermal collagen has been bonded to the nylon/silicone membrane to provide a flexible and conformable dressing with adherence properties and a hydrophilic surface.

Indications

BIOBRANE is indicated for:

- Covering clean partial thickness burn wounds
- Split thickness donor sites

It is widely used in the management of burns and other skin loss injuries.^{1,2}

Code	Product description	Items per unit
66800174	13cm x 13cm	Box/5
66800175	13cm x 38cm	Box/5
66800176	25cm x 38cm	Box/5
66800177	38cm x 50cm	Box/1
66800182	Small Glove	Box/2
66800183	Medium Glove	Box/2
66800184	Large Glove	Box/2

1. Greenwood, JE. Development of patient pathways for the surgical management of brain injury. ANZ J Surg 2006; 76: 805-811. 2. Lang EM, et al. Biobrane in the treatment of burn and scald injuries in children. Ann Plast Surg 2005; 55: 485-489.

CICA-CARE[◇]

Silicone gel sheet for scar management

CICA-CARE is a soft, self-adhesive silicone gel sheet for the treatment and prevention of scarring. CICA-CARE is effective in the treatment of scars up to 30 years old.¹

The CICA-CARE adhesive gel sheet is simple to apply and can be cut to size to fit most scars. The silicone membrane and gel combination make CICA-CARE more durable. It is washable and reusable.


Indications

CICA-CARE is designed for temporary use:

- In the management of both existing and new hypertrophic scars and keloids and
- As a prophylactic therapy on healed wounds to help to prevent hypertrophic scarring and keloids

Code	Product description	Items per unit
66250704	12cm x 6cm	Box/1
66250706	12cm x 15cm	Box/1
66250707	12cm x 15cm	Box/10

1. Smith & Nephew 1999. An Investigation of the Use of CICA-CARE in the Treatment of Hypertrophic and Keloid Scars: A Clinical, Histological and Immunohistochemical Study. Hull. CRR/00000/CT96/06.



CUTICERIN[◇]


Low Adherent Tulle Gras Dressing

Used in place of a tulle gras dressing, CUTICERIN is an impregnated, non-aqueous, non-linting dressing ideal as a wound contact layer. It reduces the risk of granulation tissue growing through the dressing, decreasing trauma on removal.

Indications

Management of superficial exuding wounds such as:

- Burns
- Radiation injuries
- Abrasions
- Split-thickness skin graft donor sites



Code	Product description	Items per unit
66045560	7.5cm x 7.5cm	Box/50
66045561	7.5cm x 20cm	Box/50
66045502	20cm x 40cm	Box/25*

*Not available in New Zealand

CUTIFILM[◇] Plus


Post-operative dressing

CUTIFILM Plus is a waterproof dressing ideal for use on surgical wounds, lacerations and abrasions.

Indications

A dressing for superficial wounds such as:

- Blisters
- Minor burns
- Abrasions and lacerations
- Post operative wounds (e.g. sutures)
- To protect the skin from friction and external contamination



Code	Product description (dressing; pad)	Items per unit
76248	5cm x 7.2cm; pad: 2.7cm x 4cm	Box/50
76244	8cm x 10cm; pad: 4cm x 6.5cm	Box/25
76245	8cm x 15cm; pad: 4cm x 11cm	Box/25
76246	10cm x 12.5cm; pad: 5.5cm x 8.5cm	Box/25
76247	10cm x 20cm; pad: 5.5cm x 16cm	Box/25
76249	10cm x 30cm; pad: 5.5cm x 24cm	Box/50
76250	10cm x 38cm; pad: 5.5cm x 32cm	Box/50*

Tan pad

76251	5cm x 7.2cm; pad: 2cm x 4cm	Box/50*
76252	8cm x 10cm; pad: 4cm x 6.5cm	Box/25*

*Not available in New Zealand

CUTILIN[◇]

Low-adherent absorbent dressing


CUTILIN is a low adherent, sterile, absorbent dressing.

Indications

- Low to moderately exuding wounds
- Can be cut to shape without linting
- It is available in sterile or non-sterile presentations

Code	Product description	Items per unit
76240	5cm x 5cm	Box/100
76242	10cm x 10cm	Box/100
76243	10cm x 20cm	Box/100
76238	45cm x 7m Roll (Non-sterile)	Roll/1*

*Not available in New Zealand



CUTIPLAST[◇]/ CUTIPLAST STERIL


Fabric Island Dressing

CUTIPLAST is a flexible fabric dressing with a low-adherent wound pad ideal for use in areas of curvature and flexion such as knees and elbows and hirsute areas.

Indications

It is indicated for dressing surgical wounds or general cuts and lacerations, where a waterproof dressing is not required.

Code	Product description (dressing; pad)	Items per unit
66001478	7.2cm x 5cm; pad: 4cm x 2.5cm	Box/100
66001473	10cm x 8cm; pad: 5.3cm x 4.5cm	Box/50
66001474	15cm x 8cm; pad: 11cm x 4cm	Box/50
66001475	20cm x 10cm; pad: 15.5cm x 5.3cm	Box/50
66001476	25cm x 10cm; pad: 20cm x 5.3cm	Box/50
66001477	30cm x 10cm; pad: 25cm x 5.3cm	Box/50
66001464	4cm x 5m Roll; pad: 1.8cm x 5m (non-sterile)	Roll/1
66001466	6cm x 5m Roll; pad: 2.5cm x 5m (non-sterile)	Roll/1
66001468	8cm x 5m Roll; pad: 4cm x 5m (non-sterile)	Roll/1



DURAFIBER^o

Absorbent gelling fibrous dressing

DURAFIBER is a hydrophilic gelling fibrous filler dressing with rapid fluid absorption and retention of fluid. It has a high level of structural integrity and wet strength for the management of medium to highly exuding wounds.¹

NB: Whilst DURAFIBER assists in the management of wounds prone to bleeding, it is not intended to be used as a surgical sponge in heavily bleeding wounds.

DURAFIBER requires a secondary or retention dressing such as ALLEVYN^o or OPSITE^o.

DURAFIBER is designed to be left in place for up to seven days, depending on the amount of exudate present and the condition of the wound.


Indications

DURAFIBER is an absorbent, gelling dressing for the management of chronic and acute, full thickness, partial thickness, or shallow granulating exuding wounds. For example:

- Leg ulcers
- Pressure injuries
- Diabetic ulcers
- Surgical wounds
- Wounds left to heal by secondary intent
- Donor sites
- Tunnelling and fistulae wounds
- Traumatic wounds

Code	Product description	Items per unit
66800559	5cm x 5cm	Box/10
66800560	10cm x 10cm	Box/10
66800561	15cm x 15cm	Box/5
66800563	2cm x 45cm (rope)	Box/5
66800546	4cm x 10cm	Box/5
66800547	4cm x 20cm	Box/5
66800548	4cm x 30cm	Box/5

1. Dowler A. DURAFIBER dressing physical properties, Data on file, 2010: report DS/10/056/R1.



DURAFIBER^o Ag

Absorbent Gelling Silver Fibrous Dressing

DURAFIBER Ag is an absorbent, non-woven, silver containing antimicrobial dressing composed of cellulose ethyl sulphonate fibres. The ionic silver in the dressing provides antimicrobial activity against a broad spectrum of common wound pathogens¹ which may help to reduce bacterial bioburden and the risk of infection.

Indications

Chronic and acute, full thickness, partial thickness, or shallow granulating exuding wounds such as:

- Leg ulcers
- Pressure injuries
- Diabetic ulcers
- Surgical wounds
- Traumatic wounds
- Donor sites
- Partial thickness burns
- Tunnelling and fistulae wounds
- May be used on infected wounds



Code	Product description	Items per unit
66800578	5cm x 5cm	Box/10
66800579	10cm x 10cm	Box/10
66800580	15cm x 15cm	Box/5
66800581	20cm x 30cm	Box/5
66800582	2cm x 45cm	Box/5
66800583	4cm x 10cm	Box/5
66800584	4cm x 20cm	Box/5
66800585	4cm x 30cm	Box/5

1. Vaughan KL, et al. Antimicrobial activity of DURAFIBER Ag against bacteria, yeast and fungi commonly found in wounds over a 7 day period. Data ion file:2010:report 1009012.


EXU-DRY[◇] Dressing



Specialised absorbent dressings with anti-shear layer

EXU-DRY is a one-piece, highly absorbent, non-adherent dressing. It has been designed to save nursing time and lower total dressing costs. Available in a wide range of sizes and shapes, EXU-DRY may be used as a wet or dry dressing on clean or contaminated wounds. EXU-DRY is compatible with most topical agents, liquids and creams.[#]

EXU-DRY is a multi-layer dressing, replacing multiple layers of gauze and abdominal or surgical pads. The outer layer is non-occlusive and permeable, while the absorbent layer wicks away excess exudate, reducing the risk of maceration.



Indications

May be used in the management of superficial to full thickness wounds.

- Ulcers (venous, arterial, diabetic)
- Pressure injuries
- Donor sites
- Surgical incisions and excisions
- Skin grafts
- Bio-engineered skin substitutes
- Burns
- Draining wounds
- Moist skin desquamation
- Fungating neoplasms
- Chemotherapy extravasation
- Kaposi's sarcoma
- Lymphoedema
- Graft vs. host disease
- Dermatological wounds
- Skin sloughing disorders



[#]Petroleum jelly based products should not be used in conjunction with EXU-DRY as they interfere with absorbency.

Code	Product description (size; absorbency)	Unit of measure	Min. Order Qty
5999003	7.5cm disc; medium	Each	50
5999034	7.5cm x 10cm; full	Car/50*	1
5999004120	10cm x 15cm; full	Car/10	12
5999006	15cm x 23cm; full	Car/12	1
5999009	23cm x 38cm; full	Each	30
5999018	38cm x 46cm; full	Each	30
5999024	38cm x 61cm; full	Each	30
5999028	50cm x 71cm; full	Each*	20
Slit Tube Dressings			
5999003S	8cm disc; medium	Each*	50
5999PTM	5cm x 8cm; full	Each*	50
5999034S	8cm x 10cm; full	Each*	50
Specialty Dressings			
5999LPG	35cm x 35cm; hand dressing (glove shaped); heavy	Each	20
5999LPL	87cm x 95cm; leg dressing (shaped); medium	Car/20*	1
5999M36	61cm x 91cm; non-permeable; heavy	Each*	15
5999L72	91cm x 183cm; non-permeable	1 Dressing**	15
5999L73	91cm x 183cm; non-permeable with straps	Each*	15

*Not available in New Zealand

**Not available in Australia

FLAMAZINE[◇] Cream

1% w/w silver sulfadiazine

FLAMAZINE is a topical antibacterial cream containing silver sulfadiazine. FLAMAZINE is available on prescription only.

PBS Information

50g Tube: restricted benefit. Prevention and treatment of infection in partial or full skin thickness loss due to burns; Prevention and treatment of infection in partial or full skin thickness loss due to epidermolysis bullosa; Stasis ulcers. 500g Jar: This size is not listed on the PBS.

Indications

- Used to prevent infection in severe burns
- Pressure injuries
- Leg ulcers

Code	Product description	Items per unit
66800531	50g Tube	Tube/1* Carton/36 Tubes
66800532	500g Jar	Jar/1

*Pharmacy-Only Medicine in New Zealand

Please review Product Information before prescribing. Product Information can be accessed from <https://www.ebs.tga.gov.au/ebs/picmi/picmirepository.nsf/PICMI?OpenForm&t=pi&q=flamazine> (Australia) and <https://www.medsafe.govt.nz/profs/Datasheet/f/flamazinecrm.pdf> (New Zealand)



INTRASITE[◇] Conformable

Conformable Hydrogel Dressing

INTRASITE Conformable is a hydrogel non-woven dressing impregnated with INTRASITE Gel. It conforms to the wound without adhering to it, for gentle packing of deep, shallow, open or undermined wounds. INTRASITE Conformable helps create a moist wound environment to speed healing.

Indications

INTRASITE Conformable is indicated for use in open wounds:

- Pressure injuries
- Surgical wounds
- Malignant wounds
- Shallow wounds
- Excoriated skin
- Radiation burns

Code	Product description	Items per unit
66000324	10cm x 10cm	Box/10
66000325	10cm x 20cm	Box/10
66000326	10cm x 40cm	Box/10*

*Not available in New Zealand



INTRASITE[◇] Gel

Hydrogel Wound Dressing

INTRASITE Gel is an amorphous hydrogel that gently and effectively debrides and desloughs necrotic and sloughy wounds. INTRASITE Gel consists of a cross-linked polymer that allows the INTRASITE product to respond to the level of moisture within the wound, and balance absorption and rehydration appropriately.

INTRASITE Gel is a sterile, clear, free-flowing gel, presented in a convenient applicator pack for controlled placement of the gel. It works in two ways to promote fast, comfortable healing:

- By loosening and absorbing slough and exudate into the gel, without damaging fragile granulation tissue
- Through rehydration of the wound as water is donated from the gel, creating a moist healing environment at the wound interface

Indications

INTRASITE Gel is indicated for the removal of non-viable tissue from shallow, undermined and deep wounds such as:

- Pressure injuries
- Leg ulcers
- Diabetic foot ulcers
- Malignant wounds
- Burns, scalds & radiation damage
- Open surgical wounds
- Lacerations & grazes
- Fungating ulcers
- Granulating cavity wounds

Code	Product description	Items per unit
7308	8g Applipak	Box/10
7311	15g Applipak	Box/10
7313	25g Applipak	Box/10



IODOSORB[◇]

0.9% Cadexomer Iodine

IODOSORB is a sterile formulation of cadexomer iodine (0.9%) indicated for the treatment and healing of chronic ulcers. The iodine in IODOSORB is bound within the cadexomer smart micro-beads. When IODOSORB comes in contact with the wound, it absorbs exudate and debris and swells, gradually releasing iodine over time.

Changes colour from brown to white when the dressing requires changing, it helps to reduce odour associated with heavily contaminated wounds.¹

IODOSORB works in four ways:

1. Disrupts biofilms.²
2. Kills micro-organisms, reducing the bioburden.³
3. Removes exudate and slough, promoting a clean wound bed and reducing the risk of infection.⁴
4. Forms a gel which creates a moist wound healing environment, protecting delicate granulation tissue.⁵

Indications

- Treatment of chronic exuding wounds
- Can be used under compression therapy
- May be used on infected wounds



IODOSORB[◇] Ointment

Code	Product description	Items per unit
66051240	10g Tube	Pack/4
66051230	20g Tube	Pack/2



IODOSORB[◇] Powder

Code	Product description	Items per unit
66051070	3g Powder	Pack/7*



IODOSORB[◇] Dressing

Code	Product description	Items per unit
66051330	5g (6cm x 4cm)	Pack/5
66051340	10g (8cm x 6cm)	Pack/3
66051360	17g (10cm x 8cm)	Pack/2



*Not available in New Zealand

1. Lindsay G, et al. A study in general practice of the efficacy of cadexomer iodine in venous leg ulcers treated on alternate days. *Acta Therapeutica* 1986; 544 (sup): 60-61. 2. Akiyama H, et al. Assessment of cadexomer iodine against *Staphylococcus aureus* biofilm in vivo and in vitro using confocal laser scanning. 3. Mertz PM, et al. Can antimicrobials be effective without impairing wound healing? The evaluation of a cadexomer iodine ointment. *Wounds* 1994; 6: 184-193. 4. Hansson C. The effects of cadexomer iodine paste in the treatment of venous leg ulcers compared with hydrocolloid dressing and paraffin gauze dressing. *Int J Derm* 1998; 37: 390-396. 5. Sundberg J, Mellor R. A retrospective review of the use of cadexomer iodine in the treatment of wounds. *Wounds* 1997; 9(3): 68-86.



Moisture Responsive Intravenous Catheter Dressing


IV3000 is a moisture responsive, transparent film dressing, specifically designed to meet the needs of catheter fixation.

IV3000 is made from moisture responsive film, which is breathable and significantly more permeable to water vapour than ordinary film. The film is waterproof and impermeable to liquids and bacteria¹, and while the dressing remains intact it protects the IV site from outside contamination. These properties result in drier catheter sites, reduced skin maceration and lower skin colonisation, thus reducing the risk of catheter related infection.

IV3000 can be applied easily and aseptically, and due to the pattern spread adhesive, pain on removal is reduced.

Four ranges are available – Ported, Non-Ported, Frame Delivery, and Orange Handle. The dressings are available in a variety of shapes and sizes to meet required applications. Ported, Non-Ported and Frame Delivery dressings include two sterile securing strips for catheter or tube fixation and a sterile label documentation.

Ported




Non-Ported



Frame Delivery



Orange Handle



Code	Product description	Items per unit	Recommended indication
Ported			
66004011	5cm x 6cm	Box/100	Paediatric
4006	7cm x 9cm	Box/100	Peripheral
66004009	9cm x 12cm (oval)	Box/50	Central
66800512	11cm x 14cm	Box/25	PICC
Non-Ported (1-Hand)			
4007	6cm x 7cm	Box/100	Peripheral
4008	10cm x 12cm	Box/50	Central
Frame Delivery			
59410082	6cm x 7cm	Box/100	Peripheral
59410882	10cm x 12cm	Box/50	Central
Orange Handle			
4923	6cm x 8cm	Box/100	Peripheral (Ported)
4924	6cm x 8.5cm	Box/100	Peripheral
4925	10cm x 14cm	Box/10	Central
4973	10cm x 14cm	Box/50	Central
4649	10cm x 20cm	Box/50	PICC/Epidural

1. Hornabrook RA. The bacterial barrier properties of OPSITE IV 3000. Data on file; 1990, report MBP131.

JELONET[◇]

Paraffin Gauze Dressing


JELONET is a soft paraffin dressing made from open-weave gauze. The gauze has interlocking threads, which minimise fraying when the dressing is cut to shape.

JELONET is soothing and low-adherent, and allows the wound to drain freely into an absorbent secondary dressing.

Indications

- Minor burns or scalds
- Leg ulcers

Code	Product description	Items per unit
7403	5cm x 5cm	Box/50
7404	10cm x 10cm	Box/10
66007450	10cm x 10cm	Box/50
7409	10cm x 10cm	Box/100
7459	10cm x 40cm	Box/10
7415	15cm x 2m roll	Box/12



LEUKOSTRIP[◇]


Skin Closure Strips

LEUKOSTRIP is a hypoallergenic skin closure strip that keeps wound edges together securely. It is made of 100% polyamide material coated with a hypoallergenic adhesive.

Indications

LEUKOSTRIP is indicated for cuts and surgical incisions, and can help minimise scar formation when used to replace sutures or staples.¹

Code	Product description	Items per unit
White		
66002876	4mm x 38mm (8 strips/pouch)	Box/50
66002878	6.4mm x 76mm (3 strips/pouch)	Box/50
66002879	6.4mm x 102mm (5 strips/pouch)	Box/50
66002880	13mm x 102mm (6 strips/pouch)	Box/50
66002881	26mm x 102mm (4 strips/pouch)	Box/25
66002903	4mm x 76mm (4 strips/pouch)	Box/50
66002952	6.4mm x 76mm (3 strips/pouch)	Box/10*
66002953	13mm x 102mm (6 strips/pouch)	Box/10*
Tan		
66002882	4mm x 38mm (4 strips/pouch)	Box/50
66002883	6.4mm x 76mm (3 strips/pouch)	Box/50



*Not available in New Zealand

1. G Lepsien et al: Wound Healing Process and Cosmetic results after early removal of clamps and the replacement thereof by adhesive strips Akt. Chir 27 (1992).

MELOLIN[◇]

Low Adherent Absorbent Wound Dressing

MELOLIN is a low-adherent, absorbent dressing that absorbs excess fluid from moderately exuding wounds. It consists of three layers: a low-adherent perforated film, a highly absorbent cotton/polyester pad and a hydrophobic backing layer. This structure allows for drainage of exudate, reducing trauma to healing tissue.¹ It is easily conformable and retains its integrity when cut.

Indications


Light-to-moderately exuding wounds including:

- Clean sutured wounds
- Abrasions
- Lacerations
- Minor burns

Code	Product description	Items per unit
66974940	5cm x 5cm	Box/100
66974932	5cm x 5cm	Box/25*
66974933	10cm x 10cm	Box/10*
66974941	10cm x 10cm	Box/100
66974939	10cm x 20cm	Box/100
66974930	50cm x 7m Roll (non-sterile)	Roll/1

*Not available in New Zealand

1. Robb WAT, Clinical Trial of MELOLIN: a new non-adherent dressing. *Brit J. Plast. Surg* 1961; 15 47-49.



MELOLITE[◇]



Low Adherent Absorbent Pads

MELOLITE dressings are constructed from an absorbent pad which helps distribute the fluid, thus avoiding the problems of immediate exudate strike-through, covered on both sides with a heat-laminated polyethylene film. This film acts as a low-adherent wound contact layer by allowing the passage of exudate through its elliptical holed construction into the absorbent pad. This minimises the risk of adherence to the wound and the related trauma on dressing removal. Either side of the dressing can be placed directly onto the wound site.

Indications

- Clean, sutured wounds
- Lacerations
- Abrasions

Code	Product description	Items per unit
4811	7.5cm x 5cm	Box/100
4812	7.5cm x 10cm	Box/100
4814	7.5cm x 20cm	Box/100



OPSITE[◇] Range

Film dressings



OPSITE products are designed specifically to treat wounds in a range of healthcare settings. The transparent, adhesive polyurethane film helps protect wounds from bacterial contamination.¹


OPSITE[◇] FLEXIFIX[◇]

Transparent Film Roll

OPSITE FLEXIFIX combines all the advantages of OPSITE film with a simple application system, to allow easy and secure fixation of all kinds of dressings. It is resistant to water and bacteria, and is ideal for application on awkward areas of the body. Non-sterile.

Indications

- Provides secure fixation of primary dressings and tubes
- To protect the skin around stoma and under leg bags
- Aids in preventing bacterial contamination of the wound



Code	Product description	Items per unit
66000040	5cm x 10m roll	Roll/1
66000041	10cm x 10m roll	Roll/1
66000375	15cm x 10m roll	Roll/1*

*Not available in Australia

1. Data on file. Bacterial barrier properties of OPSITE Post-op film against methicillin-resistant MRSA. Data on file: report 0505004

OPSITE[◇] FLEXIFIX[◇] Gentle

Transparent Film Roll with Silicone Adhesive

OPSITE FLEXIFIX Gentle combines breathability, transparency, and up to a seven day wear time to offer a versatile silicone film roll. It removes cleanly without disrupting fragile skin layers or causing patients undue pain. Non-sterile.

Indications

- Indicated to retain primary wound dressings
- To protect the skin from friction and external contamination
- For patients with fragile or at risk skin

Code	Product description	Items per unit
66801195	2.5cm x 5m roll	Roll/1
66801196	5cm x 5m roll	Roll/1
66801197	10cm x 5m roll	Roll/1



OPSITE[◇] FLEXIGRID




Transparent film dressing

OPSITE FLEXIGRID is a sterile film dressing that facilitates moist wound healing for superficial wounds such as clean grazes and minor burns. The dressing can be left in place for up to fourteen days, minimising interference at the wound site and reducing the amount of nursing time required.

Indications

- The management of superficial wounds (e.g. superficial pressure injuries, minor burns, scalds, abrasions, lacerations and leg ulcers in the final stage of healing)
- The protection of skin from friction and external contamination
- The retention of primary dressings e.g. INTRASITE[◇] Gel, ALLEVYN[◇] Non-Adhesive
- The fixation of catheters
- The protection of skin around stoma and under leg bags
- The dressing of post-operative wounds



Code	Product description	Items per unit
4628	6cm x 7cm	Box/100
4629	10cm x 12cm	Box/10
4630	10cm x 12cm	Box/50
4631	15cm x 20cm	Box/10
4632	12cm x 25cm	Box/20

OPSITE[◇] Incise

Adhesive Incision Drape

OPSITE Incise Drapes are transparent, adhesive polyurethane films which allow the skin to breathe, so preventing moisture build-up under the drape. As a result OPSITE Incise Drapes adhere throughout long surgical operations to the surrounding skin and most importantly to the wound edge.



Code	Product description	Items per unit
4967	14cm x 25cm	Box/20*
4986	28cm x 15cm	Box/10
4987	28cm x 30cm	Box/10
4988	28cm x 45cm	Box/10
4995	42cm x 40cm	Box/10
4989	55cm x 45cm	Box/10
4994	84cm x 56cm	Box/10

*Not available in New Zealand

OPSITE[◇] Post-Op

Absorbent, Waterproof and Bacteria Proof Film Dressing


OPSITE Post-Op with its REACTIC film, offers a high level bacterial barrier and moisture management in a post-op dressing. The low-allergy adhesive reduces the risk of irritation while ensuring the dressing stays in place, even over awkward contours. The pad surface leaves the wound site clean and comfortable, reducing the pain and site trauma when changing the dressing.

Indications

To address acute wounds where levels of exudate are low/moderate such as:

- Lacerations
- Cuts
- Abrasions
- Post op wounds
- Minor burns

Code	Product description (dressing; pad)	Items per unit
66000708	6.5cm x 5cm; pad: 4cm x 2.5cm	Box/100
66000709	9.5cm x 8.5cm; pad: 7.3cm x 3.7cm	Box/20
66000710	12cm x 10cm; pad: 7.5cm x 5.2cm	Box/10
66000712	15.5cm x 8.5cm; pad: 12cm x 3.7cm	Box/20
66000713	20cm x 10cm; pad: 15cm x 5.6cm	Box/20
66000714	25cm x 10cm; pad: 20cm x 5.6cm	Box/20
66000715	30cm x 10cm; pad: 25cm x 5.6cm	Box/20
66000716	35cm x 10cm; pad: 30cm x 5.6cm	Box/20



OPSITE[◇] Post-Op Visible

Waterproof, Bacteria-proof Dressing with See-through Absorbent Pad

OPSITE Post-Op Visible is a dressing with a transparent REACTIC[◇] film and a see-through absorbent pad. The highly breathable moisture responsive film is waterproof and bacteria-resistant, and the unique lattice structure foam pad enables regular wound assessment without the need to lift or remove the dressing. OPSITE Post-Op Visible is ideal for moderately exuding wounds including invasive surgical wounds, lacerations and abrasions. The low allergy adhesive reduces the risk of irritation and keeps the dressing in place, even over awkward contours.



OPSITE Visible Drain Dressing is designed specifically for dressing drains and other percutaneous catheters.

Indications

To address acute wounds where levels of exudate are low/moderate such as:

- Lacerations
- Cuts
- Abrasions
- Post op wounds
- Minor burns

Code	Product description	Items per unit
66800136	8cm x 10cm	Box/20
66800137	15cm x 10cm	Box/20
66800138	20cm x 10cm	Box/20
66800139	25cm x 10cm	Box/20
66800140	30cm x 10cm	Box/20
66800141	35cm x 10cm	Box/20
OPSITE Visible Drain Dressing		
66800842	10cm x 9cm	Box/20



OPSITE[◇] Spray

Moisture Vapour Permeable Spray

OPSITE Spray is a water-resistant and quick-drying film. It is permeable to moisture vapour and air, allowing the surrounding skin to breathe and helping prevent maceration. OPSITE Spray is not sterilised.

Indications

Provides protection for dry, minor surgical and surface wounds.

Code	Product description	Items per unit
66004978	100ml Can	Car/12




PICO[◇] 7

Single Use Negative Pressure Wound Therapy System

PICO 7 Single Use Negative Pressure Wound Therapy System (sNPWT) technology is completely portable and clinically effective in the treatment of chronic; acute; traumatic; subacute and dehisced wounds; partial-thickness burns; ulcers (such as diabetic or pressure); flaps and grafts; closed surgical incisions. Suitable for use in both a hospital and home care setting.

PICO 7 delivers a number of enhancements which benefit both the clinical effectiveness and patient experience.

- Use uninterrupted for up to 7 days, potentially saving time and money by optimising dressing changes¹⁻³
- Twice as effective at dealing with air leaks, compared to the previous version of PICO sNPWT⁴
- Can be used in hard-to-seal or awkward places⁵
- Improved patient experience with a more efficient, quieter pump and integrated belt clip for full portability⁶⁻⁷
- A new 'dressing full' change indicator, designed to help optimise dressing change frequency
- Area to write start date of therapy, helping with healthcare protocols
- The PICO Dressing, with AIRLOCK[◇] Technology, distributes pressure evenly across the entire wound surface for optimal performance⁸
- The PICO dressing consists of four layers which allows for effective fluid management. The dressing evaporates up to 80% of fluid, whilst 20% is absorbed, without the need for a bulky canister⁸
- PICO 7Y is uniquely designed to manage two incisions simultaneously on one patient via an integral Y adaptor



Code	Product description	Items per unit
Two dressings and one pump kits		
66802002	10cm x 20cm	1 Kit
66802003	10cm x 30cm	1 Kit
66802004	10cm x 40cm	1 Kit
66802005	15cm x 15cm	1 Kit
66802006	15cm x 20cm	1 Kit
66802007	15cm x 30cm	1 Kit
66802008	20cm x 20cm	1 Kit
66802009	25cm x 25cm	1 Kit
66802000	Multisite 15cm x 20cm	1 Kit
66802001	Multisite 20cm x 25cm	1 Kit
66802031	7Y Multisite 20cm x 25cm	1 Kit
One dressing and one pump kits		
66802012	10cm x 20cm	1 Kit
66802013	10cm x 30cm	1 Kit
66802014	10cm x 40 cm	1 Kit
Accessories		
66801021	Foam Filler 10cm x 12cm x 1.5cm	Pack/5
66802127	Antimicrobial Gauze Rolls 11.4cm x 3.7m	Pack/5

PICO^o 14

Single Use Negative Pressure Wound Therapy System

PICO 14 Single Use Negative Pressure Wound Therapy System (sNPWT) is portable and clinically effective in the treatment of chronic; acute; traumatic; subacute and dehisced wounds; partial-thickness burns; ulcers (such as diabetic or pressure); flaps and grafts; closed surgical incisions. Suitable for use in both a hospital and home care setting.

PICO 14 sNPWT has been designed to challenge hard-to-heal and deep wounds.⁹ PICO 14 System uses the same innovative technology as PICO 7 with features including:

- A pump duration of up to 14 days¹⁰
- Designed for use on deep wounds⁹
- Use in hard-to-seal or awkward places⁵
- Dressing full change indicator and 'start date' note area helps support efforts to reduce the cost of unnecessary dressing changes and patient visits⁵
- The PICO Dressing, with AIRLOCK^o Technology, distributes pressure evenly across the entire wound surface for optimal performance⁸
- The PICO dressing consists of four layers which allows for effective fluid management. The dressing evaporates up to 80% of fluid, whilst 20% is absorbed, without the need for a bulky canister⁸
- Multipacks of 5 dressings now available, allowing therapy to be tailored to the patients clinical needs and provides flexibility in dressing choice



Code	Dressing sizes	Items per unit	Code	Dressing sizes	Items per unit
Two dressings and one pump kits			Multipack of 5 dressings		
66802042	10cm x 20cm	1/Kit	66802022	10cm x 20cm	Pack/5
66802043	10cm x 30cm	1/Kit	66802023	10cm x 30cm	Pack/5
66802044	10cm x 40cm	1/Kit	66802024	10cm x 40cm	Pack/5
66802045	15cm x 15cm	1/Kit	66802025	15cm x 15cm	Pack/5
66802046	15cm x 20cm	1/Kit	66802026	15cm x 20cm	Pack/5
66802047	15cm x 30cm	1/Kit	66802027	15cm x 30cm	Pack/5
66802048	20cm x 20cm	1/Kit	66802028	20cm x 20cm	Pack/5
66802049	25cm x 25cm	1/Kit	66802029	25cm x 25cm	Pack/5
66802040	Multisite 15cm x 20cm	1/Kit	66802020	Multisite 15cm x 20cm	Pack/5
66802041	Multisite 20cm x 25cm	1/Kit	66802021	Multisite 20cm x 25cm	Pack/5

Code	Product description	Items per unit
Accessories		
66801021	Foam Filler 10cm x 12cm x 1.5cm	Pack/5
66802127	Antimicrobial Gauze Rolls 11.4cm x 3.7m	Pack/5

1. Dowsett C et al. Use of PICO to improve clinical and economic outcomes in hard-to-heal wounds. Wounds Int. 2017; 8(2):52-58.
2. Hampton J. Providing cost-effective treatment of hard-to-heal wounds in the community through use of NPWT. Br J Community Nurs. 2015; 20:S14-S20.
3. Nherera LM et al. Cost-effectiveness analysis of single-use negative pressure wound therapy dressings (sNPWT) to reduce surgical site complications (SSC) in routine primary hip and knee replacements. Wound Repair & Regeneration 2017; 25(3):474-482. DOI:10.1111/wrr.12530.
4. Data on file, reference DS/17/666/R2: Comparison of PICO 1.6 and 2.1 device air leak tolerance. January 2018.
5. Data on file report OR-DOF/42. A volunteer trial to assess dressing performance of Palermo (PICO) multisite vs control. March 2014.
6. WMP.11446.UEF/R3 Project Fairbanks Human Factors Summary Report Issue 2. R. Knight, August 2016.
7. Data on file, reference DS/17/701/R2: Acoustic Testing Report: Comparison of PICO 7 to PICO 1.6 devices. January 2018.
8. Malmjö M. et al. Biological effects of a disposable, canisterless Negative Pressure Wound Therapy system. Eplasty 2014; 14:e15.
9. Smith & nephew 2018. The review of evidence supporting the use of PICO in wounds ≥2cm in depth. Internal report EO AWM. PCS230.001 v2.
10. Smith & Nephew December 2018. PICO 14 Service Life Testing: 14 Day Device Lifespan. Internal Report. RD/18/132.

PRIMABAND[◇]

Elastic Conforming Bandage

A retention bandage suitable for high-movement areas such as joints. Suitable for all types of patients, especially those with fragile or sensitive skin.

- Soft
- Porous, permitting air circulation
- White, conforming retention bandage
- Lightweight
- Non-adhesive
- Highly conformable

Indications

To hold dressings and other devices in place, such as:

- Synthetic casts
- Back slabs
- IV lines



Code	Product description	Items per unit
36361487	2.5cm x 1.75m unstretched	Pkt/12 Rolls
36361488	5cm x 1.75m unstretched	Pkt/12 Rolls
36361489	7.5cm x 1.75m unstretched	Pkt/12 Rolls
36361490	10cm x 1.75m unstretched	Pkt/12 Rolls
36361491	15cm x 1.75m unstretched	Pkt/12 Rolls

PRIMACREPE[◇]

Elastic Crepe Bandage (Medium Weight)

A retention bandage where light support is required.

- White, cotton crepe bandage
- Medium weight
- Air permeable
- Easy to apply
- Good skin tolerance

Indications

To provide firm directional support to joint injuries involving ligaments.



Code	Product description	Items per unit
36361478	5cm x 1.6m unstretched (3.5m stretched)	Pkt/12 Rolls
36361479	7.5cm x 1.6m unstretched (3.5m stretched)	Pkt/12 Rolls
36361480	10cm x 1.6m unstretched (3.5m stretched)	Pkt/12 Rolls
36361481	15cm x 1.6m unstretched (3.5m stretched)	Pkt/12 Rolls

PRIMACREPE[◇]

Elastic Crepe Bandage (Heavy Weight)

A retention bandage where firm directional support for joint injuries involving ligaments is required. Ongoing support for minor sprains and strains.

- Tan, cotton crepe bandage
- Heavy weight
- Centre guide for easy application
- Washable and reusable
- Excellent skin tolerance

Indications

To provide firm directional support to joint injuries involving ligaments.



Code	Product description	Items per unit
36361483	5cm x 2.3m unstretched (4.5m stretched)	Pkt/12 Rolls
36361484	7.5cm x 2.3m unstretched (4.5m stretched)	Pkt/12 Rolls
36361485	10cm x 2.3m unstretched (4.5m stretched)	Pkt/12 Rolls
36361486	15cm x 2.3m unstretched (4.5m stretched)	Pkt/12 Rolls

PRIMAFIX[◇] Plus

Conformable Retention Tape

Secure retention of primary dressings, particularly on joints and contoured areas of the body, allowing them to be dressed without restricting movement.

- Non-woven retention tape
- Air and water vapour permeable
- Low-allergy adhesive making it skin friendly
- Highly conformable
- Can be used under taping

Indications

Wide area fixation of wound dressings, particularly on:

- Joints
- Frequently moved and contoured parts of the body
- Fixation of tubes
- Catheters
- Measuring devices
- Instruments



Code	Product description	Items per unit
36361503	2.5cm x 10m	Box/1 Roll
36361504	5cm x 10m	Box/1 Roll
36361505	10cm x 10m	Box/1 Roll
36361506	15cm x 10m	Box/1 Roll
36361509	20cm x 10m	Box/1 Roll

PRIMAGAUZE[◇]

Elastic Cohesive Bandage

A cohesive retention bandage suited to areas which are highly mobile and difficult to dress.

- White, low bulk bandage
- Sticks to itself, not to skin, hair or clothing
- Highly conformable
- Air-permeable
- Stays in place, even on limbs with varying contours
- Skin friendly, suitable for fragile skin

Indications

To hold dressings and other devices in place, such as:

- Synthetic casts
- Back slabs
- IV lines



Code	Product description	Items per unit
36361492	2.5cm x 2m unstretched (4m stretched)	Box 2 Rolls
36361493	4cm x 2m unstretched (4m stretched)	Box/1 Roll
36361494	6cm x 2m unstretched (4m stretched)	Box/1 Roll
36361495	8cm x 2m unstretched (4m stretched)	Box/1 Roll
36361496	10cm x 2m unstretched (4m stretched)	Box/1 Roll
36361497	6cm x 10m unstretched (20m stretched)	Box/1 Roll
36361498	8cm x 10m unstretched (20m stretched)	Box/1 Roll
36361499	10cm x 10m unstretched (20m stretched)	Box/1 Roll

PRIMAPORE[◇]

Adhesive Non-Woven Wound Dressing



PRIMAPORE is a soft, conformable fabric island dressing that can be used to cover a variety of wounds. The highly absorbent pad provides an effective barrier, minimising fluid strike-through. PRIMAPORE also offers a low-adherent wound contact layer, reducing trauma to underlying tissue and ensuring pain is minimised on removal. The low-allergy acrylic adhesive provides secure fixation over the entire wound site.

Indications

It is ideal for use in treating wounds in curved areas of the body, including surgical wounds, general cuts, lacerations and abrasions.

Code	Product description (dressing; pad)	Items per unit
66003634	7.2cm x 5cm; pad: 4cm x 2.5cm	Box/50*
7135	8.3cm x 6cm; pad: 5.7cm x 3.5cm	Box/50
66000317	10cm x 8cm; pad: 5.3cm x 4.5cm	Box/20*
66000318	15cm x 8cm; pad: 11cm x 4cm	Box/20
66000319	20cm x 10cm; pad: 15.5cm x 5.3cm	Box/20
66000320	25cm x 10cm; pad: 20cm x 5.3cm	Box/20
66007140	35cm x 10cm; pad: 29.5cm x 5.3cm	Box/20*

*Not available in New Zealand



PRIMASWAB[◇]

Skin Cleansing Alcohol Swab

Simple and quick cleansing.

- Two-ply skin cleansing wipes
- Individually packed
- Saturated with 70% v/v isopropyl alcohol

Code	Product description	Items per unit
36361526	PRIMASWAB Alcohol Wipe	Box/200



PRIMATAPE[◇]

Universal Tape

A universal zinc oxide tape, suitable for securing dressings, particularly on joints and contoured parts of the body.

- Tan, rigid tape
- Strong adhesion
- Breathable
- Conformable

Indications

Intended for fixation of wound dressings, particularly on joints, frequently moved and contoured parts of the body. Secure fixation of:

- Tubes
- Catheters
- Measuring devices and instruments

May also be applied to cover and exert pressure on a body part.

Code	Product description	Items per unit
36361510	1.25cm x 5m	Box/12 Rolls
36361511	2.5cm x 5m	Box/12 Rolls
36361512	5cm x 5m	Box/6 Rolls
36361513	7.5cm x 5m	Box/6 Rolls



PRIMATAPE[◇]



Elastic Tape

An elastic zinc oxide tape, suitable for securing dressings, particularly on joints, finger taping and high movement areas.

- Tan, elastic tape
- Strong adhesion
- Breathable
- Flexible

Indications

Intended for fixation of wound dressings, particularly on joints, frequently moved and contoured parts of the body. Secure fixation of:

- Tubes
- Catheters
- Measuring devices and instruments



Code	Product description	Items per unit
36361518	2.5cm x 1m	Box/12 Rolls
36361519	5cm x 1m	Box/6 Rolls
36361520	2.5cm x 2.5m	Box/12 Rolls

PRIMATAPE[◇]

Waterproof Tape

A waterproof zinc oxide tape, suitable for securing dressings which may be exposed to wet and/or dirty environments.

- White, waterproof tape
- Strong adhesion
- Breathable
- Easy to tear

Indications

Intended for fixation of wound dressings, particularly on joints, frequently moved and contoured parts of the body. Secure fixation of:

- Tubes
- Catheters
- Measuring devices and instruments



Code	Product description	Items per unit
36361523	1.25cm x 5m	Box/12 Rolls
36361524	2.5cm x 5m	Box/12 Rolls
36361525	5cm x 5m	Box/6 Rolls

PROFORE[◇]

Multi-layer Compression Bandage System

PROFORE is a convenient presentation of the clinically effective four-layer bandaging system for venous leg ulcer management. It reinforces a safe and consistent protocol to heal up to 72% of venous leg ulcers in 12 weeks.¹

Proven effective in both hospital and community treatment studies, PROFORE is designed to deliver 40mm Hg pressure at the ankle, decreasing to 17mm Hg at the knee, therefore offering effective graduated compression.²

Indications

- Specifically developed for the management of venous leg ulcers and associated conditions
- This system pack is designed to give best results on ankle circumferences of 18-25cm

Code	Product description	Items per unit
66000016	Profore 4-layer kit (18-25cm)	1 Kit
66000770	Profore LATEX FREE 4-layer kit (18-25cm)	1 Kit
66000017	Profore 4-layer kit (25-30cm)	1 Kit
66000018	Profore 4-layer kit (over 30cm)	1 Kit

1. Moffatt, C., Simon, D., Franks, P., et al. Randomised trial comparing two fourlayer bandage systems in the management of chronic leg ulceration. *Phlebology*. 1999;14(4):139-142. 2. Blair, S.D., Wright, D., Backhouse, C.M., Riddle, E., Mccollum, C.N. Sustained compression and healing of chronic venous ulcers. *BMJ*. 1988;297(6657):1159-1161.



PROFORE[◇] Lite

Multi-layer Compression Bandage System

PROFORE Lite is a multi-layer compression bandage system offering reduced compression from the original PROFORE kits.

It has been specially developed for the management of “mixed” aetiology leg ulcers where assessment of the patient has identified that there is some degree of arterial impairment which prevents the use of full compression.

Code	Product description	Items per unit
66000415	PROFORE Lite 3-layer kit (18-28cm)	1 Kit





Skin Care

Prevention and treatment of incontinence-associated dermatitis (IAD) is possible using the PROSHIELD skincare range.

PROSHIELD Foam and Spray Incontinence Cleanser

A gentle, no-rinse, pH-balanced cleanser effective in the prevention and management of moisture lesions/IAD when used with PROSHIELD PLUS.

Indications

Cleansing of intact or injured skin associated with common severe or chronic diarrhoea or incontinence. Also indicated for removing urine and faeces, dried blood, post-operative antiseptic solutions and other hard to remove debris from the patient and for removing skin barriers such as PROSHIELD Plus Skin Protectant.

PROSHIELD PLUS Skin Protectant

26% reduction in the incidence of pressure injuries and IAD was seen in a study, when used as part of the PROSHIELD system, compared to standard care.¹

Indications

For the care of intact or injured skin associated with common severe or chronic diarrhea or incontinence.



Code	Product description	Items per unit
8213015008	PROSHIELD Foam & Spray 235ml	Each
8213030004	PROSHIELD Plus Skin Protectant 115g	Each

1. Wall L, Vernon T. Launch of a Skin Care Regime to reduce the incidence of Moisture Associated Skin Damage. Doncaster and Bassetlaw NHS Foundation Trust. Presented at TVS, UK, 2016.

REMOVE[◇]

Adhesive remover

REMOVE is a gentle solvent used to remove adhesive dressings and tapes.

REMOVE is aloe-formulated to moisturise and condition tender or friable skin. It is easy to apply and cost-effective, as the single-use wipes contain just the right amount of adhesive remover, so there is no mess or waste.

Indications

It dissolves adhesives and assists in removing acrylic-based, rubber-based and hydrocolloid-based residues from the skin without irritation or trauma.



Code	Product description	Items per unit
403100	REMOVE Wipe	Box/50

RENASYS[◇]

Traditional Negative Pressure Wound Therapy



The RENASYS range of Traditional NPWT devices provides clinicians with a range of flexible, powerful and portable devices, canisters and wound fillers to manage a variety of wound types, in any care settings.

NPWT consists of a canister, a wound filler (foam or gauze), a drainage channel and a transparent film that are all connected to a suction device.

Indications

For patients who would benefit from a suction device (Negative Pressure Wound Therapy), as it may promote wound healing by removing fluid, including wound exudate and infectious materials.


Examples of appropriate wound types include:

- Chronic
- Acute
- Traumatic
- Sub-acute and dehisced wounds
- Pressure injuries
- Partial-thickness burns
- Flaps and grafts
- Diabetic foot ulcers

RENASYS[◇] TOUCH

RENASYS TOUCH combines an intuitive design with power and flexibility, facilitating effective therapy whilst supporting individualised patient care from the hospital to the home.

- **Intuitive** – Easy to navigate and a high resolution touch screen provides users with a simple and easy interface
- **Flexible** – Designed for portability, flexibility and simplicity
- **Powerful** – Manages a variety of wound types and prescribed therapies, all in one device, reducing the need for multiple devices and canisters



Code	Product description	Unit of measure	Min. order qty
66802134	RENASYS TOUCH device	Each	1
66801273	RENASYS TOUCH 300ml canister with solidifier	Each	5
66801274	RENASYS TOUCH 800ml canister with solidifier	Each	5
66801276	RENASYS TOUCH carry straps	Each	1
66801277	RENASYS TOUCH carry bag	Each	1
66801278	RENASYS TOUCH IV/Bed pole clamp	Each	1
66801279	RENASYS TOUCH Transit case	Each	1
66801286	RENASYS TOUCH Class 2 Power Supply	Each	1
66801560	RENASYS Power Cord	Each	1
66801284	RENASYS TOUCH Odour Filter	Pack/20 units	1
66801283	RENASYS TOUCH O-Ring	Pack/20 units	1

Foam Dressing Kits



RENASYS^o-F

Foam Dressing Kit – Small

Component list

- 1 - Soft Port – length 69cm applicator head 15cm x 10cm
- 1 - Foam block 10cm x 8cm x 3cm
- 1 - RENASYS Transparent Film 20cm x 30cm



RENASYS-F

Foam Dressing Kit – Medium

Component list

- 1 - Soft Port – length 69cm applicator head 15cm x 10cm
- 1 - Foam block 20cm x 13cm x 3cm
- 2 - RENASYS Transparent Film 20cm x 30cm



RENASYS-F

Foam Dressing Kit – Large

Component list

- 1 - Soft Port – length 69cm applicator head 15cm x 10cm
- 1 - Foam block 25cm x 15cm x 3cm
- 3 - RENASYS Transparent Film 20cm x 30cm



RENASYS-F

Foam Dressing Kit – X-Large

Component list

- 1 - Soft Port – length 69cm applicator head 15cm x 10cm
- 1 - Foam block 50cm x 63cm x 1.5cm
- 6 - RENASYS Transparent Film 20cm x 30cm



RENASYS-F AB

Abdominal Dressing Kit

Component list

- 1 - Soft Port – length 69cm applicator head 15cm x 10cm
- 2 - Foam blocks 43cm x 30cm x 3cm
- 6 - RENASYS Transparent Film 20cm x 30cm
- 1 - organ protection layer 89cm x 66cm



Code	Product description	Unit of measure	Min. order qty
66800794	RENASYS-F Foam Dressing Kit with Soft Port – Small	1 Dressing Kit	5
66800795	RENASYS-F Foam Dressing Kit with Soft Port – Medium	1 Dressing Kit	5
66800796	RENASYS-F Foam Dressing Kit with Soft Port – Large	1 Dressing Kit	5
66800797	RENASYS-F Foam Dressing Kit with Soft Port – X-Large	1 Dressing Kit	5
66800980	RENASYS AB Abdominal Dressing Kit with Soft Port	1 Dressing Kit	1

Gauze Dressing Kits



RENASYS[◇]-G

Gauze Dressing Kit – Small

Component list

- 1 - saline bullet 30mL
- 1 - wound ruler
- 1 - sachet of SECURA[◇] No-Sting Barrier Film
- 1 - Soft Port – length 69cm applicator head 15cm x 10cm
- 1 - Kerlix AMD[™] Gauze Dressing Pack 15cm X 17cm
- 1 - RENASYS Transparent Film 20cm x 30cm
- 1 - non-adherent layer 8cm x 30cm



RENASYS-G

Gauze Dressing Kit – Medium

Component list

- 1 - saline bullet 30mL
- 1 - wound ruler
- 1 - sachet of SECURA No-Sting Barrier Film
- 1 - Soft Port – length 69cm applicator head 15cm x 10cm
- 2 - Kerlix AMD[™] Gauze Dressing Packs 15cm X 17cm
- 1 - RENASYS Transparent Film 20cm x 30cm
- 1 - non-adherent layer 8cm x 30cm



RENASYS-G

Gauze Dressing Kit – Large

Component list

- 1 - saline bullet 30mL
- 1 - wound ruler
- 1 - sachet of SECURA No-Sting Barrier Film
- 1 - Soft Port – length 69cm applicator head 15cm x 10cm
- 1 - Kerlix AMD[™] Gauze Dressing Roll 11cm x 4m
- 2 - RENASYS Transparent Film 20cm x 30cm
- 2 - non-adherent layers 8cm x 30cm



RENASYS[◇]-G

Gauze Dressing Kit – X-Large

Component list

- 2 - saline bullet 30mL
- 1 - wound ruler
- 2 - sachet of SECURA[◇] No-Sting Barrier Film
- 1 - Soft Port – length 69cm applicator head 15cm x 10cm
- 2 - Kerlix AMD[™] Gauze Dressing Rolls 11cm x 4m
- 3 - RENASYS Transparent Film 20cm x 30cm
- 4 - non-adherent layers 8cm x 20cm



Code	Product description	Unit of measure	Min. order qty
66800933	RENASYS-G Dressing Kit with Soft Port – Small	1 Dressing Kit	5
66800934	RENASYS-G Dressing Kit with Soft Port – Medium	1 Dressing Kit	5
66800935	RENASYS-G Dressing Kit with Soft Port – Large	1 Dressing Kit	5
66800936	RENASYS-G Dressing Kit with Soft Port – X-Large	1 Dressing Kit	5

Kit Accessories



RENASYS[◇] Transparent Film

Component list

- 10 - Transparent Films per kit 20cm x 30cm



RENASYS Y-Connector

Component list

- 1 - Y-Connector per kit



RENASYS Soft Port Stand-Alone Kit

Component list

- 1 - Soft Port per kit Length 69cm applicator head 15cm x 10cm



RENASYS 10mm Flat Drain Pack

Component list

- 1 - 10mm Flat Drain
- 1 - RENASYS Quick-Click Adapter
- 2 - non-Adherent Layer 7.5x20cm
- 1 - AMD gauze roll 11cm x 4m
- 2 - RENASYS Transparent Film
- 4 - strips waterproof tape
- 1 - RENASYS Adhesive Gel Patch
- 2 - sachet of SECURA No-Sting Barrier Film
- 1 - saline bullet 30mL
- 1 - wound ruler



RENASYS[◇] 15FR Channel Drain Pack

Component list


- 1 - 15 French Channel Drain
- 1 - RENASYS Quick-Click Adapter
- 2 - non-Adherent Layer 7.5x20cm
- 1 - AMD gauze roll 11cm x 4m
- 2 - RENASYS Transparent Film
- 4 - strips waterproof tape
- 1 - RENASYS Adhesive Gel Patch
- 2 - sachet of SECURA[◇] No-Sting Barrier Film
- 1 - saline bullet 30mL
- 1 - wound ruler



RENASYS Drains

Component list


- 1 - Drain (available in 2 configurations listed)
- 1 - RENASYS Quick-Click Adapter
- 2 - Strips waterproof tape
- 1 - RENASYS Adhesive Gel Patch



RENASYS Adhesive Gel Patch

Component list

- 10 gel patches per box 10cm x 7cm
- Double sided silicone adhesive hydrogel



RENASYS Accessories

Code	Product description	Unit of measure	Min. order qty
66802127	Antimicrobial Gauze Rolls	Pack/5	1
66800394	RENASYS Transparent Film	Box/10	1
66800853	RENASYS Transparent Film X-Large	Box/5	1
66800971	RENASYS Y-Connector	Each	10
66800799	RENASYS Soft Port Stand-Alone Kit	Each	5
66802252	RENASYS 10mm Flat Drain	Each	10
66802253	RENASYS 15FR Channel Drain	Each	10
66802256	RENASYS-G 10mm Flat Drain Pack*	Each	5
66801082	RENASYS Adhesive Gel Patch	Box/10	1

*Contains dressing kit with gauze and drain

SECURA[◇]

No sting barrier film skin protectant

SECURA No Sting Barrier Film is a liquid, film-forming, alcohol free skin protectant that minimises stinging. With fast drying times of less than 30 seconds¹⁻², SECURA No Sting Barrier Film will create an effective water proof barrier for up to 96 hours (four days).³

Indications

It may be used to protect intact and damaged skin from irritation due to:

- Urinary or faecal incontinence
- Bodily fluids
- Adhesive trauma
- Friction
- The product may be used on adults, children and full term infants from one month

Code	Product description	Items per unit	Cartons per case
66800710	Spray 28ml	Box/1	12
66800711	Swab 3ml	Box/50	12
66800712	Wipe 1ml	Box/50	20

SECURA[◇]

No Rinse Cleanser

SECURA No Rinse Cleanser is a liquid solution that is an alternative to soap and water. SECURA No Rinse Cleanser should be used in order to maintain the normal pH and moisture content of the skin, provide a no-rinse and easy one step cleansing solution.


SECURA No Rinse Cleanser is pH balanced to protect the skin's acid mantle.

Indications

- Aids in removal of urine and faeces, or other foreign material without friction
- Reduces odour
- Non-irritating

Code	Product description	Items per unit
59430979	SECURA No Rinse Cleanser 236ml	Each

1. Time-To-Dry Test: No-Sting Barrier Film. SAP, 2009. 2. Tactile Characterization (Stickiness) Test: No-Sting Barrier Film. SAP, 2009. 3. Study to Compare the Wash-Off Resistance and Durability of a Barrier Film. Rpt # 09-58A, Cyberderm Clinical Studies, 2009.





Moisturisers

SECURA moisturising creams and lotions are used to keep the skin intact and supple by preventing water loss. Well-hydrated skin and good skin integrity is important as the skin provides the first line of defence for the immune system.

Indications

- SECURA Moisturising Cream restores moisture to very dry skin, conditions, soothes and rehydrates whilst being absorbed quickly
- SECURA Moisturising Lotion should be used daily to moisturise dry skin in order to maintain healthy skin/acid mantle, is fragrance free and absorbs quickly



Code	Product description	Items per unit
59432000	SECURA Moisturising Cream 184g	Each
59433400	SECURA Moisturising Lotion 236ml	Each

S+N Skin Care Lotion

Moisturising Lotion

Effective but gentle, this non-greasy penetrating Skin Care Lotion contains 1.5% dimethicone, has humectants to moisturise the skin and effective emollients to relieve dry skin. While softening the skin, it provides a silicone barrier for lasting protection against the harsh conditions encountered by skin every day. Available in a handy 500ml pack, it dispenses easily and fits neatly on walls.

Indications

Protects the hands and the skin from chapping, dryness, and the effects of detergent alkalinity.



Code	Product description	Items per unit
HSCL500	500ml Pump Pack	Box/12 Dispensers

SKIN-PREP[◇]


Protective Wipe

SKIN-PREP is a liquid preparation that when applied to intact skin, forms a protective film. Only one coat is required.

SKIN-PREP can be easily applied to even the most awkward areas. It moves naturally with the patient's skin, so it will not crack or peel. The protective coating helps to preserve proper skin integrity and prevent damage or injury during the removal of tapes and films, making adhesive removal less painful. SKIN-PREP itself is also easily removed using skin cleanser or soap and water.

Indications

Prepares the skin for attachments such as drainage tubes, external catheters and adhesive dressings.



Code	Product description	Items per unit
420200	SKIN-PREP Spray 118ml	1 Bottle*
420400	SKIN-PREP Wipe	Box/50

*Not available in New Zealand

SOLOSITE[◇] Gel

Cooling, Soothing, Hydrating Gel

SOLOSITE Gel is a non-sterile, preserved hydrogel capable of donating and absorbing moisture.

Indications

Used to create an environment that supports moist wound healing for the treatment of minor conditions such as:

- Minor burns, including sunburn
- Superficial cuts and lacerations
- Skin tears and abrasions

Code	Product description	Items per unit
36100614	20g Tube	Box/10
36361354	50g Tube	Box/12
449600	3oz/85g	Box/1*

*Not available in Australia



UNI-SOLVE[◇]

Adhesive Remover

UNI-SOLVE thoroughly dissolve adhesives reducing trauma to the skin caused by frequent or aggressive tape adhesion.

Indications

It is ideal for removing adhesive dressings, ostomy appliances, adhesive orthopaedic strapping and Montgomery straps.

Code	Product description	Items per unit
402300	UNI-SOLVE Adhesive Remover Wipes	Box/50



VERSAJET[◇] II


Hydrosurgical Debridement

The VERSAJET II Hydrosurgery System uses a razor-thin saline jet to optimise surgical debridement. The system enables the surgeon to precisely select, excise and evacuate non-viable tissue, bacteria and contaminants from wounds, burns and soft tissue injuries using a tissue-sparing technique.

The VERSAJET PLUS handset enables even faster removal of more non-viable tissue than the VERSAJET EXACT. With the increased speed and strength of the saline jet, a higher deck height and wider channel within the operating window, the PLUS handset performs more powerfully and should therefore, only be used within the operating theatre.

Indications

- Precision: Precisely removes unwanted tissue whilst sparing viable tissue³
- Control: Grasps, cuts and removes unwanted tissue with one instrument
- Versatility: Effectively excises/debrides a wide spectrum of acute and chronic wounds, burns and surgical sites




Code	Product description	Items per unit
VERSAJET II Hydrosurgery System		
66800039	Console (includes foot switch and power cord)	Box/1
66800979	Surgical cart	Each
66800472	Multi-function footswitch	Each
VERSAJET II Hydrosurgery System Handsets		
66800041	VERSAJET II EXACT disposable handset 45°/14mm	Each
66800042	VERSAJET II EXACT disposable handset 45°/8mm	Each
66800044	VERSAJET II PLUS disposable handset 45°/14mm	Each

1. Hyland EJ, D'cruz R, Menon S, et al. Prospective, randomised controlled trial comparing Versajet hydrosurgery and conventional debridement of partial thickness paediatric burns. *BURNS*. 2015;41(4):700-707. 2. Matsumura H, Nozaki M, Watanabe K, et al. The Estimation of Tissue Loss During Tangential Hydrosurgical Debridement. *Annals of Plastic Surgery*. 2012;69(5):521-525. 3. Cubison TCS, et al. Dermal preservation using the Versajet hydrosurgery system for debridement of paediatric burns. *Burns* 2006; 32: 714-720.

VERSAJET II EXACT handsets

VERSAJET II PLUS handsets



Adhesive Remover

ZOFF is a universal adhesive remover ideal for the painless removal of adhesive tapes, dressings and acrylics such as hydrocolloids. It contains no trichloroethane and is easy to use.

Use

To clean all types of adhesive residue from the skin, including hydrocolloids, acrylics and rubber based adhesives.



Code	Product description	Items per unit
36105120	ZOFF Adhesive Remover 250ml	Box/12*

*Not available in New Zealand

Training and education from Smith+Nephew

From new techniques and products, to clinical evidence and health economics, Smith+Nephew's continuing education and professional development programs keep you up-to-date with the latest in wound management.

CLOSER TO ZERO^o Online Education

CLOSER TO ZERO Online Education is an exciting new platform for healthcare professionals to access peer-to-peer educational speed learning and symposiums. This is a free subscription-based service whereby new content is released regularly and available for on-demand viewing so you can decide when and what you want to watch and listen - the choice is entirely up to you! It is optimised to work on computers, tablets and smartphones.

Subscribe today at closertozero.com.au

Wound Education Program

The Smith+Nephew Wound Education Program provides comprehensive education for healthcare professionals on best practice in wound prevention, assessment and management. This program covers the fundamentals of wound management including:

- Threats to skin integrity
- Pressure injures
- Wound healing physiology
- Holistic approaches to assessment and management
- T.I.M.E Clinical Decision Support Tool
- Chronic wounds
- Infection management
- Product application, tips and tricks
- Case studies

For further information, contact your Smith+Nephew representative.

Customised Education Sessions

Smith+Nephew provides a comprehensive range of in-service options based on customer need. This degree of flexibility is designed to optimise educational support so that resources supporting product use are congruent with the dynamics of the diverse health environment.



Educational Events

In conjunction with leading experts, Smith+Nephew develops, organises and delivers educational events focused on specific indication areas. For example:

- Forums/ Expert meetings
- Seminars
- Round-table discussions
- Webinars

Education and Evidence

Smith+Nephew's free online resource 'Education and Evidence', offers a learning environment for wound care professionals who are interested in practical resources to support their professional development. Trusted by healthcare professionals around the world, this website is the place to go for:

- On-line courses
- Video presentations from leading world experts
- Literature summaries

Subscribe today at smith-nephew.com/education*

*Education and Evidence website is designed for the education of healthcare professionals on the subject matters therein contained, and for the global exchange of professional knowledge, scientific research results and practical medical information among such professionals.

The website may contain information and discussion (including the promotion of) methods, procedures or products that may not be available in certain countries or regions, or may be available under various other trade or service marks, names or brands.





Wound Management Resources



Skin Tear Management with ISTAP* Classification Tool

I Identify

- Correctly identify on first presentation.

S Stop the bleeding and clean

- Control superficial bleeding.
- Select appropriate cleanser.
- Gently clean wound bed.

T Tissue alignment

- Align skin flap and classify according to ISTAP classification system (measure and document).

A Assess and dress

- Use TIME Clinical Decision Support Tool to set goals of treatment (refer to TIME CDST).
- Mark the dressing with an arrow to indicate correct direction of removal.

P Prevent recurrence

- Consider limb protector (CoolFAST).
- Consider limb support (CoolGRIP) if clinically indicated.
- Use peri-wound protection and moisturise 2x/day with SECURA Moisturising Lotion.
- Implement Risk Reduction tool from ISTAP guidelines.

Product Solutions

Skin Tear Kit to promote best practice.
ALGISITE[®] M to assist with control of bleeding.



Dress with **ALLEVYN[®] LIFE, Gentle Border, Gentle Border Lite or Gentle Border Ag.**
SECURA[®] No-Sting Barrier film, Skin-Prep[®] for peri-wound protection.

CoolFAST[™]
CoolGRIP[™]
SECURA Moisturising Lotion
SECURA Moisturising Cream




Type 1
No skin loss. Flap can be repositioned.



Type 2
Partial flap loss. Flap will not fully cover wound bed.











Type 3
Total flap loss. Exposing entire wound bed.



*International Skin Tear Advisory Panel. Adapted from LeBlanc K et al. Best practice recommendations for the prevention and management of skin tears in aged skin. Wounds International 2018. Used with permission. Download from www.woundsinternational.com. Products listed and outline of care are examples only. Product selection and management should always be based on comprehensive clinical assessment. Smith & Nephew Pty Ltd T +61 2 9857 3999 www.smith-nephew.com/australia Smith & Nephew Ltd +64 9 820 2840 www.smith-nephew.com/new-zealand. Trademark of Smith & Nephew SN14693 (10/19).

+ IAD differentiation, classification and treatment

	Incontinence-Associated Dermatitis ^{1,3}	Pressure Injury ^{1,3}
Cause	Moisture (urinary, faecal, urinary-faecal)	Pressure/Ischaemia
Location	 Perineal* (buttocks, natal cleft, perineum, perianal)	 Over bone
Depth	 Superficial	 Superficial to deep
Tissue colour	 Pink/red *	 Pink/red, black, yellow
Edges/Shape	 Diffuse/wandering/irregular*	 Well defined
Associated factors	Urinary and/or faecal incontinence	Reduced mobility/Sensory impairment

*Images courtesy of Professor Dimitri Beekman, Ghent University, Belgium

IAD Categorisation tool²

Category 1A	Persistent redness without signs of infection. A variety of tones may be present (could be pale or purple in dark skin)
Category 1B	Persistent redness WITH clinical signs of infection - such as white scaling of skin or satellite lesions
Category 2A	Skin loss without clinical signs of infection - such as excoriation, denudation or skin erosion
Category 2B	Skin loss WITH clinical signs of infection - white scaling skin/slough may be visible, shiny wound with excess exudate, green appearance may suggest Pseudomonas

Choosing the right product

	Healthy skin with exposure to urine/faeces	Category 1A	Category 2A	Category 1B/2B
Skin condition				
Product selection				
Intervention³	<ul style="list-style-type: none"> Use PROSHIELD[®] Foam & Spray Cleanser to gently clean soiled area. Apply PROSHIELD Plus Skin Protectant. Avoid rubbing or friction. Reapply as often as necessary. 	<ul style="list-style-type: none"> Use PROSHIELD Foam & Spray Cleanser to clean soiled area. Apply PROSHIELD Plus Skin Protectant. Re-apply after each incontinence episode or as often as necessary. 	<ul style="list-style-type: none"> Use PROSHIELD Foam & Spray Cleanser to clean soiled area. Apply PROSHIELD Plus Skin Protectant with a gentle patting motion, or SECURA[®] No-Sting Barrier Film Spray. Re-apply as often as necessary. 	<ul style="list-style-type: none"> Use PROSHIELD Foam & Spray Cleanser to clean soiled area. Apply appropriate treatment for identified infection. Once resolved treat as per 1A or 2A as appropriate.

Biofilm management protocol

Cadexomer iodine (IODOSORB[®]) is indicated for the management of sloughy and non-healing chronic ulcers such as diabetes-related foot ulcers, venous leg ulcers, pressure injuries and other wounds healing by secondary intention.


IODOSORB has several modes of action:



1. Disrupts biofilm¹⁻⁴



2. Reduces bioburden^{1,5-7}



3. Absorbs exudate^{5,8-12}



4. Removes slough & necrotic tissue^{5,8-10}

Clean wound as per local protocol with non-antimicrobial cleanser. Do not dry wound.

Apply Cadexomer Iodine (IODOSORB) to wound:

Ointment

- Apply ointment, making sure there is enough to cover whole wound with 3mm depth
- Use glove to spread out Cadexomer Iodine (IODOSORB) to ensure coverage of whole wound
- Cover with secondary dressing. If wound is heavily exuding use foam dressing



Powder

- Apply cadexomer iodine (IODOSORB) to surface of wound to depth of at least 3mm.
- Cover with secondary dressing. If wound is heavily exuding use foam dressing



Change Cadexomer Iodine (IODOSORB) every 2-3 days OR when colour has changed.

- If necessary soak dressing for several minutes with saline to loosen.
- Remove using stream of sterile water or saline



Reassess after two weeks

- If wound has started healing, but slough still present, continue Cadexomer Iodine (IODOSORB use) for a further two weeks, then reassess.

- If no apparent healing, discontinue and reassess patient to ensure correct diagnosis and all comorbidities controlled.

- If clean granulation tissue is present change to appropriate non-antimicrobial dressing (e.g. foam). If risk of further infection persists, consider the use of nanocrystalline silver (ACTICOAT[™]) to prevent biofilm reforming.

References

1. Phillips PL, et al. Antimicrobial dressing efficacy against mature Pseudomonas aeruginosa biofilm on porcine skin explants. Int Wound J 2015; 12: 469-483.
2. Akiyama H, et al. Assessment of cadexomer iodine against Staphylococcus aureus biofilm in vivo and in vitro using confocal laser scanning microscopy. J Dermatol 2004; 31:529-534.
3. Lumb H. Testing the biofilm disruption activity of IODOFLEX dressing. Data on file; 2007: report TSG015-07-001.
4. Lumb H. Testing the biofilm prevention activity of IODOFLEX dressing. Data on file; 2007: report TSG015-07-001.
5. Sundberg J, Meller R. A retrospective review of the use of cadexomer iodine in the treatment of wounds. Wounds 1997; 9(3): 68-86.
6. Salman H, Leakey A. The in vitro activity of silver sulphadiazine and cadexomer iodine against recent clinical isolates of methicillin-resistant Staphylococcus aureus, methicillin-resistant coagulase-negative Staphylococci and Pseudomonas aeruginosa. Data on file; 2001: report 194-03-01.
7. Zhou LH, et al. Slow release iodine preparation and wound healing: in vitro effects consistent with lack of in vivo toxicity in human chronic wounds. Brit J Derm 2002;146: 365-374.
8. Troeng T, et al. A randomized multicenter trial to compare the efficacy of cadexomer iodine and standard treatment in the management of chronic venous ulcers in out-patients. In: Fox JA, Fischer H, eds. Cadexomer iodine. New York: F.K. Schattauer Verlag; 1983:43-50.
9. Ormiston MC, et al. Controlled trial of Iodosorb in chronic venous ulcers. BMJ 1985; 291: 308-310.
10. Hansson C, et al. The effects of cadexomer iodine paste in the treatment of venous ulcers compared with hydrocolloid dressing and paraffin gauze dressing. Int J Dermatol 1998; 37: 390-396.
11. Lindsay G, et al. A study in general practice of the efficacy of cadexomer iodine in venous leg ulcers treated on alternative days. Acta Therapeutica 1986; 12:141-147.
12. Skog E, et al. A randomized trial comparing cadexomer iodine and standard treatment in the out-patient management of chronic venous ulcers. Brit J Derm 1983; 109: 77-83.

A route to more effective infection management

Improve patient outcomes¹ with accurate decision making, a fast response and effective treatment choices

T.I.M.E.
Clinical decision support tool

A
B
C
D
E

Start with following steps to undertake a comprehensive assessment²

A Assess patient, wellbeing and wound

B Bring in a multi-disciplinary team and informal carers to promote holistic patient assessment


C Control and treat the underlying causes and barriers to wound healing

D Decide appropriate treatment

E Evaluate and reassess the treatment and wound management outcomes

What clinical signs and symptoms of infection are present?

<div style="display: flex; justify-content: space-around;"> </div> <div style="display: flex; justify-content: space-around; margin-top: 5px;"> <div style="border: 1px solid black; padding: 5px; width: 45%;"> <p style="text-align: center; background-color: #008080; color: white; margin: -5px -5px 5px -5px;">Biofilm^{1,3-5}</p> <ul style="list-style-type: none"> Antibiotic/antimicrobial treatment failure Recurrence of delayed healing on cessation of antibiotic treatment Delayed healing despite optimal wound/patient management Low level chronic inflammation Low level erythema Friable granulation Covert (subtle) signs of infection </div> <div style="border: 1px solid black; padding: 5px; width: 45%;"> <p style="text-align: center; background-color: #008080; color: white; margin: -5px -5px 5px -5px;">Covert (subtle)^{1,3}</p> <ul style="list-style-type: none"> Delayed wound healing Serous drainage with concurrent inflammation Hypergranulation Bleeding, friable granulation Epithelial bridging and pocketing in granulation tissue Wound breakdown & enlargement New or increasing pain Increasing malodour </div> </div>	<div style="display: flex; justify-content: space-around;"> </div> <div style="display: flex; justify-content: space-around; margin-top: 5px;"> <div style="border: 1px solid black; padding: 5px; width: 45%;"> <p style="text-align: center; background-color: #800080; color: white; margin: -5px -5px 5px -5px;">Overt (classic)^{1,3}</p> <ul style="list-style-type: none"> Erythema Warmth Oedema/swelling Purulent discharge Pain Increasing malodour Delayed wound healing </div> <div style="border: 1px solid black; padding: 5px; width: 45%;"> <p style="text-align: center; background-color: #800080; color: white; margin: -5px -5px 5px -5px;">Spreading or systemic infection^{1,3}</p> <ul style="list-style-type: none"> Spreading erythema, warmth May include cellulitis, crepitus Wound breakdown/dehiscence with or without satellite lesions Malaise/lethargy Loss of appetite Systemic inflammatory response Sepsis Organ dysfunction </div> </div>
--	--



Use standard wound care (i.e. non-antimicrobial dressings) or advanced therapies until healing (follow local protocol)⁵

TWO-WEEK CHALLENGE^{1,6,9}

Antimicrobial dressings are recommended to be used for a minimum of two weeks' duration. After two weeks, re-evaluate and either:

1. discontinue if signs and symptoms of infection have resolved,
2. continue with antimicrobial if wound is progressing but there are still signs and symptoms, or
3. consider an alternative antimicrobial and refer to an appropriate specialist if no improvement.



* No one sign or symptom can reliably confirm the presence of infection, and those with immunosuppression may not exhibit signs and symptoms of clinical infection.
 † Cleanse wound and periwound skin thoroughly. Should an antiseptic cleanser be selected, the product's Instructions for Use (IFU) and soak time should be followed.
 ‡ Consider the use of DURAFIBER[®] Ag Silver Gelling Fibre Dressing for deep infected wounds.
 Ω Unless iodine contraindicated.
 ∞ For very-high risk patients and wounds (e.g. osteomyelitis), it may be appropriate to use antimicrobial treatment for longer than the two-week challenge.
 For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's Instructions for Use (IFU).

References 1. International Wound Infection Institute (IWII) Wound infection in clinical practice. Wounds International (2016). 2. Moore Z, et al. J Wound Care 28(3):154-161 (2019). 3. Weir D, Schultz G. Assessment and Management of Wound-Related Infections. In Doughty D & McNichol L (Eds.). Wound, Ostomy and Continence Nurses Society Core Curriculum: Wound Management (p. 156-180). 2016. Philadelphia: Wolters-Kluwer. 4. Wolcott RD, et al. J Wound Care 19(2): 45-53 (2010). 5. Schultz G, et al. Wound Repair Regen 25(5): 744-757 (2017). 6. Ayello EA, et al. Wounds Int 1-24 (2012). 7. Roche ED, et al. Int Wound J 1-10 (2019). 8. Malone M, et al. J Antimicrob Chemother 72, 2093-2101 (2017). 9. Schwarzer S, et al J Infect 80(3):261-270 (2020). 10. Gago M, Garcia F, Gaztelu V, Verdu J, Lopez P, Nolasco A. A comparison of three silver-containing dressings in the treatment of infected, chronic wounds. Wounds. 2008;20(10):273-278.
 Photographs (from left to right) courtesy of Kerlyn Carville, Kevin Woo, and Henri Post.
 ©Trademark of Smith+Nephew. All Trademarks acknowledged. ©June 2020 Smith+Nephew. 25993 | GMC1067

PICO sNPWT hard-to-heal wounds pathway⁵

Patient selection


Identify appropriate patients



*Wounds with overt signs of clinical infection (e.g. increased pain, levels of exudate, cellulitis etc.) should be excluded from the evaluation. Colonised/critically colonised wounds are not excluded from the evaluation. Site standard protocol should be implemented to address bacterial burden; †Wounds that have healed by <10% but have shown significant improvement in granulation tissue quality/ quantity may be considered for further PICO treatment based on clinician judgement.

ABPI: Ankle-brachial pressure index **VLU:** Venous leg ulcers

Smith & Nephew does not provide medical advice. The information presented is not, and is not intended to serve as, medical advice. It is the responsibility of healthcare professionals to determine and utilise the appropriate products and techniques according to their own clinical judgment for each of their patients. For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.



Smith & Nephew Pty Ltd
Australia
T +61 2 9857 3999
F +61 2 9857 3900

Smith & Nephew Ltd
New Zealand
T +64 9 820 2840
F +64 9 820 2841

www.smith-nephew.com/australia
www.smith-nephew.com/new-zealand

[®]Trademark of Smith+Nephew
All Trademarks acknowledged
©2019 Smith+Nephew
18456 SN14817 REVO 03/20

1. Troxler M, Vowden K, Vowden P. Integrating adjunctive therapy into practice: the importance of recognising 'hard-to-heal' wounds. World wide wounds 2006. Available from: <http://www.worldwidewounds.com/2006/december/Troxler/Integrating-Adjunctive-Therapy-Into-Practice.html> 2. Forssgren, et al. Leg Ulcer Point Prevalence can be Decreased by Broad-scale Intervention: a Follow-up Cross-sectional Study of a Defined Geographical Population. Acta Derm Venereol 2008; 88: 252–256. 3. Dowsett C et al. Reducing the burden of chronic wounds in the Community using single use NPWT. JCN supplement 2015; 29(5):1-20. 4. European Wound Management Association EWMA. Position Document: Hard-to-heal wounds: a holistic approach. London: MEP Ltd; 2008. 5. Dowsett C, et al. Use of PICO[®] to improve clinical and economic outcomes in hard-to-heal wounds. Wounds International. 2017;3, 53–58. 6. Moore Z, Dowsett C, Smith G, et al. TIME CDST: an updated tool to address the current challenges in wound care. J Wound Care. 2019;28(3):154-161.



Shaping

what's possible
in wound care

Smith & Nephew Pty Ltd
Australia
T +61 2 9857 3999
F +61 2 9857 3900
smith-nephew.com/en-au

Smith & Nephew Ltd
New Zealand
T +64 9 820 2840
F +64 9 820 2841
smith-nephew.com/en-nz

®Trademark of Smith+Nephew
All trademarks acknowledged
©July 2023 Smith+Nephew
30501-anz V2 07/23

This material is intended for healthcare professionals. For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use. Products listed and outline of care are examples only. Product selection and management should always be based on comprehensive clinical assessment. The results of in vitro simulation testing have not been proven to predict clinical performance.