








Fitting*

1 Select cannula size

Select appropriate cannula size; A clear gap must be visible around each prong. The sizing chart on the following page should only be used as a guide.

Ensure gap around prongs*



						
	XS	S	M	L	XL	XXL
	OJR410	OJR412	OJR414	OJR416	OJR418	OJR520

2 Prepare skin

Prepare the patient's skin according to hospital protocol.

3 Connect to gas source

Connect the cannula to the gas source.

TIP: Place hand close to nasal prongs to ensure that there is gas flow exiting the prongs.

4 Remove the first F&P Wigglepad™ 2 tabs

Remove the first backing tabs from the Wigglepads 2 and avoid touching the adhesive.

5 Insert cannula

Insert the cannula into the nares.

Ensure the cannula bridge rests close to the nose without touching the septum.

DO NOT stretch the cannula during application. Stick the Wigglepads 2 to the patient's cheeks.

6 Secure cannula

Remove the second backing tabs and stick the Wigglepads 2 onto the cheeks.



* IMPORTANT: Always refer to the user instructions supplied with the product for full set up instructions, warnings, cautions and contraindications.

Remove Cannula

- 1 Place fingertip**
Place fingertip on the outside edge of the Wigglepads 2 and gently peel the cannula away from the Wigglepads 2.



Wigglepads 2 removal and replacement

- 1 Lift the edge**
Lift the edge of the Wigglepads 2. Use a damp cloth to wipe the patient's skin and the underside of the Wigglepads 2 while gently peeling away from the patient's face.



- 2 Adhere**
Adhere the replacement Wigglepads 2 to the cannula, remove the first backing tabs and stick the Wigglepads 2 onto the patient's cheeks.



- 3 Remove**
Remove the second backing tabs and stick onto the patient's cheeks.

Sizing Information

NOTE:

Prongs must not create a seal in the nares. A clear gap must be visible around each prong. Patient weight should only be used as a guide.

Weight (Kg)	0.5	1	2	2.5	0.9	1	3.5	4	1	1.5	8	10
Correlated age*	23 wkGA	29 wkGA	33.5 wkGA	35 wkGA	27 wkGA	28 wkGA	40 wkGA	42.5 wkGA	28 wkGA	31 wkGA	6.6 mo	15.4 mo
Interface size	OJR410 (XS)				OJR412 (S)				OJR414 (M)			

Weight (Kg)	3	3.5	18	20	5	7	25	30	8	12.5	30	40
Correlated age*	37.5 wkGA	40 wkGA	4.9 yr	5.6 yr	47.5 wkGA	4.7 mo	7.6 yr	9.5 yr	6.4 mo	2.2 yr	9.5 yr	12 yr
Interface size	OJR416 (L)				OJR418 (XL)				OJR520 (XXL)			

wkGA = weeks of gestation; mo = months; yr = years
 * Weight data is based on F&P product validation studies.
 ** Age data is a correlation to weight data based on a combination of Fenton, WHO and CDC growth charts.