

The logo consists of a dark blue square with the word "SIGVARIS" written in white, uppercase, sans-serif font in the center.

INSTRUCTIONS AND CARE

- Wash daily by hand or in a laundry bag in a washing machine (delicate wash cycle) at max. 40 degrees. For hand-washing we recommend using our mild hand-washing detergent SIGVARIS WASHING SOLUTION. Do not use fabric softener!
- Please do not put the stockings in the dryer as this will cause them to overheat and generate static electricity. Press the compression stockings between two dry towels and lay them out to dry in a shady place with plenty of ventilation. Do not place on heaters/radiators.
- Do not iron.
- Do not dry clean. Do not use chemicals.

Thigh stocking with grip-top:

Human skin releases sweat, oil and skin cells. Thus, the silicone surface of the grip-top must be treated with particular care. Clean the silicone surfaces with SIGVARIS CLEAN PAD from time to time for better grip.

Thigh stocking with waist attachment:

Before washing, please match and fasten the Velcro strip carefully so that the entire area is protected.

Caution:

Sharp fingernails and rings with sharp edges can quickly result in defects. Pay due attention to these risks when donning and removing your stockings, and always wear rubber gloves!

Do not cut off any loose threads inside the stockings.