

HARTMANN



Continence Management Catalogue

2023



MoliCare[®]

MoliCare®

Protection. Every day.

HARTMANN has been providing healthcare products for over 200 years and incontinence products for over 50 years. We take a holistic approach to healthcare, placing emphasis on the tangible benefits for patients and those who care for them, aiming to enhance their well-being.

The MoliCare Difference

MoliCare is dedicated to meeting the most stringent standards of quality, dependability, and effectiveness by utilising our advanced technology to develop innovative solutions for managing incontinence. The aim is to address the needs of individuals living with incontinence, enhance the daily duties of skilled caregivers, and satisfy the expectations of healthcare professionals.

Incontinence can disrupt the skin's natural barrier and acid protection mantle, leaving it more susceptible to irritation and damage from urine and faeces.¹ In order to prevent Incontinence Associated Dermatitis (IAD), which can cause skin inflammation, blisters, itchiness or skin erosion²⁻⁴, it is important to take measures to protect the skin from the harmful effects of incontinence.



MoliCare provides absorbent products that are specifically designed to keep the skin dry and maintain its acid protection layer. This can help alleviate the negative consequences of stress factors associated with incontinence and prevent the development of IAD. This is in line with expert recommendations, which emphasises the importance of preventing IAD as a means of improving quality of life for individuals living with incontinence.⁵

pH of Skin on Incontinence Briefs⁶



MoliCare Absorbent Products:


- ✓ Significantly lower skin pH values compared to standard briefs, when skin is exposed to alkaline solution.^{6,7}
- ✓ Lower the prevalence of erythema^{6,7} which is skin redness, an early sign of skin irritation.

Graph data from Bliss Study Absorbent Briefs Containing Curly Fiber Lower (Acidify) Skin pH Reducing Risk for Incontinence Associated Skin Damage in Older Nursing Home Residents.⁶ ANOVA F=8.83 for thigh and F=10.12 for forearm pH differences, p < .001 for both; in post-hoc comparisons, skin pH values after wet brief using MoliCare are lower than all other conditions (p=.001 for all).⁶




✓ MoliCare controls odour and promotes a healthy skin environment.

- Buffers urine to a skin-friendly pH
- Maintains a skin-healthy pH balance
- Preserves the acid protection mantle
- Facilitates rapid fluid distribution within the brief
- Keeps the skin feeling dry



View MoliCare's optimal skin compatibility video to understand how it works.

Scan the code or visit:
lead.me/ph-balance



1. Shigeta, Yoshie, et al. "Exploring the relationship between skin property and absorbent pad environment." Journal of Clinical Nursing 18.11 (2009): 1607-1616. 2. Gray, Mikel. "Optimal management of incontinence-associated dermatitis in the elderly." American journal of clinical dermatology 11.3 (2010): 201-210. 3. Gray, Mikel, et al. "Incontinence-associated dermatitis: a consensus." Journal of Wound Ostomy & Continence Nursing 34.1 (2007): 45-54. 4. Junkin, Joan, and Joan Lerner Selekof. "Beyond "diaper rash": Incontinence-associated dermatitis Does it have you seeing RED?" Nursing2020 38.11 (2008): 56h1-56h10. Warner RR, Stone KJ, Boissy YL. Hydration disrupts human stratum corneum ultrastructure. J Invest Dermatol. 2003;120(2):275-84. 5. Gray, Mikel, et al. "Incontinence-associated dermatitis: a comprehensive review and update." Journal of Wound Ostomy & Continence Nursing 39.1 (2012): 61-74. 6. Bliss et al., An economic evaluation of four skin damage prevention regimens in nursing home residents with incontinence; economics of skin damage prevention; J Wound Ostomy Continence Nurse. 2007. 7. Data on file. The protective effect of the new diaper design on the epidermal barrier. 2022.



The MoliCare Skintegrity System

MoliCare offers a comprehensive approach to incontinence management by combining MoliCare absorbent products with MoliCare Skin care products. This combination, known as the MoliCare Skintegrity System, helps to maintain a healthy skin pH balance and help prevent the occurrence of Incontinence Associated Dermatitis (IAD).

MoliCare absorbent products work in conjunction with MoliCare Skin care products to maintain the acid protection mantle and provide a skin-neutral pH balance.¹

Our goal is to provide innovative solutions that improve the management of incontinence, allowing individuals to live with greater comfort and dignity.



Health Economics of MoliCare

Incidence

In an acute care study of 5,432 patients¹:

- 45.7% of patients had a prevalence of IAD
- 73% of IAD patients developed IAD following admission

IAD accompanied by unrelieved pressure results in a 5-fold increase in pressure injury risk²


Cost - IAD¹

The average cost to treat one mild to moderate IAD episode in one resident is \$220.

IAD occurrence can be up to 52%, resulting in costs of \$3,294 per month.

(This example assumes 20% in a 100 bed long term care facility with 75% incontinence at April 2018 rate¹)


Cost - Pressure Injuries^{2,3}



¹ Gray, M., Giuliano, K.K., JWOCN, Incontinence Associated Dermatitis, Characteristics and Relationship to Pressure Injury A Multisite Epidemiologic Analysis, January/February 2018. ² Melter, Pressure Ulcer Assessment and Treatment; Wild Iris Medical Education 2014; 17. ³ Bliss et al., An economic evaluation of four skin damage prevention regimens in nursing home residents with incontinence; economics of skin damage prevention; J Wound Ostomy Continence Nurse. 2007;34(2): 143-152.

¹ Except MoliCare® Premium Men Pad 3 Drops, MoliCare® Skin Body Lotion, MoliCare® Hand Cream, MoliCare® Protection Foam

3 easy steps to select the right pad



STYLE

Various styles suit different levels of mobility.

Active


Most styles suitable. Selection based on personal preference.



MoliCare Premium Pads and Pants

Active with restricted mobility and/or cognitive difficulties (e.g. Dementia)

Pull-up recommended. Looks and feels like normal underwear.



MoliCare Premium Mobile

Restricted/supported mobility or bedbound

All-in-one product recommended. Easiest option for the carer.



MoliCare Premium Elastic, Premium Form and Fixpants

ABSORBENCY

Giving the right level of protection – Choose an absorbency that suits your level of incontinence and keeps you dry.



Less Absorbent




More Absorbent

Look for the drop rating on pack
More drops = Higher absorbency

SIZE

To ensure maximum comfort and leakage protection:

- Measure waist at the largest width between the waist and hips.
- Select the product size based on the size guide printed on the product packaging.



Leakage may occur if the product is too large or too small.

MoliCare Advanced Features



Active Skin Protection

Skin-friendly, pH balance ensures maximum skin compatibility. Keeps skin dry with a 3-layer absorbent core and advanced fast dry effect.



Odour neutraliser

Neutralises unpleasant odours for enhanced well-being.



Dermatologically tested

Helps to maintain healthy skin.



Adjusted fit

Safety & comfort through optimal fit; elasticised side panels ensure better adaptation to body shape and movements, delivering improved comfort.



Air active

Skin-friendly due to high quality air active material with an extremely low potential to trigger allergic reactions.



All-round protection

Provided by fast absorption, anti-leakage system and odour neutraliser.



Anatomically shaped fit

MoliCare Premium MEN Pants: Feels like regular underwear. Front shifted absorbent core for extra protection where men need it most.



Anti-leakage system

For an improved feeling of safety.



Hygienic care

Single-use secure protection.



Prefolded & compressed

For easy use and storage.



Re-closable tapes

2 re-closable tapes (one for each side) for easy and secure closing.



Secure fixation

Designed to be worn with regular underwear.



Secure fixation

MoliCare Premium Fixpants: Ensure a safe and hygienic fixation of anatomically shaped large incontinence pads.



Soft & skin-friendly

Soft, breathable material for healthy skin.



Underwear like

Ideal body fit with cross-elasticated threads and extra long leg joints to ensure high wearing comfort.



Washable

Washable up to the number of times indicated on the label in a laundry net and can be tumble dried at low heat.



Wetness indicator on the back sheet

Shows when it is time to change.



3-layer absorbent core

Improved absorption core technology that helps to prevent skin irritations.



Comfortable fit

The soft waistband, the skin-friendly, air active material and the cross-elasticated threads ensure high wearing comfort.



Discreet underwear feeling

Discreet underwear-free through anatomical shape and soft, stretchy material for a comfortable fit.



Elasticated side panels

Elasticised side panels in the hip area ensure better adaptation to body shape and movements, delivering improved comfort for the patients.



Ergonomic use

Approximately 20% faster application as well as easier and more ergonomic handling for caregivers.¹



Extra protection

MoliCare Premium Bed Mat: Absorbs liquid and provides extra protection for various surfaces.



Advanced skin protection

Thanks to a balanced pH.



Anti-leakage system

For an improved feeling of safety.



Dermatologically tested

Helps to maintain healthy skin.



Extra odour neutralisation

Superabsorbent polymers neutralise unpleasant odours



Active skin protection

Up to 64% faster urine absorption compared to the previous MoliCare® Premium Form.






¹ W. Schneider, H. Mutschler, R. Haller, C. le Mathe: Test for ergonomic evaluation of incontinence products, Oct 5-8, 2015 Dijon. The reduction from 4 to 2 tapes (one for each side) combined with no side panels in the front simplifies the application by reducing the steps from 8 to 4 compared to a traditional slip with 4 side panels and 4 tapes. Additionally, due to the elasticated side panels, there is less need of readjusting the tapes.

PADS

1-5 

MoliCare® Premium Lady Pad



1 Drop	Code	Fit	ISO Max Absorbency in mL*	Units Per Pack	Pack Per Carton
	168132	Women	181	14	12
2 Drops					
	168634	Women	331	14	18
3 Drops					
	168644	Women	504	14	12
4 Drops					
	168682	Women	782	14	12
5 Drops					
	168670	Women	1029	14	12

* Measurement in accordance with ISO 11948/1 (Rothwell method)



MoliCare® Premium Lady Pad 1 Drop



MoliCare® Premium Lady Pad 2, 3, 4, 5 Drops

PADS



MoliCare® Premium Men Pad



ANATOMICAL FIT
all-round
protection

2 Drops	Code	Fit	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	168600	Men	330	14	12
3 Drops					
	168603	Men	406	14	8
4 Drops					
	168705	Men	546	14	12
5 Drops					
	168801	Men	852	14	12

Except MoliCare® Premium Men Pad 3 drops

MoliCare® Premium Men Pad 2 Drops




MoliCare® Premium Men Pad 3, 4, 5 Drops

PADS

2 - 4 

MoliCare® Pad




2 Drops	Code	Fit	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	168100	Unisex	321	28	12
3 Drops					
	168101	Unisex	481	28	6
4 Drops					
	168102	Unisex	865	28	6

SUPPORT PADS

3.5 

MoliCare® Rectangular



3.5 Drops	Code	Dimensions	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	148530	11 x 35 cm	691	56	4


PULL-UPS

5 7 M-L 

MoliCare® Premium Lady Pants




5 Drops	Code	Size	Waist Measurement	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	915867	M	80 - 120 cm	843	8	8
	915868	L	100 - 150 cm	844	7	8


7 Drops	Code	Size	Waist Measurement	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	915877	M	80 - 120 cm	1442	8	4
	915870	L	100 - 150 cm	1373	7	4

 Available in Outpatient markets

MoliCare® Premium Men Pants




5 Drops	Code	Size	Waist Measurement	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	915817	M	80 - 120 cm	974	8	4
	915818	L	100 - 150 cm	1001	7	4

7 Drops	Code	Size	Waist Measurement	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	915827	M	80 - 120 cm	1455	8	4
	915828	L	100 - 150 cm	1412	7	4

 Available in Outpatient markets

PULL-UPS





5 - 10 XS - XL 

MoliCare® Premium Mobile



DISCREET & COMFORTABLE



5 Drops	Code	Size	Waist Measurement	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	915861	S	60 - 90 cm	941	14	4
	915852	M	80 - 120 cm	1017	14	3
	915853	L	100 - 150 cm	1198	14	4
	915854	XL	130 - 170 cm	1183	14	4
6 Drops						
	915840	XS	45 - 70 cm	1361	14	4
	915831	S	60 - 90 cm	1475	14	4
	915832	M	80 - 120 cm	1662	14	3
	915833	L	100 - 150 cm	1963	14	4
	915834	XL	130 - 170 cm	2140	14	4
8 Drops						
	915871	S	60 - 90 cm	1791	14	4
	915872	M	80 - 120 cm	2015	14	3
	915873	L	100 - 150 cm	2279	14	4
	915874	XL	130 - 170 cm	2421	14	4
10 Drops						
	915878	M	80 - 120 cm	2626	14	3
	915879	L	100 - 150 cm	2613	14	4
	915880	XL	130 - 170 cm	2757	14	4

LARGE PADS

4 – 9

MoliCare® Premium Form



Recommended with MoliCare® Premium Fixpants

4 Drops	Code	Fit	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	168404	Unisex	1493	32	4
5 Drops					
	168405	Unisex	1662	32	4
6 Drops					
	168406	Unisex	2353	32	4
6 Drops					
	168401	Men	2441	32	4
7 Drops					
	168407	Unisex	2598	32	4
8 Drops					
	168408	Unisex	3017	32	4
9 Drops					
	168409	Unisex	3497	16	4

LARGE PADS

8 +SIZE


MoliCare® Premium Form + SIZE



Recommended with MoliCare® Premium Fixpants

8 Drops	Code	Fit	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	168519	Unisex	3347	18	3

MoliCare® Premium Form STOOL




Recommended with MoliCare® Premium Fixpants

STOOL	Code	Fit	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	168419	Unisex	1295	32	4

FIXPANTS

S - 5XL

MoliCare® Premium Fixpants Long Leg



	Code	Size	Waist Measurement	Colour	Units Per Pack	Pack Per Carton
	947790	S	45 - 80 cm	Yellow	25	1
	947791	M	60 - 100 cm	Blue	25	1
	947792	L	80 - 120 cm	Brown	25	1
	947793	XL	100 - 160 cm	Green	25	1
	947794	XXL	140 - 180 cm	Grey	25	1
	947789	XXXL	160 - 200 cm	Orange	25	1
	947825	5XL	200 - 240 cm	Purple	25	1

MoliCare® Premium Fixpants Short Leg



	Code	Size	Waist Measurement	Colour	Units Per Pack	Pack Per Carton
	947711	S	45 - 80 cm	Yellow	25	1
	947712	M	60 - 100 cm	Blue	25	1
	947713	L	80 - 120 cm	Brown	25	1
	947714	XL	100 - 160 cm	Green	25	1
	947715	XXL	140 - 180 cm	Grey	25	1

ALL-IN-ONE SLIPS

6 - 10 S - XL

MoliCare® Premium Elastic



COMFORTABLE & EASY USE



Drops	Code	Size	Waist Measurement	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
6	165271	S	70 - 90 cm	1667	30	3
	165272	M	85 - 120 cm	2159	30	3
	165273	L	115 - 145 cm	2349	30	3
	165274	XL	140 - 175 cm	2786	14	4
7	165371	S	70 - 90 cm	1998	30	3
	165372	M	85 - 120 cm	2414	30	3
	165373	L	115 - 145 cm	2719	30	3
	165374	XL	140 - 175 cm	3180	14	4
8	165471	S	70 - 90 cm	2356	26	3
	165472	M	85 - 120 cm	3144	26	3
	165473	L	115 - 145 cm	3220	24	3
	165474	XL	140 - 175 cm	3591	14	4
9	165571	S	70 - 90 cm	2546	26	3
	165572	M	85 - 120 cm	3521	26	3
	165573	L	115 - 145 cm	3856	24	3
	165574	XL	140 - 175 cm	3951	14	4
10	165671	S	70 - 90 cm	2771	22	3
	165672	M	85 - 120 cm	3696	14	4
	165673	L	115 - 145 cm	4279	14	4
	165674	XL	140 - 175 cm	4380	14	4

ALL-IN-ONE SLIPS


6 9 XS-XL

MoliCare® Slip XS




6 Drops	Code	Size	Waist measurement	ISO max absorbency in mL	Units per pack	Pack per carton
	169190	XS	40-70 cm	1367	30	4

MoliCare® Slip



9 Drops	Code	Size	Waist Measurement	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	165531	S	60 - 90 cm	2153	14	4
	165532	M	90 - 120 cm	3196	14	4
	165533	L	120 - 150 cm	3815	14	4

Dignity® Bariatric




Code	Size	Waist Measurement	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
916680	Bariatric	160 - 239 cm	2445	8	4



BED MATS



MoliCare® Premium Bed Mats



5 Drops




Code	Dimensions	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
161061	40 x 60 cm	449	30	6
161065	60 x 90 cm	971	30	4

7 Drops



161070	60 x 90 cm	1241	25	4
--------	------------	------	----	---

 40 x 60 cm variant not compressed



Order Samples

Scan the QR code or visit:
<http://l.ead.me/molicare-samples>



MoliCare Premium Lady Pad

Absorbency	Code	Size	Units Per Sample Pack
1 Drop	168430	-	2
2 Drops	168433	-	2
3 Drops	168434	-	2
4 Drops	168437	-	2
5 Drops	168438	-	2



MoliCare Premium Lady Pants

Absorbency	Code	Size	Unit Per Sample Pack
5 Drops	915810	M	1
	915811	L	1
7 Drops	915815	M	1
	915816	L	1

7 Drops Available in outpatient markets.



MoliCare Premium Form

Absorbency	Code	Size	Units Per Sample Pack
4 Drops	168204	-	2
5 Drops	168205	-	2
6 Drops	168206	-	2
6 Drops for Men	168201	-	2
7 Drops	168207	-	2
8 Drops	168208	-	2
9 Drops	168209	-	2
+ SIZE	168974	-	2
STOOL	168874	-	2



MoliCare Premium Elastic

Absorbency	Code	Size	Units Per Sample Pack
6 Drops	165291	S	2
	165292	M	2
	165293	L	2
7 Drops	165294	XL	2
	165391	S	2
	165392	M	2
	165393	L	2
8 Drops	165394	XL	2
	165491	S	2
	165492	M	2
9 Drops	165493	L	2
	165494	XL	2
	165591	S	2
	165592	M	2
10 Drops	165593	L	2
	165594	XL	2
	165691	S	2
	165692	M	2
	165693	L	2
	165694	XL	2



MoliCare Premium Men Pad

Absorbency	Code	Size	Units Per Sample Pack
2 Drops	168503	-	2
3 Drops	168604	-	2
4 Drops	168501	-	2
5 Drops	168609	-	2



MoliCare Premium Men Pants

Absorbency	Code	Size	Unit Per Sample Pack
5 Drops	915802	M	1
	915803	L	1
7 Drops	915807	M	1
	915808	L	1

7 Drops Available in outpatient markets.



MoliCare Premium Fixpant Long Leg

Size	Code	Unit Per Sample Pack
Small	947800	1
Medium	947801	1
Large	947802	1
XL	947803	1
XXL	947804	1



MoliCare Premium Mobile

Absorbency	Code	Size	Units Per Sample Pack
5 Drops	915615	S	2
	915616	M	2
	915617	L	2
	915648	XL	2
6 Drops	915618	XS	2
	915619	S	2
	915620	M	2
	915621	L	2
	915622	XL	2
8 Drops	915623	S	2
	915624	M	2
	915625	L	2
	915656	XL	2
10 Drops	915825	M	2
	915826	L	2
	915829	XL	2



MoliCare Pad

Absorbency	Code	Size	Units Per Sample Pack
2 Drops	168728	-	2
3 Drops	168828	-	2
4 Drops	168928	-	2



MoliCare Slip XS

Absorbency	Code	Size	Units Per Sample Pack
6 Drops	169160	XS	2



MoliCare Slip

Absorbency	Code	Size	Units Per Sample Pack
9 Drops	165751	S	2
	165752	M	2
	165753	L	2



Helps. Cares. Protects.

For further information please contact
your HARTMANN representative,
or call HARTMANN customer care
toll free 1800 805 839

www.molicare.au

PAUL HARTMANN Pty Ltd. ABN 35 000 099 589.
Level 5, 1 Thomas Holt Drive, Macquarie Park NSW 2113
PO Box 943, North Ryde BC NSW 1670
© Registered trademark PAUL HARTMANN AG.

9166444_042023