

danuButton® Gastrostomy Button



Safety meets Comfort:

- German Manufacturer ✓
- Reliable Supply Chain ✓
- 84 References ✓
- Accessories in the Set ✓


- Enhanced Comfort ✓
- Long-lasting ✓
- Secure Hold ✓
- Apple-shaped ✓
- Easy Cleaning ✓

• Latex-free • PVC-free • DEHP-free • BPA-free

danuButton® Gastrostomy Button


The danuButton® is a low-profile feeding tube made of high-quality medical silicone.

- Small, low in profile and discreet for more mobility and quality of life
- Large selection of sizes and lengths for individual and perfect fit
- Soft, flexible silicone for high patient comfort
- Better patient safety due to high-quality materials, free of harmful substances
- Apple-shaped inner retention balloon with high balloon stability for long usage, secure fixation and to avoid leakage




The danuButton® is supplied in three different sets:

The **danuButton® Premium Set** and the **danuButton® Compact Set** include a comprehensive range of ENFit® accessories for convenient and easy replacement of the gastrostomy button. The **danuButton® Compact XS** is a set with reduced accessories.



Recommended balloon filling volume	
Shaft size in Fr	Balloon volume in ml
12	3.0
14	4.0
16	5.0
18	5.0
20	8.0
24	8.0



The danumed® Extension Set


The Extension Set is the connecting piece between the danuButton® and usually an enteral syringe for the administration of enteral nutrition, liquids or medication.

The Extension Sets included in the danuButton® Premium Set are made of kink-resistant, biocompatible polyurethane. They are equipped with ENFit® connectors according to DIN EN ISO 80369-3 and danumed's innovative rotating ring, the ENSwivel®.



Set content of the danuButton® Premium Set:

1. danuButton® Gastrostomy Button
2. danumed® Extension Set for Gastrostomy Button, 30 cm, Bolus-ENFit®, ENSwivel®, straight button connector
3. danumed® Extension Set for Gastrostomy Button, 30 cm, Y-ENFit®, ENSwivel®, angled button connector
4. danumed® Introducer
5. danumed® Balloon Syringe 10 ml
6. danumed® Enteral Syringe, 60 ml, ENFit®
7. optiLube Lubricant
8. 4 Gauze pads



danuButton® Premium Set Product Codes:

Shaft size 12 French	
Length in mm	REF
8	DA 1041 1208
10	DA 1041 1210
12	DA 1041 1212
15	DA 1041 1215
17	DA 1041 1217
20	DA 1041 1220
23	DA 1041 1223
25	DA 1041 1225
27	DA 1041 1227
30	DA 1041 1230
35	DA 1041 1235
40	DA 1041 1240
45	DA 1041 1245
50	DA 1041 1250

Shaft size 14 French	
Length in mm	REF
8	DA 1041 1408
10	DA 1041 1410
12	DA 1041 1412
15	DA 1041 1415
17	DA 1041 1417
20	DA 1041 1420
23	DA 1041 1423
25	DA 1041 1425
27	DA 1041 1427
30	DA 1041 1430
35	DA 1041 1435
40	DA 1041 1440
45	DA 1041 1445
50	DA 1041 1450

Shaft size 16 French	
Length in mm	REF
8	DA 1041 1608
10	DA 1041 1610
12	DA 1041 1612
15	DA 1041 1615
17	DA 1041 1617
20	DA 1041 1620
23	DA 1041 1623
25	DA 1041 1625
27	DA 1041 1627
30	DA 1041 1630
35	DA 1041 1635
40	DA 1041 1640
45	DA 1041 1645
50	DA 1041 1650

Shaft size 18 French	
Length in mm	REF
8	DA 1041 1808
10	DA 1041 1810
12	DA 1041 1812
15	DA 1041 1815
17	DA 1041 1817
20	DA 1041 1820
23	DA 1041 1823
25	DA 1041 1825
27	DA 1041 1827
30	DA 1041 1830
35	DA 1041 1835
40	DA 1041 1840
45	DA 1041 1845
50	DA 1041 1850

Shaft size 20 French	
Length in mm	REF
8	DA 1041 2008
10	DA 1041 2010
12	DA 1041 2012
15	DA 1041 2015
17	DA 1041 2017
20	DA 1041 2020
23	DA 1041 2023
25	DA 1041 2025
27	DA 1041 2027
30	DA 1041 2030
35	DA 1041 2035
40	DA 1041 2040
45	DA 1041 2045
50	DA 1041 2050

Shaft size 24 French	
Length in mm	REF
8	DA 1041 2408
10	DA 1041 2410
12	DA 1041 2412
15	DA 1041 2415
17	DA 1041 2417
20	DA 1041 2420
23	DA 1041 2423
25	DA 1041 2425
27	DA 1041 2427
30	DA 1041 2430
35	DA 1041 2435
40	DA 1041 2440
45	DA 1041 2445
50	DA 1041 2450



danumed Medizintechnik GmbH

Individual solutions are our standard.

We are a family-run, innovative, German medical technology company with high quality standards, many years of expertise in enteral nutrition and a clear customer focus.

Individual customer requests are met with flexibility and reliability.
Personal, competent, cooperative.

For more information, please contact us on **(02) 9882 3666** or at **amsl@amsl.com.au**

amsl.com.au

ALWAYS READ THE LABEL AND FOLLOW THE DIRECTIONS FOR USE.
Manufactured by danumed. Distributed in Australia by AMSL.
AMSL is a Dexcom company. PR-100-763 July 2023

