

Be Independent.™



# ASSIST-A-TRAY™

BECAUSE YOUR TRAY CAN SWIVEL AND HELP YOU STAND



PRODUCT #2050



WATCH A DEMO





Easy Stand Safety Handle



Accommodates Any Lift Chair



Extends in Height and Depth to Fit Any Couch or Chair



Arm Pivots and Tray Swivels

## FEATURES & BENEFITS

- **Ergonomic Safety Handle** - Makes standing and sitting easy and safe.
- **Swivel Tray** - Multi-use laptop/TV tray with cup holder and utensil compartment. Tray swivels 360°.
- **Adaptable** - Adjusts in length and height to fit either side of your favorite couch or chair.
- **Super Grip Rubber Pads** - Ensure stability and protect floor from damage.

## SPECIFICATIONS

Height of Tray: 26" to 32"  
 Height Adjustment Handle: 34" to 40"  
 Dimensions of Handle: 6" x 6" 1.25 Diameter  
 Base Length Adjustment: 20" to 36"  
 Weight: 25 lbs  
 Size of Tray: 16" x 18"  
 Base Feet Measurements: 4" x 7"  
 Weight Limit of Tray: 30 lbs  
 Weight Limit of Handle: 250 lbs  
 Case of 2 Dimensions: 29" L x 19" W x 8" H; 50 lbs  
 Package Dimensions: 29" L x 19" W x 4" H; 25 lbs

## EASY ASSEMBLY (No Additional Tools Required)



Slide Two Base Pieces Together



Insert Handle Into Base



Attach Tray to Handle



Adjust Depth of Base and Handle to Desired Height



Place Under Your Couch or Chair