


EZ ADJUST BED RAIL

Features:

- extends in length after installation; allowing you to have a "half" bed rail to help you get out of bed, and full-sized bed rail at night for maximum fall protection.
- secures to any bed with included safety strap
- accommodates either side of bed
- installs in minutes with 4 bolts and an Allen key (included)
- Fits mattress thickness of 30 - 45 cm



BEDSIDE MIGHTY RAIL	METRIC		IMPERIAL	
Bed Type	Attaches to any home or hospital bed with included safety strap			
Rail height from base of mattress	58.42	cm	23	inch
Width of rail	Locks into place at 66, 86 & 106 after installed on bed			
Base Dimensions	46.99 x 63.5	cm	18.5 x 25	inch
Gap between Rails	7.62	cm	3	inch
Tube Diameter	2.54	cm	1	inch
Product weight	6.35	kg	14	lbs
Made of	Powder coated steel			
Weight Capacity	136.08	kg	300	lbs