

Be Independent™

# PT BEDCANE

BECAUSE EVERYONE CAN USE A HELPING HAND TO GET OUT OF BED

NEW



INCLUDES  
ORGANIZER POUCH

PRODUCT #2100





Ergonomic Handle Allows For 2 Handed Grip



Includes Organizer Pouch



Secures To Bed Frame with Included Safety Strap

### FEATURES & BENEFITS

- **Ergonomic Handle** – Allows for 2 handed grip
- **A Pocket Organizer** – Provides storage space for handy items.
- **Height Adjustment** – Fits 10"-17" mattress
- **Easy Installation** – Installs in seconds with no additional tools required.

### SPECIFICATIONS

Height Adjustment From Mattress Base to Top of Handle: 19.25" to 25.25"  
 Handle Width: 6.5"  
 Depth of Horizontal Support Structure Under Mattress: 19"  
 Weight Capacity of Rail: 300 lbs  
 Weight of Product: 5.6 lbs  
 Made of: Steel frame with durable powder coat finish  
 Package Dimensions: 2.625" L x 8.25" W x 21.25" H; 6 lbs  
 Case of 4 Dimensions: 11" L x 8.75" W x 21.75" H; 24.5 lbs

### WHAT IS IN THE BOX



NO ADDITIONAL TOOLS REQUIRED



WEIGHT CAPACITY