



## BEDSIDE ECONORAIL

### Features:

- bedside handle for easy standing and sitting
- secures to any bed with included safety strap
- installs in seconds with no tools required
- includes convenient organiser pouch
- easy four step installation



BEDSIDE MIGHTY RAIL	METRIC		Imperial	
	Rail height from base of mattress	54.61	cm	21.5
Width of rail from base to top	40.64 - 30.48	cm	16 - 12	inch
Base Dimensions	50.8 x 40.64	cm	20 x 16	inch
Tube Diameter	2.54	cm	1	inch
Product Dimensions	54.61 x 50.8 x 40.64	cm	21.5 x 20 x 16	inch
Made of	Anodised Aluminium			
Weight Capacity	136.08	kg	300	lbs