

IV Pole Nesting



4Healthcare



4Healthcare

IV Nesting Pole

PRODUCT CODE: 4H164

25mm stainless steel upright can accept 1x syringe pump

2 or 4 hook options


Adjustable pole with pneumatic damping mechanism to retard adjustment speed

Dimensions
Base: 520mm diameter
Height adjustment: 1400mm - 2250mm

Nestable soft polyurethane easy clean base

Diameter 50mm swivel rubber castors

Option



Push Ring
4H160PH



6 Welsh Place Landsdale WA 6065
Phone: (08) 9302 4757 Fax: (08) 9302 4524
Website: www.4healthcare.com.au
Email: sales@4healthcare.com.au

*Due to 4Healthcare's commitment to continuous improvement, our products are subject to change without notice.