



**ABLE LIFE**  
Comfortably Independent



**SURGICAL HOUSE**  
Medical & Surgical Supply Specialists

46 KING EDWARD ROAD, OSBORNE PARK 6017 TEL. +61 8 9381 4199  
PO BOX 1537, OSBORNE PARK DC WA 6916 FAX. +61 8 9382 3009  
surgical@surgicalhouse.com.au www.surgicalhouse.com.au

## FOLD-N-GO WALKER



### Features:

- Portable - quickly & easily folds up for storage in small spaces
- 6 inch swivel wheels with locking pin and rear Easy-Glide feet - allows for trouble-free maneuvering
- height adjustable - and supports up to 140kg
- contemporary look - three stylish colours to choose from



Black Walnut



Cobalt Blue



Regal Rose



FOLD-N-GO WALKER	Metric		Imperial	
Collapsed Dimensions	91.44 x 13.3 x 20.32	cm	36 x 5.25 x 8	inch
Wheels	15.24	cm	6	inch
Height Adjustment	11.43 x 12.7	cm	4.5 x 5	inch
User Height	147.32 to 203.2	cm	4'10" to 6'8"	
Open Depth	54.61		21.5	inch
Open Width	59.69		23.5	inch
Open Handle Width	45.72		18	inch
Weight	3.2	kg	7	lbs
Weight Capacity	140	kg	287	lbs

