

816 BREEZY BASIX TRANSIT WHEELCHAIR



alphacare



- Swingaway removable legrests and swingback removable armrests
- Push handles, seat depth, seat height and footplate angle are adjustable to suit individual client requirements
- Breathable padded nylon upholstery
- Tension adjustable backrest
- Seat to back connecting flap
- Heel straps
- Transport tie down points
- Anti tip bars and rollers help prevent tipping
- Lightweight aluminium folding frame

Optional:

- Drum brakes
- Folding back rest

Models Available:

- 81635 35cm seat width
- 81640 40cm seat width
- 81646 46cm seat width
- 81651 51cm seat width

Seat Width:	35/40/46/51cm
Seat Depth:	40-46cm
Seat Height:	48cm
Back Height:	41-43.5cm
Arm Height:	20.5-23cm
Overall Width:	54/59/65/70cm
Gross Weight:	14.9kg
Net Weight:	8.4kg
Max User Weight:	125 kg
Max User Weight:	100kg

