


Y-Suction Catheters Open Tip (Control Vent)

Designed for easy passage through endotracheal and tracheostomy tubing our Y-suction catheters offer an open distal end with an elliptical eye for easy and effective suctioning, incremental measurement markings commencing at 4cm from distal tip.

- ⊗ Polished open tip with bilateral side eyelets
- ⊗ Stepped connector
- ⊗ Soft clear tubing
- ⊗ Colour coded to support product selection
- ⊗ Available in a wide range of sizes
- ⊗ Inner sleeve with perforations to assist with maintaining aseptic technique
- ⊗ Incremental measurement markings



Y-Suction Catheters				STERILE
RE-ORDER	DESCRIPTION	GAUGE	COLOUR	CARTON QTY
AN042000	Suction Catheter with 380mm Y Type Control Vent	5Fr	Grey	50
AN042001	Suction Catheter with 380mm Y Type Control Vent	6Fr	Light Green	50

Table continued on next page

Y-Suction Catheters

Y-Suction Catheters				STERILE
RE-ORDER	DESCRIPTION	GAUGE	COLOUR	CARTON QTY
AN042002	Suction Catheter with 560mm Y Type Control Vent	8Fr	Light Blue	50
AN042003	Suction Catheter with 560mm Y Type Control Vent	10Fr	Black	50
AN042004	Suction Catheter with 560mm Y Type Control Vent	12Fr	White	50
AN042005	Suction Catheter with 560mm Y Type Control Vent	14Fr	Green	50
AN042006	Suction Catheter with 560mm Y Type Control Vent	16Fr	Orange	50
AN042007	Suction Catheter with 560mm Y Type Control Vent	18Fr	Red	50