

Air Lift Chair

COBALT
HEALTH



Proven pressure relief suitable for all-day use, while offering the mobility support of powered lift.

Product Specs

Air Lift 125

Safe working load: 125kg
Overall width: 800mm
Overall height: 1100mm
Overall depth: 1000mm
Seat width: 600mm
Seat depth: 500mm
Seat height: 500mm
Single motor seat recline: 45-105 °
Dual motor seat recline: 10-105 °
Foot rest angle: 5-90 °
Unit weight: 63kg

Order codes:

Single Motor - C09-SA-PU-BLU

Dual Motor - C09-DA-PU-BLU

Features

- Air cushion system for pressure relief
- Powered adjustment with remote control
- Adjustable backrest and leg rest
- Chair raises upwards to assist user in standing
- Hidden castors for portability
- Four-way stretch PU cover - spill resistant and fire retardant
- Ergonomic push handles
- Dual motor model offers near-to-flat recline

Air Lift 180

Safe working load: 180kg
Overall width: 800mm
Overall height: 1100mm
Overall depth: 1000mm
Seat width: 600mm
Seat depth: 500mm
Seat height: 500mm
Single motor seat recline: 45-105 °
Dual motor seat recline: 10-105 °
Foot rest angle: 5-90 °
Unit weight: 65kg

Order codes:

Single Motor - C09-SA-180-BLU

Dual Motor - C09-DA-180-BLU

Features

- Increased SWL of 180kg
- Built in battery back up for safety
- Air cushion system for pressure relief
- Powered adjustment with remote control
- Adjustable backrest and leg rest
- Chair raises upwards to assist user in standing
- Hidden castors for portability
- Four-way stretch PU cover - spill resistant and fire retardant
- Ergonomic push handles
- Dual motor model offers near-to-flat recline