

Delivering F&P Optiflow™ High Flow Therapy



Introducing F&P Airvo™ 2, delivering respiratory support and airway hydration.

The Airvo 2 is an integrated solution for delivering Optiflow high flow therapy. Fisher & Paykel Healthcare's leading humidification technology allows the Airvo 2 to comfortably provide high flows of humidified air (with supplemental oxygen if required) to spontaneously breathing patients, through the unique series of Optiflow interfaces.

Adjustable temperature and flow settings

- Three temperature settings, including 37°C (body temperature pressure saturated for a normothermic patient) i.e. optimal humidity
- Integrated flow generator delivers a wide flow range 10 - 60 L/min for adult patients and 2 - 25 L/min for pediatric patients
- No wall air supply required

Optiflow+ Nasal Interface

- F&P Evaqua™ technology, which reduces the formation of mobile condensate
- Soft, flexible prongs
- Prongs contoured to the patient's nose



Supplementary oxygen (when required)

- Oxygen can be added from a wall supply or cylinder
- Inbuilt ultrasonic oxygen analyzer requires no calibration, service or replacement

Designed for simple setup, use and cleaning


- Helpful onscreen animations assist with setup and troubleshooting

F&P AirSpiral™ heated breathing tube

- No separate temperature probes or heater-wire adapters required

An integrated system for delivering Optiflow high flow therapy

- The Airvo 2, with its wide range of accessories, provides versatility, mobility and convenience



Product codes (some accessories may not be available in all countries)		
	PT101xx	Airvo 2 humidifier
Breathing circuit	900PT561	Airvo 2 AirSpiral tube and chamber kit (10-pack)
	900PT562,563	Airvo 2 AirSpiral tube and chamber kit with nebulizer adapter (10-pack)
Adult interface	OPT1042, 1044, 1046	Optiflow 3S S/M/L nasal interface (20-pack)
	OPT942, 944, 946	Optiflow+ S/M/L interface (20-pack)
	OPT970	Optiflow+ Tracheostomy interface (20-pack)
	OPT971	Sputum guard spare part (20-pack)
	OPT980	Optiflow+ Mask adapter interface (20-pack)
Junior interface	OPT 316, 318	Optiflow Junior L/XL (20-pack)
	OJR416, 418	Optiflow Junior 2 Nasal interface L/XL (20-pack)
	OJR520	Optiflow Junior 2+ Nasal interface XXL (10-pack)
Hardware	900PT421	Mobbile pole stand
	900PT405	Pole mounting tray
	900PT426	Plastic basket
	900PT427	Oxygen bottle holder
	900PT411	UPS mounting kit

Find out more at www.fphcare.com/airvo

F&P, AirSpiral, Airvo, Evaqua, and Optiflow are trademarks of Fisher & Paykel Healthcare Limited. For patent information, see www.fphcare.com/ip

185048296 REV E © 2022 Fisher & Paykel Healthcare Limited

Fisher & Paykel
HEALTHCARE