

## New

### Apollo 4 Plus

### Dynamic Alternating Cushion


Product Code - APH055

12 Hours Lithium-ion battery powered  
Moulded / Contoured visco foam base / surround  
Maximum weight limit 24 stone / 152kg  
6 Alternating Cells  
Cycle - 1in2 - 6 minutes  
Ultra Light 1 kg pump unit




The Apollo 4 Plus Cushion is a high specification alternating cushion that can be used in both the hospital and community environments and is ideal for wheel chair users giving back freedom to go anywhere without compromising pressure relief and comfort.

Six alternating cells sit in a moulded visco foam base that not only ensure patient comfort and the best possible pressure relieving qualities, but also the design reduces the risk of the patient slipping forward.



Cycle 1



Cycle 2

Pressure Mapping Results Male 67 years 1.8m height 85kg weight