



OVER BED TABLE RANGE



axishealth

OVER BED TABLES

The Axis Health range of Over Bed Tables offers contemporary style combined with functionality.

A range of models are available to fit seamlessly with patient requirements. Colour options also available as a stylish solution to coordinate with the surrounding decor.



Vacuum formed top with lip surround



Laminated MDF timber top



Assisted lift height adjustment



Manual height adjustment



Split Table Top — Assisted Lift Height

- Sturdy mild steel construction powder coated off white
- Vacuum formed top with lip surround to retain liquids 700 x 425 and 195 x 400 (mm)
- Lightweight assisted height mechanism
- Locks in compression no release required for raising
- Modern ABS plastic base cover
- 55mm Diameter twin wheeled castors
- Height adjustment range 745 to 1050 (mm)
- Tilt section angles 24°, 31°, 37°, 42°, 47°
- Safe Working Load: 25kg evenly distributed

Order Code	Colour
AXST210-AW	Ash White (Standard)
AXST210-YB	Young Beech
AXST210-SB	Select Beech
AXST210-MW	Milano Walnut



Vacuum formed top with lip surround



Fixed Table Top — Assisted Lift Height

- Sturdy mild steel construction powder coated off white
- Vacuum formed top with lip surround to retain liquids 850 x 400 (mm)
- Lightweight assisted height mechanism
- Locks in compression no release required for raising
- Modern ABS plastic base cover
- 55mm Diameter twin wheeled castors
- Height adjustment range 745 to 1050 (mm)
- Safe Working Load: 25kg evenly distributed

Order Code	Colour
AXMT112TH-AW	Ash White (Standard)
AXMT112TH-YB	Young Beech
AXMT112TH-SB	Select Beech
AXMT112TH-MW	Milano Walnut



Vacuum formed top with lip surround

Fixed Table Top — Assisted Lift Height



- Sturdy mild steel construction powder coated off white
- Laminated MDF (15mm) timber top 850 x 400 (mm)
- Lightweight assisted height mechanism
- Locks in compression no release required for raising
- 55mm Diameter twin wheeled castors
- Height adjustment range 745 to 1050 (mm)
- Safe Working Load: 25kg evenly distributed

Order Code	Colour
AXKT110-AW	Ash White (Standard)
AXKT110-YB	Young Beech
AXKT110-SB	Select Beech
AXKT110-MW	Milano Walnut



Laminated MDF timber top

Fixed Table Top — Manual Height Adjustment



- Mild steel construction powder coated off white
- Laminated MDF (15mm) timber top 850 x 400 (mm)
- Locking hand wheel compression spring assisted raising
- 55mm Diameter twin wheeled castors
- Height adjustment range 745 to 1055 (mm)
- Safe Working Load: 25kg evenly distributed

Order Code	Colour
AXBT110-AW	Ash White (Standard)
AXBT110-YB	Young Beech
AXBT110-SB	Select Beech
AXBT110-MW	Milano Walnut



Laminated MDF timber top

Fixed Table Top — Manual Height Adjustment



- Mild steel construction powder coated off white
- Vacuum formed top with lip surround to retain liquids 850 x 400 (mm)
- Locking hand wheel compression spring assisted raising
- 55mm Diameter twin wheeled castors
- Height adjustment range 745 to 1055 (mm)
- Safe Working Load: 25kg evenly distributed

Order Code	Colour
AXBT112TH-AW	Ash White (Standard)
AXBT112TH-YB	Young Beech
AXBT112TH-SB	Select Beech
AXBT112TH-MW	Milano Walnut



Vacuum formed top with lip surround

Colours Available:



Ash White
(Standard)



Young Beech



Select Beech



Milano Walnut



405 – 409 Ferntree Gully Rd,
Mt Waverley, VIC, 3149, Australia

T +61 (3) 9922 0300
F +61 (3) 9544 5392

W www.axishealth.com.au