

Toilet Frames



Days Over Toilet Aid

The Days Over Toilet Aid assists transfers on and off the toilet. This toilet surround frame is height adjustable, with a push button adjustment. Powder coated steel frame with non slip rubber feet. Comes complete with a seat and splashguard.

Specifications

Seat Width	370mm
Width Between the Arms	450mm
Height Adjustable	440mm - 590mm
Maximum User Weight	125kg

Complies with AS/NZS 3973:2009

Code	Description	Qty/Box
JAN-AC20	Days Over Toilet Aid	4
JAN-AC20SP-008	Toilet Seat	
JAN-AC20SP-009	Bowl / Splashguard	
JAN-AC20SP-004	Rubber Tip	

Bariatric All In One

Redesigned with the carer in mind, the Bariatric All In One is an over toilet aid, shower chair and commode. Comes with square bowl with easy front access to remove bowl. The square bowl does not have a lid. Optional splash guard available. The seat lifts up and comes with a seat cover. The removable backrest adds stability when used as a bedside commode and can be removed when using over a toilet. The moulded armrests assist gripping for safe and easy transferring.

Dimensions

Width Between Arms	600mm
Height Adjustable	430 - 570mm
Safe Working Limit	295kg

Complies with AS/NZS 3973:2009

Code	Description	Qty/Box
MAX-AUSBAT70282	Bariatric All in One	2
Spare Parts		
MAX-70282-01	Toilet Seat with Hole	
MAX-70282-02	Toilet Seat Cover without Hole	
MAX-70282-03	Square Bowl (no lid)	
MAX-70282-05	Rubber Tips	

