



Perimeter Safety Cover

Perimeter Safety Cover is an easy-to-use coverlet that can be placed on top of any anti-decubitus mattress systems. Without the need to use side rails, Perimeter Safety Cover is able to safeguard patient's position in bed, and prevent patient fall and entrapment between the side rail and the mattress. It is an ideal coverlet that provides patient comforts, as well as securing patient safety at the same time.

PERIMETER Safety Cover

Product feature

- Smooth and slippery material will prevent friction, shear, and damage to patient skin.
- Simple installation that can attach the coverlet to bed frame with fastening straps on each sides.
- Fits most standard size mattresses.
- Easy to clean with removable side foam insert.
- Removable heel protection that helps positioning patient's feet.
- High quality water proof and high MVP transmission rate.
- Breathable top cover that can be use over Low Air Loss and Alternating pressure mattress systems.
- Easy patient transfer with mid-section openings.
- Effectively prevent patient fall without concerning about the dangers of using the side rails.



- 1 Easy patient transfer with mid-section openings
- 2 Removable heel protection that helps positioning patient's feet
- 3 Effectively prevent patient fall without concerning about the dangers of using the side rails
- 4 Buckle fasten to the bed frame to prevent the mattress slippage

Product Specification

Individual Model NO.

Mattress: SR607-GM11-11

Material

Coverlet: 100 % nylon with PU backing

Dimension

200L x 95W x 33H(Elastic foam height included) cm



SURGICAL HOUSE
Medical & Surgical Supply Specialists

46 KING EDWARD ROAD, OSBORNE PARK 6017

PO BOX 1537, OSBORNE PARK DC WA 6916

surgical@surgicalhouse.com.au

TEL. +61 8 9381 4199

FAX. +61 8 9382 3009

www.surgicalhouse.com.au