

# BVT02

Bladder Scanner



**KAI XIN**

# BVT02 Bladder Scanner

## Features

- ① Multi-language customized
- ② 3.5" color LCD screen, good interface style
- ③ Three patients figure mode (Standard / Obesity/ Children)
- ④ Scan time <4S
- ⑤ Machine beeps when the measurement exceeds the capacity of the alarm set value
- ⑥ Select different image contrast, can improve the measurement accuracy
- ⑦ Probe preset time automatically scan within specified time
- ⑧ Work time≥9 hours



## Application & Effects

### Home health

Control the micturition time of patients by measuring post-void residual volume , avoiding the irrevocable damage of bladder function caused by long-term open catheterization

### Geriatrics

Test urine volume can evaluate urinary retention , uracratia and other urinary system diseases

### Rehabilitation care

Prevent urinary tract infection caused by unnecessary catheterization and reduces patient pain

### Urology

Evaluate the therapeutic effect of certain drugs for urinary system disease by continuously measuring post-void residual volume



## Standard configuration

- Main unit 1pc
- 2.5MHz 3D probe 1pc
- Power adapter 1pc
- Internal batter 1pc

## Optional configuration

- calibrmtion tool



Ph: 1300 099 663 Fax: 1300 099 300  
Email: [info@zonemedical.com.au](mailto:info@zonemedical.com.au)  
7/22 Mavis Court, Yatala, Qld, 4207  
[www.zonemedical.com.au](http://www.zonemedical.com.au)