

ELECTRIC SUCTION UNIT with rechargeable battery

ASKIR 36BR



Antibacterial and Hydrophobic Filter

Silicone Autoclavable Tubes

LED indicators of battery status

Suction Regulator

Vacuum Indicator (kPa and Bar)

Collection Jar in Autoclavable PC with Overflow Valve System

available with disposable liners

Internal Rechargeable Battery

AC/DC Universal Adapter with Detachable Power Cord

12V Car Adapter



ASKIR 36BR is a desk-type electric suction unit for the aspiration of body liquids, oral, nasal and tracheal aspiration in adults or children. The appliance is designed for easy transport and almost continuous use. The unit is operated by means of an electronic power management system, also monitoring the internal battery operation and status. A LED system positioned on the front panel indicates appliance start-up and battery charge status. The unit is equipped with suction regulator on the front panel and polycarbonate autoclavable jar complete with overflow valve. Thanks to these characteristics and to its functions, this device is particularly suitable for different applications: utilization in hospital wards, for tracheotomy, for suction of body liquids and for minor surgery.

MAIN APPLICATIONS:

Patient Transfer

Tracheostomy

Emergency Dept.

CE 0123



SURGICAL HOUSE
Medical & Surgical Supply Specialists

CA-MI

TECHNICAL FEATURES

Motor _____	Oiless and maintenance-free piston pump
Power Feeding _____	<ul style="list-style-type: none"> ■ 14V --- 4A with AC/DC Universal Adapter mod. UE60-140429SPA3 (100-240V~ - 50/60Hz – 100VA) ■ Internal Rechargeable Pb Battery (12V --- 4A) ■ Car Adapter (12V --- 4A)
ISO 10079-1 Classification _____	HIGH VACUUM / HIGH FLOW
Max Vacuum (adjustable) _____	-0.80 Bar
	-80 kPa
	-600 mmHg
Power consumption _____	4A
Noise Level _____	65,5 dB (with jar) / 68,5 dB (without jar)
Max free air flow rate _____	36 l/m (without jar)
Weight _____	Kg 4.39
Size _____	cm 35 x 21 x 18
Years of Warranty _____	2
Shipping carton _____	4
Place of Manufacturing _____	Italy

STANDARD CONFIGURATIONS

code

Askir 36BR with 1000 ml autoclavable jar	RE 410200
Askir 36BR with 2000 ml autoclavable jar	RE 410200/01

STANDARD ACCESSORIES

code

1000 ml Autoclavable Jar + Cover (with overflow valve)	RE 210001
Antibacterial and hydrophobic Filter	SP 0121
8x14 mm Silicone Tubes	RE 210355/03
Conical Connector ø 10-11-12 mm	RE 210420
AC/DC Universal Adapter (100-240V~ - 50/60Hz)	SP 0208
Power cord with Shuko plug for universal adapter	SP 0020
12V Car Adapter	SP 0007/02
Internal Rechargeable Battery	NOT REPLACEABLE

OPTIONAL

code



2000 ml Autoclavable Jar + Cover (with overflow valve) RE 210351



Carrying bag for Unit and Consumables SP 0207



Stand with 5 castors 27730