

Compact

Aerosol delivery systems



COMPACT is a piston compressor nebulizer for aerosoltherapy. Indicated for home-care and intensive use, its design includes a convenient handle for easy transport of unit and a large compartment for storage of accessories and medications.

AVAILABLE CONFIGURATIONS

● RE300200

Compact

INCLUDED ACCESSORIES

HI-FLO Jet Nebulizer
Adult Masks with elastic strings
Child Masks with elastic strings
Mouthpiece
Nosepiece (non invasive)
Air Tube
Air Filter

FEATURES

Compressor	Oiless and maintenance-free piston pump
Power Feeding	230 V-50 Hz
Power consumption	170 VA
Max Pressure	3,0 bar
Operating Pressure	1,3 bar 18,9 psi 130 Kpa
Max Air Flow	16 L/min
Operating Air Flow	5,2 L/min
Neb-rate with 4 ml solution	0,5 ml/min
Mass Median Aerodynamic Diameter	3,25 μ m (MMAD)
Respirable Fraction (< 5 μ m)	> 70 %
Noise Level (EN 13544-1 Standard)	60 dB
Weight	1,60 Kg
Size (Main Unit)	cm 13,8x11,7x23,5
Years of Warranty	4
Place of Manufacturing	Italy

CA-MI srl - Via Ugo La Malfa 13 - Frazione Pilastro, 43013 Langhirano (PARMA) - Italy - P.IVA - CF - N° di Iscrizione c/o Registro Imprese di Parma: 00977090349



Le informazioni riportate nel sito www.ca-mi.it e relative a dispositivi medici non hanno alcuna natura pubblicitaria. Tutti i contenuti, in qualunque forma realizzati, (testi, immagini, anche fotografiche, descrizioni tecniche e non, ecc.), hanno natura esclusivamente informativa, e sono rivolte esclusivamente a operatori professionali. Circolare Ministero della Salute. Nuove linee guida della pubblicità sanitaria concernente i dispositivi medici, dispositivi medico-diagnostici in vitro e presidi medico chirurgici. 28 Marzo 2013

[Privacy Policy](#) - [Cookie Policy](#)