



BEDROOM DAY CHAIRS & COMMODOES

aspire for...

Independence
comfort
style



RANGE

High Back Classic Day Chair

Low Back Classic Day Chair

Waterfall Day Chair

Classic Bedside Commode

Deluxe Bedside Commode



www.aspirecare.com.au
1300 133 120



aspire HIGH BACK CLASSIC DAY CHAIR

aspire LOW BACK CLASSIC DAY CHAIR

OPTIONAL ACCESSORIES

- DVA Contracted
- AS / NZS Approved
- High Back Classic Day Chair is registered in the ARTG, AUST R ARTG 230087



Transport Legs with Wheels
CHP208090

A high back orthopaedic chair that offers adjustability, comfort and support



MODELS

- CHP208050 High Back Classic Day Chair - Champagne Vinyl
- CHP208055 High Back Classic Day Chair - Slate Vinyl

| | | | |
|---|--|---|--|
| 1 | Contoured Backrest – High back design provides increased postural support | 5 | Sturdy Design - Generous 160kg weight capacity, yet lightweight for easy manoeuvrability |
| 2 | Softly padded armrests, seat and back upholstery are comfortable and supportive | 6 | Powder Coated Steel Frame – Easy to clean, strong and durable |
| 3 | Easy clean vinyl surfaces | 7 | Non-Slip Tips – Provide stability whilst sitting or transferring |
| 4 | Height Adjustable – Easy push button adjustment for comfortable positioning and changing needs | 8 | Available in either Slate or Champagne Vinyl |



| SEAT DEPTH | OVERALL DEPTH | SEAT WIDTH | OVERALL WIDTH | SEAT HEIGHT | BACKREST HEIGHT | OVERALL HEIGHT | UNIT WEIGHT | S.W.L |
|------------|---------------|------------|---------------|-------------|-----------------|----------------|-------------|-------|
| 490mm | 600mm | 490mm | 645mm | 425-610mm | 720mm | 1155-1315mm | 16kg | 160kg |

- DVA Contracted
- AS / NZS Approved
- Low Back Classic Day Chair is registered in the ARTG, AUST R ARTG 230087

A versatile, low back orthopaedic chair



MODELS

- CHP208000 Low Back Classic Day Chair - Champagne Vinyl
- CHP208005 Low Back Classic Day Chair - Slate Vinyl

| | | | |
|---|---|---|--|
| 1 | Contoured Backrest – Provides postural support | 5 | Height Adjustable – Easy push button adjustment for comfortable positioning and changing needs |
| 2 | Softly padded armrests, seat and back upholstery are comfortable and supportive | 6 | Lightweight Design - lightweight yet sturdy design for easy manoeuvrability |
| 3 | Non-Slip Rubber Tips – Provide stability whilst sitting or transferring | 7 | Powder Coated Steel Frame – Easy to clean, strong and durable |
| 4 | Easy clean vinyl surfaces | 8 | Available in either Slate or Champagne Vinyl |



| SEAT DEPTH | OVERALL DEPTH | SEAT WIDTH | OVERALL WIDTH | SEAT HEIGHT | BACKREST HEIGHT | OVERALL HEIGHT | UNIT WEIGHT | S.W.L |
|------------|---------------|------------|---------------|-------------|-----------------|----------------|-------------|-------|
| 440mm | 600mm | 460mm | 580mm | 430-605mm | 440mm | 845-1015mm | 9kg | 160kg |

OPTIONAL ACCESSORIES

✓ DVA Contracted
✓ AS / NZS Approved
✓ Waterfall Day Chair is registered in the ARTG, AUSTRALIA ARTG 230097



MODELS

| | |
|-----------|--------------------------------------|
| CHP208070 | Waterfall Day Chair - Ink Vinyl |
| CHP208071 | Waterfall Day Chair - Latte Vinyl |
| CHP208075 | Waterfall Day Chair - Oatmeal Fabric |
| CHP208076 | Waterfall Day Chair - Mocha Fabric |

| | | | |
|---|--|---|--|
| 1 | Contoured, Waterfall Backrest – Segmented and adjustable backrest cushion provides improved ergonomics, comfort and postural support | 5 | Height Adjustable – Easy push button adjustment for comfortable positioning and changing needs |
| 2 | Padded Seat and Armrests - Generous support surfaces for better pressure distribution and sitting comfort | 6 | Non-Slip Tips - Provide stability whilst sitting or transferring |
| 3 | Removable Cushion Sections - The seat cushion and backrest overlay can be easily removed for cleaning | 7 | Powder Coated Steel Frame – Easy to clean, strong and durable |
| 4 | Sturdy Design - High safe working load of 200kg, yet lightweight for easy manoeuvrability | 8 | Available in either Vinyl (Ink or Latte) or Fabric (Oatmeal or Mocha) |

1 Year Warranty

| SEAT DEPTH | OVERALL DEPTH | SEAT WIDTH | OVERALL WIDTH | SEAT HEIGHT | BACKREST HEIGHT | OVERALL HEIGHT | UNIT WEIGHT | S.W.L |
|------------|---------------|------------|---------------|-------------|-----------------|----------------|-------------|-------|
| 490mm | 720mm | 535mm | 720mm | 420-580mm | 685mm | 1045-1200mm | 16.5kg | 200kg |

✓ DVA Contracted
✓ AS / NZS Approved
✓ Classic Bedside Commode is registered in the ARTG, AUSTRALIA ARTG 230171



MODELS

| | |
|-----------|---------------------------------------|
| BEC045005 | Classic Bedside Commode - Ink Vinyl |
| BEC045006 | Classic Bedside Commode - Latte Vinyl |

| | | | |
|---|---|----|---|
| 1 | Removable, Padded Seat – with lift-up toilet seat and fixed bowl carrier beneath | 6 | Easy clean vinyl surfaces |
| 2 | Raised Padded Armrests – Facilitate easier transfers and offer improved comfort | 7 | Non-Slip Rubber Tips – Provide stability whilst sitting or transferring |
| 3 | Removable Bowl with Lid and Carry Handle | 8 | Lightweight Design – Lightweight yet sturdy steel design for easy manoeuvrability (5.5kg) |
| 4 | Removable Back Rest – Enables easy storage or transportation | 9 | Three-sided Shroud – Discreetly conceals bowl |
| 5 | Height Adjustable – Easy push button adjustment allows for comfortable positioning and changing needs | 10 | Powder Coated Steel Frame – Strong, durable and easy to clean |
| | | 11 | Available in either Slate or Champagne Vinyl |

1 Year Warranty

| SEAT DEPTH | OVERALL DEPTH | SEAT WIDTH | OVERALL WIDTH | SEAT HEIGHT | BACKREST HEIGHT | OVERALL HEIGHT | UNIT WEIGHT | S.W.L |
|------------|---------------|------------|---------------|-------------|-----------------|----------------|-------------|-------|
| 440mm | 600mm | 450mm | 560mm | 475-615mm | 450mm | 900-1040mm | 5.5kg | 125kg |

✓ DVA Contracted
✓ AS / NZS Approved
✓ Deluxe Bedside Commode is registered in the ARTG, AUST R ARTG 230171



MODELS

- BEC045115 Deluxe Bedside Commode - Champagne Vinyl
- BEC045116 Deluxe Bedside Commode - Slate Vinyl

| | | | |
|---|--|----|--|
| 1 | Fold-Away, Padded Seat – Comfortable seat with lift-up toilet seat and fixed bowl carrier beneath | 6 | Easy clean vinyl surfaces |
| 2 | Raised Padded Armrests - Facilitate easier transfers and offer improved comfort | 7 | Non-Slip Tips - Provide stability whilst sitting or transferring |
| 3 | Removable Bowl with Lid and Carry Handle | 8 | Three-Sided Shroud – Discreetly conceals bowl |
| 4 | Contoured, Padded Back Rest – Provides postural support | 9 | Powder Coated Steel Frame – Strong, durable and easy to clean |
| 5 | Height Adjustable – Easy push button adjustment allows for comfortable positioning and changing needs. | 10 | Available in either Slate or Champagne Vinyl |

1
Year
Warranty

| SEAT DEPTH | OVERALL DEPTH | SEAT WIDTH | OVERALL WIDTH | SEAT HEIGHT | BACKREST HEIGHT | OVERALL HEIGHT | UNIT WEIGHT | S.W.L |
|------------|---------------|------------|---------------|---------------|-----------------|----------------|-------------|-----------|
| 450mm | 600mm | 460mm | 600mm | 410-565mm | 440mm | 865-1015mm | 12kg | 160kg |

BEDROOM RANGE SNAPSHOT



Bed Sticks & Poles

Lifters & Slings

Overbed Tables

Bedside Cabinets

AC Beds

Recline Chairs



ASSOCIATED PRODUCTS



LIFTERS & SLINGS



TRANSFER AIDS



DAILY LIVING AIDS



THERAPY & EXERCISE



CRUTCHES,
FRAMES & STICKS



SCOOTERS



OVERBED TABLES



BED STICKS & POLES



FOAM MATTRESSES



AIR MATTRESSES



URINALS & BED PANS



OVER TOILET AIDS



SHOWER STOOLS & CHAIRS



RAISED TOILET SEATS




SHOWER COMMODES



BATH BOARDS & BENCHES



SEAT WALKERS



WHEELCHAIRS



STD POWER WHEELCHAIRS



PRESSURE CUSHIONS




HOME CARE BEDS



LIFT & RECLINE CHAIRS

Distributed by:



ASBR08.17.100517

