

# Cordless Under-Mattress Pressure Pad



PART NUMBER: 3515

## FEATURES

- Under mattress for optimal patient comfort.
- Internal reinforcing for extra durability.
- Simple operation and maintenance
- For use with cura<sup>1</sup> Cordless Pad monitors
- No more cords between monitor and pads
- Easy to clean
- Virtually no false alarms over the lifetime of the pad
- Vent holes in the neck reduces the risk of moisture ingress

## TECHNICAL SPECIFICATIONS

- Pad Dimensions: 800mm x 290mm
- Recommended weight range: 40kg - 200kg
- Material: Non allergenic soft PVC
- Requires Cordless Pad Transmitter (2719)
- 1 Year guarantee

## ACCESSORIES

**2719** Cordless Pad Transmitter, **2502** LRCS Zone Repeater, **2501** LRCS Care Station, **2700** Six Channel Cordless Falls Monitor, **2608** Wireless Bed Monitor

cura<sup>1</sup>.com

cura<sup>1</sup> products are designed to augment well designed care programs and practices, cura<sup>1</sup> does not claim that any of its products are suitable as a substitute, or as a total solution for individuals that require care and constant supervision, and are at risk.