

# ActiveCare Pager



PART NUMBER: 3501

## FEATURES

- Large easy-to-read 7-digit LCD screen displays device type and ID.
- Can be programmed to pair with up to 9 wireless devices or up to 60 devices if used with the 3500 ActiveCare Console.
- Up to 5 simultaneous alarms can be stored in the pager's "Call Review" memory.
- Vibration and tone selection for hard of hearing.
- Out of range check.
- Wireless reset capability.
- Direct signalling from wireless devices for use in home care environments.
- Easily cleanable buttons for a hospital setting.

## TECHNICAL SPECIFICATIONS

- Operating frequency: 433.92MHz
- Typical operating range: 15-25m\*
- Power requirements: 2 x AA batteries
- Dimensions (LxWxD): 83mm x 55mm x 26mm
- Housing material: Durable ABS plastic
- Net weight: 65g (without batteries)
- Operating temp.: 5°C to 40°C
- Rel. humidity: <90%

## ACCESSORIES

**2719** Cordless Pad Transmitter, **3500** ActiveCare Console, **2588** PET, **2558** Waterproof Call Button, **2552** DoorAlert Door Strip, **3502** Wireless Handset, **3503** Magnetic Reed Switch

cura<sup>1</sup>.com

\*The range of wireless devices may vary and be dependent on your specific environmental conditions. cura<sup>1</sup> products are designed to augment well designed care programs and practices, cura<sup>1</sup> does not claim that any of its products are suitable as a substitute, or as a total solution for individuals that require care and constant supervision, and are at risk.