

Wireless Call Buttons with Caregiver Pager



PART NUMBER: 4024

FEATURES

- Ideal for independent living and residential care
- Creates independence for busy carers and security for those in their care
- Super-simple plug and play setup and operation
- Low cost battery operated (batteries included)
- Kit contains two transmitters and one pager
- Pager emits an alarm sound when the button on a transmitter is pressed
- Transmitters can be wall mounted or worn as a pendant around the neck
- Pager includes belt clip

TECHNICAL SPECIFICATIONS

- Transmitters use 1 x A23 12 V battery (included)
- Pager uses 2 x AA 1.5 V batteries (included)
- Operating frequency 433 MHz
- Range adequate for most conventional home environments*
- Durable ABS plastic construction
- 2 selectable alarm tones and off switch on the pager
- Green LED alarm indicator light

cura¹.com

*The range of wireless devices may vary and be dependent on your specific environmental conditions. cura¹ products are designed to augment well designed care programs and practices, cura¹ does not claim that any of its products are suitable as a substitute, or as a total solution for individuals that require care and constant supervision, and are at risk.