



Pressure Care Cushions

by Apollo Healthcare Technologies

CLASSIC - MED

The Classic-Med cushion is designed for comfort and pressure reduction. The castellated foam cells conform to the body, redistributing pressure across the surface area. Ideally recommended for those At Risk to Medium Risk of pressure ulcer development, including light weight patients.

A soft multi stretch, water-resistant cover provides protection from fluid ingress. Welded seams and a protective zip flap support infection control requirements. Easy to wipe down or remove and launder.

Features	Benefits
Caters for Patients up to Medium Risk	At Risk to Medium Risk
Comfort for Light Weight Patients	Soft Density Foam Redistributes Pressure from Bony Prominence Allows Further Immersion and Contour to Body Shape
Castellated Foam	Allows Further Immersion and Contour to Body Shape
Multi-Stretch Soft Cover	Soft on the Skin Conforms with the Foam
High Frequency Welded Seams and Water-Resistant Cover	Reduce Risk of Fluid Ingress and Easy to Wipe Clean
Rear Zip with Welded Protective Flap	Removable Cover, Launder up to 95 Degrees
Acti-Fresh Biocide Anti-Bacterial Agent	Integral Anti-Bacterial Protection



Foam Specification			
Foam Component	36/100 Beige	Design	Single Castellated Foam
Nominal Density	34-37 kg/m ³	Nominal Density	110-140N

Product Information			
Code	APH010/APH011	Risk Level	Medium Risk
		Weight Capacity	108Kg
Cushion Size	APH010 - 43 x 43 x 10cm	Product Weight	2Kg
	APH011 - 51 x 51 x 10cm	Warranty	2 Years
Safety Standards	BS EN 597-1 Cig, BS EN 597-2 Match, BS7177: 1996 Crib Source 5		

CLASSIC-MED CUSHION