


SIGVARIS

DOFF N' DONNER

The revolutionary donning
and removal device for medical
compression stockings.



SURGICAL HOUSE

Medical & Surgical Supply Specialists

Properties and benefits

- Facilitates the donning and removal of compression stockings and promotes treatment compliance.
- Less effort spent on donning and removal.
- Allows compression stockings to slide painlessly over wounds, scars and open legs.
- Suitable for use by the wearers themselves, as well as caregivers.
- Is extremely durable and suitable for different leg shapes.
- Extends the life of the compression stockings.
- Suitable for models A–D calf stockings and A–G thigh stockings with/without grip top in all sizes and compression classes as well as with open or closed toe.

A: DOFF N' DONNER

The DOFF N' DONNER (DND) is a very soft device that allows you to quickly and painlessly slide compression stockings over an arm or a leg. It is suitable for use by the wearers themselves, as well as caregivers.

B: CONE

The CONE simplifies the rolling of a compression stocking onto the DND. With the aid of the CONE, the medical compression stocking can also be easily unrolled from the DND. The CONE has a suction cup built into the base so that it attaches well to flat surfaces.

- simplest
- fastest
- most effortless donning aid

A



B



SURGICAL HOUSE
Medical & Surgical Supply Specialists

46 KING EDWARD ROAD, OSBORNE PARK 6017

PO BOX 1537, OSBORNE PARK DC WA 6916

surgical@surgicalhouse.com.au

TEL. +61 8 9381 4199

FAX. +61 8 9382 3009

www.surgicalhouse.com.au

How to Use the Doff N Donner- [Click here for video](#)