

3 Section Couches

3 section couch medical

Product Description

- Electric Hi Lo operation by foot switch
- Sleek design with rounded frame legs for safety
- Minimum height 50cm to maximum height 90cm in 30 seconds
- Adjustable head section featuring double gas strut, angle adjustable from -45° to +70°
- Adjustable backrest featuring double gas struts, angle adjustable from -45° to +70°
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm castors for safety and mobility
- Retractable dual locking castors on all 4 corners
- T Motion motor
- Padding thickness 5cm
- Foam Density: 36g/m²

Certification: CE, FDA, ISO9001-2008, ISO13485-2007

Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant

Specification:

- Total Length: 190cm
- Width: 71cm
- Head section length: 62cm
- Middle section length: 62cm
- Leg section length: 62cm
- Height: adjustable from 50cm to 90cm

Loading weight: 250Kgs

ORDER CODE:

ET33NB - NAVY BLUE

ET33BL - BLACK

