

HESTIA CONTOUR

The Hestia Contour cushion is designed to contour to the patient body, effectively redistributing pressure. This unique design improves posture and provides a stable pressure reducing support surface. The combination of improved posture and increased support surface area will contribute significantly to pressure ulcer prevention. Ideally recommended for those At Risk to High Risk of pressure ulcer development.



A soft, water-resistant multi stretch cover provides protection from fluid ingress. The cover is easy to wipe down or remove and launder.

FEATURES		BENEFITS	
Caters for Patients up to High Risk		At Risk to High Risk	
Moulded Foam Cushion		Maintains Positioning and Stability	
Castellated Foam		Allows Further Immersion and Contour to Body Shape	
Integral Sag Compensator		Supports Sacral Area Maintaining Comfort and Position	
Anti-Slip Base		Reduces Shear and Friction Forces	
High Frequency Welded Seams and Water-Resistant Cover		Reduce Risk of Fluid Ingress, Easy to Wipe Clean	
Rear Zip with Welded Protective Flap		Removable Cover, Launder up to 95 Degrees	
Acti-Fresh Biocide Anti-Bacterial Agent		Integral Anti-Bacterial Protection	
FOAM SPECIFICATION			
Foam Component	Visco Elastic Foam (PU)	Design	Single Mould
Nominal Density	86kg/m3	Nominal Density	
PRODUCT INFORMATION			
Code		RISK LEVEL	HIGH RISK
Cushion Size (Sizes available on request)	45 x 45 x 10 cm	Weight Capacity	40 - 120kg
Warranty	2 years	Product Weight	1.5kg
Safety Standards			