

# INSTILLAGEL®

STERILE LUBRICANT available with and without chlorhexidine

SETTING BEST  
PRACTICE STANDARDS  
FOR STERILE MEDICAL  
LUBRICATION

MATCHING THE NEEDS  
OF THE PATIENT WITH A  
LUBRICANT THAT PROTECTS  
WHILST PERFORMING


## INSTILLAGEL® LIDO

- Pain relief through the anaesthetising effect of 2% Lignocaine Hydrochloride.
- Excellent lubrication properties of hydroxymethylcellulose base adheres especially well to the urethra.
- Long-lasting lubrication due to the non-drying effects of added propylene glycol.
- 11ml syringe for that extra bit of lubrication.
- Now available as a certified and registered medical device.

## INSTILLAGEL® ALMED


- All the benefits of Instillagel® Lido PLUS 0.05% Chlorhexidine for patients without Chlorhexidine sensitivity.
- Two added hydroxybenzoates clinically analysed to have extra efficacy against bacteria<sup>2</sup>.
- Superior adhesion to the mucosa provides an antiseptic shield before, during and after the procedure.
- Clinically evaluated effective against MRSA<sup>1</sup>.
- Closely monitored low allergy history over 48 years.

## ACCESSORIES TO ENHANCE YOUR PRACTICE



### PENILE CLAMP

- The best solution for allowing anaesthetic lubricants to take effect – recommended 5 min for peak anaesthetic performance.
- Sterile, single use.
- Minimise patient discomfort – allowing patient discretion and dignity.
- Optimise work flow.



### INSTILLAQUILL®

- An extension quill facilitating the use of Instillagel® in gynaecological investigations.
- Allows delivery of anaesthetic gel into the cervix and uterus.
- Excellent for use in hysteroscopy and insertion of intra-uterine devices.
- An additional benefit of hysterosonography - once gel has been instilled in the endometrial cavity for 10–15 minutes it facilitates instrumentation of the uterus<sup>3</sup>.

1. Efficacy of antiseptic catheter lubricants against methicillin-resistant staphylococcus aureus [MRSA]. Der Urologe [A] 44 (2005) 282-285 Published online 1st March 2005. P. Faul.  
2. Local antiseptic therapy for non-gonococcal urethritis: phase II study with Instillagel®. World J Urol 2002. 20: 45-47. Published online 23rd November. H. Sperling.  
3. Experience with Instillagel® for hysterosonography and analgesia in a complex contraception clinic: a QIPP initiative. JFPRHC Published online 23 Oct 2011. Mary Pillai, Julia Shefras

# INSTILLAGEL®

**STERILE LUBRICANT**

available with and without chlorhexidine

## WHAT ARE THE DIFFERENCES WITH WHAT I ALREADY USE?



### BETTER FORMULA TO SUIT YOUR USE

Whether it's the added hydroxybenzoates clinically analysed to have extra efficacy against bacteria with the Chlorhexidine Almed version or the added propylene glycol to maintain non-drying effects in our standard Lido syringe hydroxymethylcellulose base, we have the sterile anaesthetic lubricant suited to your patient.

In the simple test of thickness, our formula is less viscous; designed to spread more easily up the urethra and coat the entire surface. It is designed to adhere especially well to the mucosa and last through a standard transurethral resection of over 2000 resectoscope movements. It contains no fragment particles or micro bubbles and makes an excellent ultrasound medium.

### NOT A PRESCRIPTION MEDICATION

Instillagel® is registered as a medical device, NOT a pharmaceutical and is therefore not required to be bought, stored or controlled by a Pharmacist or supplied through a pharmacy.

Instillagel® can be ordered from wholesale channels and supply departments. It does not need special licensing to store, supply or use the product as an S2 pharmaceutical syringe does, nor does it need to be prescribed. This means less paperwork and no liability for the safety requirements of handling, storing and dispensing a medication.



### STERILE PACKAGING THAT CONFORMS TO ASEPTIC TECHNIQUE

Instillagel® comes in an autoclaved blister pack for aseptic technique procedures without any assembly of the syringe prior to use. This conforms to ACORN guidelines for maintaining a sterile field:

*"All items introduced onto a sterile field should be opened, dispensed and transferred by methods that maintain sterility and integrity. Edges of all sterile wrappers and packages are considered nonsterile. Therefore, when opening an item, be sure all the edges are secured so that they do not contaminate the field or the sterile packaged item."*



### SYRINGE TIP SPECIFICALLY DESIGNED FOR LESS TRAUMA ON ENTRY

The Instillagel® tip is designed to fit easily into the urethra with the base of the syringe acting as a stopper against backflow during delivery. The Instillagel® tip is highly polished and smooth unlike other product tips which can feel sharp and rough due to being produced using a mould. It won't cause damage to the meatus or cause unnecessary discomfort when the syringe is inserted and positioned correctly.

For patient safety, it's important to know that the tip of the Instillagel® syringe has also been designed not to fit a needle so there is no risk of injecting the solution by mistake.

### PRODUCT ORDERING DETAILS

INSTILLAGEL® ALMED	10ml, 25 per box	FARCO 343367
INSTILLAGEL® LIDO	11ml, 10 per box	FARCO 350914
PENILE CLAMP	Sterile, single use, 10 per box	FARCO 337007
INSTILLAQUILL®	Sterile, single use, 10 per box	CLINI 40-200



Phone: AU 1300 857 230 Email: [sales@mdsanz.com](mailto:sales@mdsanz.com)  
Phone: NZ 0800 994 951 [www.instillagel.com.au](http://www.instillagel.com.au)

