

# Disposable Resuscitator Liberty



Disposable Bag Mask Resuscitators are made with a lightly-textured PVC bag for sure grip and comfort and a transparent re-breathing valve assembly allows for clear viewing of valve function. Available in adult, pediatric, and infant sizes



## Features

- Meets ISO-8382 and ASTM0920 standards
- Meets ISO-5356-2 standard
- Textured surface for slip free operation in wet conditions
- Faster rebound and less memory material applied for ease of operation
- Ergo designed shape provides ease of handling and less access force
- Clear resuscitation bag for ease of observation
- High efficiency silicone duckbill patient valve minimizes back leakage
- Provides low inspiratory and expiratory resistance
- Swivel 15/22mm patient connector
- 40cmH<sub>2</sub>O pressure relief for child and infant models
- 60cmH<sub>2</sub>O pressure relief available with or without lock function for adult model
- 32mm reservoir connection for lower intake resistance
- High efficiency reservoir system provides up to 99% oxygen concentration at low flow
- Oxygen tubing connection for ease of access
- Inflatable clear air cushion mask provides better seal, ease of observation during access and also patient comfort



Code	Description	Mask	Reservoir	Pop-Off
LRES3403RI	LIBERTY ADULT DISPOSABLE RESUSCITATOR	#5	1000 ML	60CM H2O
LRES3403EI	LIBERTY ADULT DISPOSABLE RESUSCITATOR	#5	1000 ML	NO POP-OFF
LRES0484QI	LIBERTY CHILD DISPOSABLE RESUSCITATOR	#3	1000 ML	40CM H2O
LRES0486QI	LIBERTY INFANT DISPOSABLE RESUSCITATOR	#0	1000 ML	40CM H2O