



# DISPOSABLE MANUAL RESUSCITATOR ADULT



## Features

- Meets ISO-8382 and ASTM0920 standards.
- Arrow shape texture surface for slip free operation in wet conditions.
- Faster rebound and less memory material applied for ease of operation.
- Ergo designed shape provides ease of handling and less access force.
- Clear resuscitation bag for ease of observation.
- High efficiency silicone duckbill patient valve minimises back leakage.
- Provides low inspiratory and expiratory resistance.
- Swivel 15/22mm patient connector.

- Meets ISO-5356-2 standard.
- 60 cm H<sub>2</sub>O pressure relief with lock function for adult model.
- High efficiency reservoir system provides up to 99% oxygen concentration at low flow.
- Swivel elbow oxygen tubing connection for ease of access. Fits-all oxygen connector can be attached to flow meter directly.
- Clear air cushion mask provides better seal, ease of observation during access and also patient comfort.



0043

CODE	DESCRIPTION
2268	DISPOSABLE MANUAL RESUSCITATOR ADULT
0043	SEE-THROUGH CARRY CASE ADULT

MASK	RESERVOIR	POP-OFF
5305	1000ML	60 CM H <sub>2</sub> O

E&OE First Published: 30 June 2009  
 Last Modified: 13 October 2011  
 File Name: info\_2268\_ADULT-DISPOSABLE-RESUSCITATOR

