


Technical Specifications

	150F	A205	200
Image			
Product Code	LSS390600	LSS390610	LSS390620
Safe Working Load (S.W.L)	150kg	205kg	200kg
Base Length	1280mm (1220mm folded)	1135mm	1280mm
External Base Width	660-1040mm (615mm folded)	620-1095mm	665-980mm
Internal Base Width	580-950mm	547-1015mm	580-900mm
Base Width Adjustment	Manual	Electric	Electric
Base Height / Under Bed Clearance	133mm	115mm	110mm
Minimum Height of Spreader Bar **	428mm	395-495mm	535mm
Maximum Height of Spreader Bar**	1640mm	1705-1805mm	1725mm
Lifter Weight	32.9kg	38.1kg	43.9kg
Under Base Height	72mm	50mm	50mm
Maximum Overall Height	1873mm (395mm folded)	2045mm	1940mm
Front Castor Diameter	75mm	75mm	75mm
Rear Locking Castor Diameter	100mm	100mm	100mm
Turning Radius	1365mm	1235mm	1385mm
Weight Of Heaviest Part (Base Assembly)	32.9kg	19.6kg	23.0kg
Handle Height	1000mm	1050-1150mm	1050mm
Warranty - Frame	5 Years	5 Years	5 Years
Warranty - Electronics	3 Years	3 Years	3 Years
Warranty - Battery*	1 Year	1 Year	1 Year
ACCESSORIES			
Digital Weigh Scale (S.W.L 300kg)	LSS390675	LSS390675	LSS390675
Pivot Frame (S.W.L 200kg)	LSS390670	LSS390670	LSS390670

* Includes Battery Pack and Handset

** Three most height settings are possible hence the range of minimum and maximums

A320	A200S	200S
		
LSS390630	LSS390640	LSS390650
320kg	200kg	200kg
1505mm	985mm	1150mm
765-1190mm	590-915mm	655-1016mm
640-1065mm	490-815mm	575-936mm
Electric	Electric	Electric
155mm	110mm	110mm
575-675mm	890mm	890mm
1860-1960mm	1395mm	1740mm
63kg	30.4kg	44.9kg
85mm	50mm	50mm
2165mm	1455mm	1790mm
100mm	75mm	75mm
125mm	75mm	100mm
1590mm	1100mm	1255mm
38.9kg	16.6kg	25.1kg
1200-1300mm	1000mm	1030mm
5 Years	5 Years	5 Years
3 Years	3 Years	3 Years
1 Year	1 Year	1 Year
LSS390675	Not Available	Not Available
LSS390670	Not Available	Not Available