

PATIENT STANDER - ALUMINIUM

The Aidacare Aspire A200S is a compact, lightweight and manoeuvrable standing lifter ideal for home, community and facility environments.



ADJUSTABLE DEPTH YOKE

Yoke is depth adjustable for optimal patient fit



ERGONOMIC HANDLES

Provide multiple grip points to suit carers of different heights and strengths



ERGONOMIC HANDSET

The latest slimline technology from Linak fits neatly in the hand with an intuitive control panel. Battery and service indicator lights assist with service scheduling




POWERED LEG FUNCTION

Opens and closes legs for lifting stability or accessibility




Need a Sling?

Find the correct sling for your manual handling scenario on page 16



Quick Tech Specs



Lifter Weight: 30.4kg

S.W.L: 200kg

LSS390640 See page 15 for detailed specs

BUILD QUALITY

Strong powder coated aluminium frame with industry leading Linak electronics and warranty provision



SAFETY

Critical features include emergency stop button, emergency lowering and anti-entrapment technology, with both patient and caregiver safety at heart



ADJUSTABLE PADDED LEGRESTS

Provide comfort and support during transfers and can be adjusted to suit the patient



INTELLIGENT DIAGNOSTICS

State of the art Linak Jumbo Care control system. Tracks usage and service requirements for safety and compliance



QUALITY CASTORS

2 x 75mm front castors for maximum under bed clearance and accessibility.
2 x 75mm rear locking castors for safety and manoeuvrability

