

# Treatment Tables - ECONOMY

## MedBed Professional Treatment Table

The MedBed 3-section treatment table is designed to provide a comprehensive range of features and offer exceptional value for money. The design is extremely robust, long lasting and provides a safe and ergonomic design for treating patients. Easy operation via footswitch to raise or lower the table. Retractable wheels for easy mobility.

Available in two styles:

- **Medical** - no face hole
- **Physiotherapy** - with face hole and plug


Specifications	Medical	Physiotherapy
Minimum Height	370mm	370mm
Maximum Height	773mm	773mm
Cushion Width	700mm	700mm
Head Section Length	740mm	740mm
Mid Section Length	550mm	550mm
Leg Section Length	640mm	640mm
Total Length	1997mm	1997mm
Head Tilt Angle	-38°/80°	-38°/80°
Leg Tilt Angle	-20°/30°	-20°/30°
Lift Mechanism	Electric	Electric
Breathing Hole	No	Yes
Gas Structure head and foot	WDF gas struct	WDF gas struct
Retractable Wheels	Standard x 4	Standard x 4
User Weight Rating	250kg	250kg

### Warranty

Frame	5 Years
Motor	3 Years
Controller	2 Years
Strut, Footswitch, & Upholstery	1 Year




### Dimensions:



Physiotherapy Medical

### Table Positions:



Code	Description
------	-------------

<b>MED-ECO70NFH</b>	3-Section Medical No Face Hole
<b>MED-ECO70FH</b>	3-Section Physiotherapy With Face Hole and Plug