

Noninvasive, Accurate, Reliable, Easy to Use





Palm Bladder Scanner

Noninvasive, Accurate,

PBS

Palm Bladder Scanner is a handheld non invasive device with advance 3-D Ultrasound technology. With the help of its advance feature it can easily differentiate between bladder and seated area. It is quite easy to measure the bladder volume in the bladder with ease and comfort.

PBS-V3.I

Feature

- · Light weight, easy to carry
- 8.0 color LCD TFT screen
- · Use Friendly software
- · Simple operation
- Store and review the imagine and patient data
- · Measurement data accurate
- · With USB interface can send the report to PC
- · Software can upgrade through internet
- · Built-in lithium battery

Interface Language:

English, French, Germany, Spanish and Chinese

Usb Ports:

for printing and upgrade the software



Weight: 1.8kg (N/W) 2.5kg (G/W)

Reliable, Easy to Use



Specification

Transducer freq: 2.6MHz

Scan method: 3D Mechanical sector scanning

Scan depth: ≥140mm

Sway angle: 120°±2°

Rotating angle: 180°±2°

• Display screen: 8.0" Color LCD TFT screen

Power supply: Battery

• Patient Case storage: 100pcs

Patient ID: 12 digits

• Max. Measurement: 999ml

· Scan mode: male; female; PRV



PBS 2



Application

PBS is widely used in hospital
Urology, Obstetrics and Gynecology,
Anorectic
Rehabilitation, Emergency
Department, General Surgery

Pediatrics, ICU

PBS's latent application area:

Appertain area, Uric incontinence guide, SCI/MS

Mentality healing institution, Rest home

Family tend:

SCI/MS, Community Medical Center, Family and Community Doctor.





Mianyang Meike Electronic Equipment Co., Itd

Add: No.238, Chuangye Road, kechuangyuan Zone, Mianyang City, Sichuan, China

Noninvasive, Accurate, Reliable, Easy to Use

