

Memo-Form Plus™ **Static Pressure Reducing Mattress**

Product Code - APH005



The **Memo-Form Plus** pressure reducing mattress features a high density visco-elastic foam surface inlay that is sat in a “U” shaped core, ideal for delivering exceptional patient comfort and a high degree of pressure reduction. The temperature sensitive foam conforms to the individual’s shape, creating a greater body contact area and effective pressure redistribution. Designed to meet the demands of the acute healthcare environment the **Memo-Form Plus** is suitable for those considered at ‘**High Risk**’ of developing pressure ulcers.



- High specification foam contributes to an effective pressure ulcer prevention plan
- Internal air channels in the foam's surface help to reduce the build up of heat and moisture on the skin
- Smooth multi-stretch cover allows the patient to glide over the surface of the mattress
- Individually castellated cells conform independently
- The waterproof, vapour permeable cover can be easily wiped clean
- High frequency welded seams reduce the risk of fluid ingress
- Full length double sided flap cover can be removed and laundered at 95°C

Specifications

Dimensions: 198cm x 88cm x 15.6cm (bespoke sizes available on request)

Maximum Weight Limit: 39 stone / 247.7 Kg

Risk Factor: High Risk

Warranty: 5 years cover / 5 years foam

Mattress Weight: 16 Kg

Safety Standards: BS EN 597-1 Cig, BS EN 597-2 Match, BS7177: 1996 Crib Source 5

