


## MEMO-PAD BARI

The Memo-Pad Bari cushion is designed to provide a high degree of comfort and pressure care management for heavier patients. Utilising memory foam technology, this cushion gently conforms to the body distributing pressure across the surface area. Ideally recommended for those At Risk to High Risk of pressure ulcer development.

A soft multi stretch, water-resistant cover provides protection from fluid ingress. Welded seams and a protective zip flap support infection control requirements. Easy to wipe down or remove and launder.



FEATURES		BENEFITS	
Bariatric Cushion for Patients up to High Risk		At Risk to High Risk	
Low Resistant Memory Foam		Conforms to the Body, Evenly Redistributing Pressure	
High Specification Foam		Provides Stability and Support	
Soft and Smooth Protective Cover		Soft on the Skin, Reducing Shear Force During Transfer	
Water-Resistant, Vapour Permeable Cover		Easy to Wipe Clean and Reduced Moisture Build Up	
High Frequency Welded Seams and Water-Resistant Cover		Reduce Risk of Fluid Ingress, Easy to Wipe Clean	
Rear Zip with Welded Protective Flap		Removable Cover, Launder up to 95 Degrees	
Acti-Fresh Biocide Anti-Bacterial Agent		Integral Anti-Bacterial Protection	
FOAM SPECIFICATION			
Layers	Top	Middle	Bottom
Foam Component			
Nominal Density	kg/m3	kg/m3	kg/m3
Nominal Hardness	N	N	N
			
PRODUCT INFORMATION			
Code	APH073	<b>RISK LEVEL</b>	<b>HIGH RISK</b>
Cushion Size (Sizes available on request)	Made to Order	Weight Capacity	254.5kg
Warranty	2 years	Product Weight	4kg
Safety Standards	BS EN 597-1 Cig, BS EN 597-2 Match, BS7177: 1996 Crib Source 5		