



Pressure Care Cushions by Apollo Healthcare Technologies

Memo-Pad

The Memo-Pad cushion is designed to provide a high degree of comfort and pressure care management. Utilising memory foam technology, this cushion gently conforms to the body distributing pressure across the surface area. Ideally recommended for those At Risk to High Risk of pressure ulcer development, including light weight patients.

A soft multi stretch, water-resistant cover provides protection from fluid ingress. Welded seams and a protective zip flap support infection control requirements. Easy to wipe down or remove and launder

Features	Benefits
Caters for Patients up to High Risk	At Risk to High Risk
Low Resistant Memory Foam	Conforms to the Body, Evenly Redistributing Pressure
High Specification Foam	Provides Stability and Support
Soft and Smooth Protective Cover	Soft on the Skin, Reducing Shear Force During Transfer
High Frequency Welded Seams and Water-Resistant Cover	Reduce Risk of Fluid Ingress and Easy to Wipe Clean
Rear Zip with Welded Protective Flap	Removable Cover, Launder up to 95 Degrees
Acti-Fresh Biocide Anti-Bacterial Agent	Integral Anti-Bacterial Protection

Layers	Top	Middle	Bottom
Foam Component	V60 Yellow Visco	V60 Blue Visco	36/100 White
Nominal Density	58-62 kg/m ³	58-62 kg/m ³	34-37 kg/m ³
Nominal Hardness	70-100 N	100-130 N	90-120 N

Product Information			
Product Codes and Sizes	APH014: 43 x 43 x 5cm	Risk Level	High Risk
	APH012: 43 x 43 x 7.6cm	Weight Capacity	152Kg
	APH015: 51 x 51 x 5cm	Product Weight	3Kg
	APH013: 51 x 51 x 7.6cm	Warranty	2 Years
Safety Standards	BS EN 597-1 Cig, BS EN 597-2 Match, BS7177: 1996 Crib Source 5		

MEMO-PAD CUSHION