

# MÖLNLYCKE - WOUND CARE PRODUCTS

## PRODUCT DESCRIPTION

## AREAS OF USE

## SIZE

## CODE

## PCS/BOX

### BORDERED FOAM DRESSINGS



#### Mepilex® Border Flex / Oval

Absorbent self-adherent foam dressing with 360° stretch Flex Technology



SafetaC  
TECHNOLOGY

- Moderate to heavily exuding wounds
- Oval shape for challenging body contours
- Diabetic foot ulcers
- Leg and foot ulcers
- Traumatic wounds, such as skin tears
- Pressure injury prevention\* and treatment

\* When a repositioning/offloading protocol is used

**Square**  
7.5cm x 7.5cm  
10cm x 10cm  
12.5cm x 12.5cm  
15cm x 15cm  
15cm x 20cm

**Square**  
595211  
595311  
595011  
595411  
595611

**Square**  
10  
10  
10  
10  
10

**Oval**  
7.8cm x 10cm  
13cm x 16cm  
15cm x 19cm

**Oval**  
583500  
583300  
583400

**Oval**  
5  
5  
5



#### Mepilex® Border Heel

Absorbent self-adherent foam dressing for the heel area

Added benefits:

- Thicker borders to minimise rolling and to allow efficient handling
- Optimised seal for enhanced protection of the high-risk heel area
- Enhanced coverage of high risk areas
- Additional handling tabs for easier skin checks




SafetaC  
TECHNOLOGY

- Moderate to heavily exuding wounds for the heel area
- Heel pressure injury and treatment (all stages)
- Diabetic foot ulcers
- Leg ulcers
- Traumatic wounds

22cm x 23cm

282750

6



#### Mepilex® Border Sacrum

Absorbent self-adherent foam dressing for the sacral area

Added benefits:

- Thicker borders to minimise rolling and to allow efficient handling
- Optimised gluteal seal for enhanced protection of the high-risk sacral area
- Enhanced coverage of high risk areas
- Additional handling tabs for easier skin checks




SafetaC  
TECHNOLOGY

- Moderate to heavily exuding wounds for the sacral area
- Sacral pressure injury prevention and treatment (all stages)

16cm x 20cm  
22cm x 25cm

282050  
282450

5  
5



#### Mepilex® Border Flex Lite

Thin absorbent self-adherent foam dressing with 360° stretch Flex Technology



SafetaC  
TECHNOLOGY

- Non- to moderately exuding wounds
- Leg and foot ulcers
- Pressure injuries
- Surgical and traumatic wounds, e.g. abrasions, blisters and skin tears
- Protection of compromised and/or fragile skin

4cm x 5cm  
5cm x 12.5cm  
7.5cm x 7.5cm  
10cm x 10cm  
15cm x 15cm

581011  
581100  
581200  
581300  
581500

10  
5  
5  
5  
5

# MÖLNLYCKE - WOUND CARE PRODUCTS

## PRODUCT DESCRIPTION

## AREAS OF USE

## SIZE

## CODE

## PCS/BOX

### NON-BORDERED FOAM DRESSINGS

	<p><b>Mepilex® XT</b> Highly conformable absorbent foam dressing, with specially designed exudate channels to assist in the absorption of thick exudate</p>  	<ul style="list-style-type: none"> <li>• Low to moderately exuding wounds in all healing stages</li> <li>• Chronic wounds with thick exudate</li> <li>• Leg ulcers and diabetic foot ulcers</li> <li>• Pressure injuries</li> <li>• Traumatic wounds</li> <li>• Can be used under compression</li> </ul>	<p>5cm x 5cm 10cm x 10cm 10cm x 20cm 15cm x 15cm 20cm x 20cm 20cm x 50cm</p>	<p>211015 211100 211200 211300 211400 211500</p>	<p>5 5 5 5 5 2</p>
	<p><b>Mepilex®</b> Highly conformable and absorbent foam dressing</p>  	<ul style="list-style-type: none"> <li>• Low to moderately exuding wounds</li> <li>• Leg ulcers and diabetic foot ulcers</li> <li>• Pressure injuries</li> <li>• Traumatic wounds</li> <li>• Partial thickness burns</li> <li>• Can be used under compression</li> </ul>	<p>5cm x 5cm 10cm x 10cm 10cm x 20cm 15cm x 15cm 20cm x 20cm 20cm x 50cm</p>	<p>294015 294100 294200 294300 294400 294500</p>	<p>5 5 5 5 5 2</p>
	<p><b>Mepilex® Heel</b> Absorbent foam dressing for the heel area</p>  	<ul style="list-style-type: none"> <li>• A wide range of exuding wounds on the heel or malleolus</li> <li>• Foot ulcers / diabetic foot ulcers</li> </ul>	<p>13cm x 20cm</p>	<p>288100</p>	<p>5</p>
	<p><b>Mepilex® Lite</b> Thin absorbent foam dressing</p>  	<ul style="list-style-type: none"> <li>• Non- to low exuding wounds</li> <li>• Diabetic foot ulcers</li> <li>• Leg ulcers</li> <li>• Radiation skin reactions</li> <li>• Epidermolysis Bullosa</li> </ul>	<p>6cm x 8.5cm 10cm x 10cm 15cm x 15cm 20cm x 50cm</p>	<p>284000 284100 284300 284500</p>	<p>5 5 5 4</p>

### WOUND CONTACT LAYERS

	<p><b>Mepitel® One</b> Wound contact layer dressing with Safetac® on one side</p>  	<ul style="list-style-type: none"> <li>• Skin tears</li> <li>• Lacerations</li> <li>• Surgical incisions</li> <li>• Fixation of partial and full thickness grafts</li> <li>• Partial thickness burns</li> <li>• Diabetic foot ulcers</li> <li>• Venous and arterial leg ulcers</li> </ul> <p>** Sizes with pre-printed symbols</p>	<p>5cm x 7.5cm 6cm x 7cm** 7.5cm x 10cm 9cm x 10cm** 13cm x 15cm** 10cm x 18cm 17cm x 25cm 24cm x 27.5cm 27.5cm x 50cm</p>	<p>289100 289170 289300 289270 289470 289500 289700 289670 289750</p>	<p>10 5 10 5 5 10 5 5 5</p>
	<p><b>Mepitel®</b> Wound contact layer dressing with Safetac® on both sides</p>  	<ul style="list-style-type: none"> <li>• Skin tears</li> <li>• Surgical incisions</li> <li>• Partial thickness burns</li> <li>• Fixation of partial and full thickness grafts</li> <li>• Diabetic foot ulcers, venous and arterial leg ulcers</li> </ul>	<p>5cm x 7.5cm 7.5cm x 10cm 10cm x 18cm 20cm x 30cm</p>	<p>290510 290710 291010 292005</p>	<p>10 10 10 5</p>
	<p><b>Mepilex® Transfer</b> Exudate transfer foam dressing</p>  	<ul style="list-style-type: none"> <li>• Highly exuding wounds</li> <li>• Difficult-to-dress areas</li> <li>• Radiation skin reactions</li> <li>• Epidermolysis Bullosa</li> </ul>	<p>7.5cm x 8.5cm 10cm x 12cm 15cm x 20cm 20cm x 50cm</p>	<p>294600 294700 294800 294502</p>	<p>5 5 5 4</p>

# MÖLNLYCKE - WOUND CARE PRODUCTS

## PRODUCT DESCRIPTION





## AREAS OF USE

## SIZE

## CODE

## PCS/BOX



### POST-OPERATIVE DRESSINGS

	<p><b>Mepilex® Border Post-Op</b> Self-adherent post-operative dressing with 360° stretch Flex Technology</p> 	<ul style="list-style-type: none"> <li>• Low to heavily exuding wounds</li> <li>• Optimised for post-operative use and blood absorption</li> <li>• Surgical wounds, cuts and abrasions</li> <li>• Ideal for use over joints, such as hips and knees</li> </ul>	6cm x 8cm	496100	10
			9cm x 10cm	496200	10
			10cm x 15cm	496300	10
			10cm x 20cm	496400	10
			10cm x 25cm	496450	10
			10cm x 30cm	496600	10
10cm x 35cm	496650	5			
	<p><b>Mepore® Film&amp;Pad</b> Transparent, shower-proof absorbent island dressing</p> 	<ul style="list-style-type: none"> <li>• For a wide variety of wounds which have low to moderate exudate levels</li> <li>• Post-operative wounds</li> <li>• Excisions and sutures</li> </ul>	9cm x 15cm	275500	30
			9cm x 20cm	275600	30
			9cm x 25cm	275700	30
			9cm x 35cm	275900	25

### FIBRE / ALGINATE DRESSINGS

	<p><b>Exufiber®</b> Gelling fibre dressing with Hydrolock® Technology</p> 	<ul style="list-style-type: none"> <li>• Leg ulcers and diabetic foot ulcers</li> <li>• Pressure injuries</li> <li>• Partial thickness burns</li> <li>• Surgical wounds</li> <li>• Donor sites</li> <li>• Malignant wounds</li> </ul>	5cm x 5cm	709900	10
			10cm x 10cm	709901	10
			15cm x 15cm	709903	10
			20cm x 30cm	709904	5
			4.5cm x 20cm	709906	10
			1cm x 45cm (Cavity)	709908	5
			2cm x 45cm (Cavity)	709909	5
	<p><b>Melgisorb® Plus</b> Highly absorbent alginate dressing</p> 	<ul style="list-style-type: none"> <li>• Moderate to heavily exuding wounds</li> <li>• Pressure injuries</li> <li>• Venous and arterial leg ulcers</li> <li>• Diabetic foot ulcers</li> <li>• Donor sites</li> <li>• Post-operative wounds</li> <li>• Dermal lesions and other external wounds inflicted by trauma</li> </ul>	5cm x 5cm	252000	10
			10cm x 10cm	252200	10
			10cm x 20cm	252500	10
			3cm x 45cm (Cavity)	253500	5

### ACTIVE DRESSING

	<p><b>Mesalt®</b> Sodium chloride impregnated dressing</p> 	<ul style="list-style-type: none"> <li>• For the management of heavily draining and infected wounds in the inflammatory phase</li> <li>• Deep cavity wounds</li> </ul>	5cm x 5cm	285580	30
			7.5cm x 7.5cm	285780	30
			10cm x 10cm	286080	30
			2cm x 1m (Cavity)	285280	10

# MÖLNLYCKE - WOUND CARE PRODUCTS

## PRODUCT DESCRIPTION

## AREAS OF USE

## SIZE

## CODE

## PCS/BOX

### ANTIMICROBIAL DRESSINGS

	<p><b>Mepilex® Border Ag</b> Antimicrobial absorbent self-adherent foam dressing</p>  	<ul style="list-style-type: none"> <li>Moderate to heavily exuding wounds</li> <li>Infected wounds</li> <li>Leg ulcers</li> <li>Pressure injuries</li> <li>Malignant wounds</li> <li>Surgical incisions</li> </ul>	<p>7.5cm x 7.5cm 10cm x 10cm 15cm x 15cm 15cm x 20cm <b>Surgical Sizing</b> 10cm x 20cm 10cm x 25cm 10cm x 30cm</p>	<p>395200 395300 395400 395600 395800 395700 395900</p>	<p>5 5 5 5 5 5 5</p>
	<p><b>Mepilex® Border Sacrum Ag</b> Antimicrobial absorbent self-adherent foam dressing for the sacral area</p>  	<ul style="list-style-type: none"> <li>Moderate to heavily exuding wounds for the sacral area where a self-adherent antimicrobial dressing is required</li> </ul>	<p>18cm x 18cm 23cm x 23cm</p>	<p>382000 382400</p>	<p>5 5</p>
	<p><b>Mepilex® Ag</b> Antimicrobial absorbent foam dressing</p>  	<ul style="list-style-type: none"> <li>Low to moderately exuding wounds</li> <li>Partial thickness burns</li> <li>Infected wounds</li> <li>Leg ulcers and foot ulcers</li> <li>Pressure injuries</li> <li>Can be used under compression</li> </ul>	<p>10cm x 10cm 10cm x 20cm 15cm x 15cm 20cm x 20cm 20cm x 50cm</p>	<p>287110 287210 287310 287410 287510</p>	<p>5 5 5 5 2</p>
	<p><b>Mepilex® Heel Ag</b> Antimicrobial absorbent foam dressing for the heel area</p>  	<ul style="list-style-type: none"> <li>Low to moderately exuding wounds</li> <li>Infected wounds on the foot or heel</li> <li>Diabetic foot ulcers</li> <li>Pressure injuries</li> <li>Can be used under compression</li> </ul>	<p>13cm x 20cm</p>	<p>388100</p>	<p>5</p>
	<p><b>Mepilex® Transfer Ag</b> Antimicrobial compressed foam, exudate transfer dressing</p>  	<ul style="list-style-type: none"> <li>Heavily exuding wounds</li> <li>Infected or malignant wounds</li> <li>Skin grafts and donor sites</li> <li>Partial thickness burns</li> <li>Leg ulcers and diabetic foot ulcers</li> <li>Difficult-to-dress areas</li> </ul>	<p>7.5cm x 8.5cm 10cm x 12.5cm 12.5cm x 12.5cm 15cm x 20cm 20cm x 50cm</p>	<p>394000 394100 394700 394800 394500</p>	<p>10 5 5 10 2</p>
	<p><b>Melgisorb® Ag</b> Highly absorbent antimicrobial alginate dressing</p> 	<ul style="list-style-type: none"> <li>Moderate to heavily exuding wounds</li> <li>Post-operative wounds</li> <li>Leg ulcers and diabetic foot ulcers</li> <li>Trauma wounds (dermal lesions, trauma injuries or incisions)</li> <li>Graft and donor sites</li> <li>Pressure injuries</li> </ul>	<p>5cm x 5cm 10cm x 10cm 15cm x 15cm 20cm x 30cm 3cm x 44cm (Cavity)</p>	<p>256055 256105 256155 256205 256605</p>	<p>10 10 10 5 10</p>
	<p><b>Exufiber® Ag+</b> Antimicrobial gelling fibre dressing with Hydrolock® Technology</p> 	<ul style="list-style-type: none"> <li>Moderate to heavily exuding wounds</li> <li>Venous leg ulcers</li> <li>Diabetic foot ulcers</li> <li>Pressure injuries</li> </ul>	<p>5cm x 5cm 10cm x 10cm 15cm x 15cm 20cm x 30cm 4.5cm x 20cm 2cm x 45cm (Cavity)</p>	<p>603401 603402 603403 603407 603405 603400</p>	<p>10 10 10 5 10 5</p>

# MÖLNLYCKE - WOUND CARE PRODUCTS

## PRODUCT DESCRIPTION





## AREAS OF USE

## SIZE




## CODE

## PCS/BOX

### SUPERABSORBENTS

	<b>Mextra® Superabsorbent</b> Superabsorbent dressing with a fluid-repellent backing 	<ul style="list-style-type: none"> <li>Moderate to very heavily exuding wounds</li> <li>Can be used as a primary or secondary dressing where high absorbency is required, e.g. together with Mepitel®, Mepitel One or Mepilex® Transfer</li> <li>Use with Tubifast® for light retention</li> </ul>	12.5cm x 12.5cm	610000	10
			12.5cm x 17.5cm	610100	10
			12.5cm x 22.5cm	610200	10
			17.5cm x 22.5cm	610300	10
			22.5cm x 27.5cm	610400	10
			22.5cm x 32.5cm	610500	10
22.5cm x 42.5cm	610600	10			
	<b>Mesorb®</b> Highly absorbent dressing with a fluid-repellent backing 	<ul style="list-style-type: none"> <li>Moderate to very heavily exuding wounds</li> <li>Can be used as a primary or secondary dressing where high absorbency is required, e.g. together with Mepitel, Mepitel One or Mepilex Transfer</li> <li>Use with Tubifast for light retention</li> </ul>	10cm x 13cm	677001	50
			13cm x 15cm	677201	50
			10cm x 23cm	677401	50
			15cm x 23cm	677601	50
			23cm x 25cm	677701	30
			23cm x 30cm	677801	30

### FILM DRESSINGS

	<b>Mepitel® Film</b> Breathable, transparent film dressing with Safetac® Technology 	<ul style="list-style-type: none"> <li>Protects fragile and sensitive skin</li> <li>For fixation of primary dressings</li> <li>For use under aggressive tapes</li> <li>For use in radiotherapy</li> <li>Can be used as a protective dressing over IV sites</li> </ul>	6cm x 7cm	296170	10
			10cm x 12cm	296270	10
			10cm x 25cm	296470	10
			15cm x 20cm	296670	10
	<b>Mepore® Film</b> Breathable, transparent film dressing	<ul style="list-style-type: none"> <li>Designed as a secondary fixation dressing and as a primary dressing for a wide range of clean wounds</li> </ul>	6cm x 7cm	270600	100
			10cm x 12cm	271500	70
			10cm x 25cm	272500	20
			15cm x 20cm	273000	10

### ISLAND DRESSINGS

	<b>Mepore®</b> Absorbent island dressing 	<ul style="list-style-type: none"> <li>For a wide variety of superficial wounds which have low to moderate levels of exudate</li> <li>Post-operative wounds</li> <li>Cuts and grazes</li> </ul>	5cm x 7cm	672000	60
			6cm x 7cm	670800	60
			8cm x 10cm	672100	50
			9cm x 10cm	670900	50
			8cm x 15cm	672200	50
			9cm x 15cm	671000	50
			9cm x 20cm	671100	30
			9cm x 25cm	671200	30
			9cm x 30cm	671300	30
9cm x 35cm	671400	30			
	<b>Mepore® Pro</b> Shower-proof, absorbent island dressing 	<ul style="list-style-type: none"> <li>For a wide variety of superficial wounds which have low to moderate levels of exudate</li> <li>Where a shower-proof dressing is required</li> <li>Post-operative wounds</li> <li>Cuts and grazes</li> </ul>	6cm x 7cm	670820	60
			9cm x 10cm	670920	40
			9cm x 15cm	671020	40
			9cm x 20cm	671120	30
			9cm x 25cm	671220	30
9cm x 30cm	671320	30			
	<b>Mepore® Film&amp;Pad</b> Transparent, shower-proof absorbent island dressing 	<ul style="list-style-type: none"> <li>For a wide variety of wounds which have low to moderate exudate levels</li> <li>Where a shower-proof dressing is required</li> <li>Cuts and abrasions</li> <li>Post-operative wounds</li> </ul>	4cm x 5cm	275100	85
			5cm x 7cm (Oval)	275200	85
			5cm x 7cm (Square)	275300	85
			9cm x 10cm	275400	30
			9cm x 15cm	275500	30
			9cm x 20cm	275600	30
			9cm x 25cm	275700	30
			9cm x 35cm	275900	25

# MÖLNLYCKE - WOUND CARE PRODUCTS

## PRODUCT DESCRIPTION

## AREAS OF USE

## SIZE

## CODE

## PCS/BOX

### BANDAGES

	<p><b>Tubifast®</b> WITH 2-WAY STRETCH TECHNOLOGY Lightweight tubular retention bandage</p>	<ul style="list-style-type: none"> <li>• Lightweight tubular dressing retention bandage that is ideal for any part of the body</li> <li>• Can be used for patch wrapping and as an undercast stockinette</li> <li>• Can be used under compression</li> <li>• For use as skin covering</li> <li>• Epidermolysis Bullosa</li> </ul>	<p><b>10m roll</b> <b>RED</b> 3.5cm x 10m <b>GREEN</b> 5cm x 10m <b>BLUE</b> 7.5cm x 10m <b>YELLOW</b> 10.75cm x 10m <b>PURPLE</b> 20cm x 10m</p>	<p>2434 2436 2438 2440 2444</p>	<p>1 1 1 1 1</p>
	<p><b>Tubifast® Garments</b> Ready-to-wear garments</p>	<ul style="list-style-type: none"> <li>• Ready-to-wear garments for dressing retention and wet or dry wrapping, are available in full-sleeved vests, tights, leggings, socks and gloves</li> <li>• For children of six months to 14 years of age</li> <li>• With no tapes, pins or ties, Tubifast Garments are quick and easy to apply, and their 2-Way Stretch construction ensures complete freedom of movement, so aiding compliance.</li> <li>• Comfortable to wear under nightwear and ordinary clothes</li> </ul>	<p><b>Vests</b> 6-24 months 992007 2-5 yr 992008 5-8 yr 992009 8-11 yr 992010 11-14 yr 992011 <b>Tights</b> 6-24 months 992012 <b>Leggings</b> 2-5 yr 992013 5-8 yr 992014 8-11 yr 992015 11-14 yr 992016 <b>Socks</b> 1 size: 2-14yr 992017 <b>Gloves</b> Small child 5923</p>	<p>6 6 6 6 6 6 6 6 6 6 6</p>	<p>6 6 6 6 6 6 6 6 6 6 6</p>
	<p><b>Tubigrip®</b> Multi-purpose, elasticated tubular support bandage that provides firm support for sprains, strains and swelling</p>	<p><b>Size / Limb Circumference (cm)</b> A 10 - 14.5 B 14.5 - 29 C 22 - 35 D 26 - 40 E 29 - 46 F 35 - 53 G 40 - 62 J 54 - 82 K 70 - 113 L 106 - 149 M 119 - 167</p>	<p><b>Natural (N)</b> 1435 1436 1443 1437 1434 1438 1439 1440 1441 1442 1444</p>	<p><b>Beige (B)</b> - 1449 1450 1451 1448 1452 1453 - - - - -</p>	<p><b>N B</b> 1 - 1 1 1 1 1 1 1 1 1 1 1 1 1 - 1 - 1 - 1 - 1 -</p>
	<p><b>Elset® 'S'</b> Light compression bandage</p>	<ul style="list-style-type: none"> <li>• Used to provide compression for the leg following vascular surgery and for general surgical and orthopaedic use</li> <li>• Specifically recommended for above and below knee stump bandaging</li> </ul>	<p><b>Elset 'S'</b> 10cm x 12m 15cm x 12m</p>	<p>993486 993487</p>	<p>12 12</p>
	<p><b>Setopress®</b> Lightweight high compression bandage</p>	<ul style="list-style-type: none"> <li>• Used in the treatment of venous leg ulcers and associated conditions requiring graduated compression of the leg</li> </ul>	<p>10cm x 3.5m</p>	<p>3505</p>	<p>12</p>

# MÖLNLYCKE - WOUND CARE PRODUCTS

## PRODUCT DESCRIPTION

## AREAS OF USE

## SIZE

## CODE

## PCS/BOX

### FIXATION



#### Mepitac®

Polyurethane tape with Safetac® Technology



- For the fixation of medical devices such as drains, tubes and probes
- Skin protection under devices

2cm x 3m  
4cm x 1.5m

298300  
298400

1  
1



#### Mefix®

Self-adhesive fabric fixation tape

- For fixation of a wide variety of dressings, tubes and cannulas

**10m roll**  
x 2.5cm  
x 5cm  
x 10cm  
x 15cm  
x 20cm  
x 30cm

310299  
310599  
311099  
311599  
312000  
313000

1  
1  
1  
1  
1  
1

### SCAR THERAPY & DERMATOLOGY



#### Mepiform®

Scar dressing with Safetac Technology



- Medical grade silicone dressing for the prevention and treatment of all scars, in particular keloid scars and/or hypertrophic scars

4cm x 30cm  
5cm x 7.5cm  
10cm x 18cm

293100  
293200  
293400

5  
5  
5



#### Epaderm® Cream

An emollient for dry skin conditions

- A 2-in-1 emollient that can be used as a skin emollient, as well as a skin cleanser
- Used for eczema, psoriasis and all other dry skin conditions
- Suitable for all ages, including babies

25g  
50g  
500g

99400822  
99400851  
99400850

12  
12  
6



#### Epaderm® Ointment

An emollient for dry skin conditions

- A 3-in-1 emollient that can be used as a skin emollient, a skin cleanser and a bath additive
- Used for eczema, psoriasis and all other dry skin conditions
- Suitable for all ages, including babies

25g  
125g  
500g  
1kg

400861  
99400846  
99400847  
99400818

12  
6  
6  
6

### WOUND IRRIGATION



#### Granudacyn®

Wound irrigation solution and gel for cleansing and moisturising acute, chronic and contaminated wounds

- All chronic wounds of any depths, such as diabetic foot ulcers, pressure injuries, venous leg ulcers, etc.
- All acute wounds, such as cuts, bite wounds, lacerations, abrasions, etc.
- Surgical wounds (intra-operative and post-operative)
- Wounds with exposed cartilage, tendons, ligaments and/or bones
- Superficial and partial thickness burns
- Radiation ulcers
- Fistulas and abscesses

**Solution**  
50ml spray  
250ml spray  
500ml  
1000ml  
**NPWT Solution**  
500ml  
1000ml  
**Gel**  
50g  
100g spray

360150  
360100  
360101  
360102  
360103  
360104  
360107  
360108

1  
1  
1  
1  
1  
1  
1  
1

# MÖLNLYCKE - WOUND CARE PRODUCTS

PRODUCT DESCRIPTION

AREAS OF USE

SIZE

CODE

PCS/BOX

## TURNING & POSITIONING SYSTEMS



### Mölnlycke® Z-Flex™ Heel Boot

To maintain offloading of heels overtime. The combination of positive air displacement and fluidised positioning helps to offload the heel and protects the vulnerable Achilles area.

- Fully offloads the heel
- Positions feet in neutral alignment
- Accommodates use of compression devices
- Built-in gate allows quick access for heel inspection
- Thermal regulating material
- One size fits most

2 boots per box  
8 boots per box

1400122  
1400123

2  
8



### Mölnlycke® Z-Flo™ Fluidized Positioner

To assist with patient positioning and offloading of at risk anatomical sites

- Redistributes pressure or offloads
  - Keeps patients in position
  - Mouldable to protect vulnerable areas
- \* More sizes available for Paediatric and Neonatal

**With Tab**  
40cm x 76cm  
**Without Tab**  
64cm x 91cm  
30cm x 51cm  
(Occipital)

1401003  
1401005  
1401007

1  
1  
8



### Mölnlycke® Tortoise® Turning & Positioning System

To assist with turning and repositioning of patients

- Facilitates patient turning, positioning, boosting and lateral transfers
- Helps maintain the intended degree of turn and protects bony prominences
- Reduces risk of caregiver injury


Standard  
Bariatric

1400800  
1400801

1  
1

## NEGATIVE PRESSURE WOUND THERAPY

PCS/BOX PCS/TRP



### Avance® Solo

A portable single patient use negative pressure system for closed surgical incisions and open wounds

- -125mmHg negative pressure for up to 14 days
- Visual and audible alarms for blockage, leakage and low battery
- Excess fluid transferred into a 50ml replaceable canister
- CFM (R) Technology ensures continuous negative pressure
- Absorbent multilayer dressing with Safetac® Technology
- Wide variety of dressing sizes

10 x 20 cm\*  
10 x 30 cm\*  
10 x 35 cm\*  
Pump & canister kit  
10 x 20 cm  
10 x 30 cm  
10 x 35 cm  
15 x 15 cm  
15 x 20 cm  
15 x 30 cm  
20 x 20 cm  
25 x 25 cm  
50ml  
1.5 x 10 x 12cm  
Nightcase

8881020  
8881030  
8881035  
8880050  
881020  
881030  
881035  
881515  
881520  
881530  
882020  
882525  
880050  
882000  
5105100-00

1  
1  
1  
1  
2  
2  
2  
2  
2  
2  
2  
2  
2  
4  
4  
1



### Avance® Solo Adapt

A portable single patient use NPWT system for pressure injuries

- -125mmHg negative pressure for up to 14 days
- Visual and audible alarms for blockage, leakage and low battery
- Excess fluid transferred into a 50ml replaceable canister
- CFM (R) Technology ensures continuous negative pressure
- single use polyurethane foam wound filler which can be cut to the shape of the wound
- Avance Solo Adpat Film with Safetac® Technology
- Wide variety of dressing sizes

'Starter Set'  
Film Dressing Kit  
Film Dressing  
Night-case

9994152  
994152  
994000  
5105100-00

1  
1  
2  
1

NB: Please refer to the Instructions of Use prior to usage.



SCAN ME

Visit **Mölnlycke ANZ YouTube channel** or **scan the QR** for our product application guides.

Find out more at [www.molnlycke.com.au](http://www.molnlycke.com.au)

You will be required to input the following code **1FtyrV** or your **AHPRA registration number** to access our websites.

The Mölnlycke, Avance, Elset, Epaderm, Exufiber, Granudacyn, Hydrolock, Mefix, Melgisorb, Mepiform, Mepilex, Mepitac, Mepitel, Mepore, Mesalt, Mesorb, Mextra, Safetac, Setopress, Tubifast, Tubigrip, 2-Way Stretch, Tortoise names and logos are registered globally to one or more of the Mölnlycke Health Care Group of Companies. Z-Flex and Z-Flo are trademarks used by Mölnlycke Health Care Group of Companies. ©2023 Mölnlycke Health Care. All rights reserved. [www.molnlycke.com.au](http://www.molnlycke.com.au). Mölnlycke Health Care, 12 Narabang Way, Level 4, Belrose, NSW, 2085. Phone 1800 005 231. New Zealand Orders & Enquiries 0800 005 231. [www.molnlycke.co.nz](http://www.molnlycke.co.nz). AUWC\_000899

