

Personal Emergency Transmitter (PET)



PART NUMBER: **2588**

FEATURES

- Large rubber button for easy use by elderly patients.
- Standard Lanyard for patient to wear around neck.
- Works with single use hospital wristbands.
- Changable battery with projected 1-year life.
- Splash resistant.

TECHNICAL SPECIFICATIONS

- Operating frequency 433.92MHz
- Operating Range: 30m line of sight
- Power Supply: 12V 27A Alkaline Battery
- Battery Life: >1 year (depends on usage)
- Dimensions: 42mm diameter, 12mm thick
- Housing Material: ABS
- Operating Temp: 20-50°C
- Rel. Humidity: <95%

ACCESSORIES

3500 ActiveCare Console, **3501** ActiveCare Pager, **2700** Cordless 6 Channel Monitor, **2544** Wireless Nurse-Call Upgrade Kit, **2589** Wrist Strap