




PROVOX® Swab

High quality swab for easy and hygienic tube cleaning



Atos
Atos Medical Your voice



- Ideal to clean tube surface and prevent mucus or obstructive build-up within the tube
- Strong, sturdy, smooth swab head, designed to clean the surface without scratching
- Robust flexible handle for ease of movement inside the tube
- Disposable swabs in 5 packs of 10 to help reduce contamination risks, which is important in the hospital environment

Provox Swab	Ref. no	Qty.	Recommended for tube sizes (I.D.)
Provox Swab Small, 5 mm	8250	50	4-5
Provox Swab Medium, 8 mm	8251	50	5-7
Provox Swab Large, 11 mm	8252	50	8-10
Provox Swab XtraLarge, 13 mm	8258	50	11+, Provox LaryTubes, Provox LaryButtons

I.D. = inner diameter

Provox Swab has been designed for ex-situ cleaning of tracheostomy tubes, laryngectomy tubes and stoma buttons. It consists of a head made from a smooth polyurethane foam and a plastic (polypropylene) handle.