

STANDARD PRODUCT PERFORMANCE

based on a medium sized glove

| | NITRILE | LATEX | VINYL |
|-------------------------|-------------|-------------|-------------|
| | ASTM D 6319 | ASTM D 3578 | ASTM D 5250 |
| Measurements | | | |
| Palm thickness | ≥ 0.05 mm | ≥ 0.08 mm | ≥ 0.08 mm |
| Finger tip thickness | ≥ 0.05 mm | ≥ 0.08 mm | ≥ 0.05 mm |
| Length | ≥ 230 mm | ≥ 230 mm | ≥ 230 mm |
| Tensile strength | | | |
| Before aging | ≥ 14 MPa | ≥ 18 MPa | ≥ 11 MPa |
| After aging | ≥ 14 MPa | ≥ 14 MPa | ≥ 11 MPa |
| Elongation | | | |
| Before aging | ≥ 500% | ≥ 650% | ≥ 300% |
| After aging | ≥ 400% | ≥ 500% | ≥ 300% |



Medicom®

SafeBasics™

**NEW
GLOVE
LINE**

The essentials you need. The quality you expect.
Value that can't be beat.

Healthcare professionals will not compromise on quality when it comes to protection. There are many options available, products with added benefits, with varying levels of protection and with many features. However, there are times when a more basic approach to disposables and every day consumables is all that is required. For Healthcare professionals who seek great quality in their essentials from a brand name they trust, Medicom introduces a new line that provides the protection you expect, at a value that can't be beat.

Our commitment is to help to make the world safer and healthier by offering consistent, reliable protection. Medicom is proud to introduce SafeBasics, a line of products with the protection you insist on. The first step in a complete line of product essentials, Medicom is proud to introduce a new family of gloves focused on delivering the best fit at the best price. With options available in latex, nitrile and vinyl, we have a glove to fit all hands and meet all needs, from a brand name you trust.

As always Medicom is driven by a commitment to quality, so Medicom SafeBasics gloves are engineered to provide the ultimate fit, are made from high quality materials and are designed to ensure protection that you can always count on.

Medicom SafeBasics, the perfect combination of protection, fit and price.



Medicom® Trust in Quality

Like you, we never compromise. Our commitment is to help to make the world safer and healthier by offering consistent, reliable protection. We use carefully selected materials and engineer our products to offer protection you can count on.

As the maker of Safe brands, we take pride in the protection we offer, so you can take pride in yours.

www.medicom.com

Medicom® Pride in Protection.

The perfect combination of protection, fit and price



Think Green with economy style packaging, available in 200 and 300 count boxes, which reduces transport costs and environmental waste to increase your savings

Medicom® SafeBasics™ True Fit™



NITRILE for best protection

For peace of mind. The unrivaled strength of nitrile delivers enhanced barrier protection. Designed to fit like a glove should, while providing the tactile sensitivity needed for the most complex procedures.

1185 Medicom SafeBasics True Fit Thin BLUE Textured Nitrile Medical Examination Gloves

- | | |
|--------------------|--------------------|
| 1185-A Extra Small | 1185-D Large |
| 1185-B Small | 1185-E Extra Large |
| 1185-C Medium | |

Maximum savings with
300 gloves per box

Packaging:
300 gloves/box; 8 boxes/case



1186 Medicom SafeBasics True Fit Thin PINK Textured Nitrile Medical Examination Gloves

- | | |
|--------------------|---------------|
| 1186-A Extra Small | 1186-C Medium |
| 1186-B Small | 1186-D Large |

Maximum savings with
300 gloves per box

Packaging:
300 gloves/box; 8 boxes/case



1187 Medicom SafeBasics True Fit BLACK Textured Nitrile Medical Examination Gloves

- | | |
|--------------------|--------------------|
| 1187-A Extra Small | 1187-D Large |
| 1187-B Small | 1187-E Extra Large |
| 1187-C Medium | |

Save more with
200 gloves per box

Packaging:
200 gloves/box; 10 boxes/case

Product specifications compliant with current standard ASTM D 6319

Medicom® SafeBasics™ Easy Fit™



LATEX for best fit

Made from all-natural rubber latex, these gloves conform perfectly to the hand for a fit like no other. Natural rubber latex offers superior barrier protection to protect you and your patients. Textured fingertips enhance grip and improve control.

1188 Medicom SafeBasics Easy Fit NATURAL Powder Free, Textured Latex Medical Examination Gloves

- | | |
|--------------------|--------------------|
| 1188-A Extra Small | 1188-D Large |
| 1188-B Small | 1188-E Extra Large |
| 1188-C Medium | |

Packaging:
100 gloves/box; 10 boxes/case



1189 Medicom SafeBasics Easy Fit NATURAL Powdered, Latex Medical Examination Gloves

- | | |
|--------------------|--------------------|
| 1189-A Extra Small | 1189-D Large |
| 1189-B Small | 1189-E Extra Large |
| 1189-C Medium | |

Save more with
200 gloves per box

Packaging:
200 gloves/box; 10 boxes/case

Product specifications compliant with current standard ASTM D 3578

Medicom® SafeBasics™ Smart Fit™



VINYL for best value

The right fit at the right price. 100% Latex free, these gloves are ideal for users with latex allergies. Delivers proven protection for all your healthcare needs. Smooth finish for increased tactile sensitivity and wearer comfort.

1194 Medicom SafeBasics Smart Fit CLEAR Powder-free, Vinyl Medical Examination Gloves

- | | |
|---------------|--------------------|
| 1194-B Small | 1194-D Large |
| 1194-C Medium | 1194-E Extra Large |

Save more with
200 gloves per box

Packaging:
200 gloves/box; 10 boxes/case



1195 Medicom SafeBasics Smart Fit CLEAR Powdered, Vinyl Medical Examination Gloves

- | | |
|---------------|--------------------|
| 1195-B Small | 1195-D Large |
| 1195-C Medium | 1195-E Extra Large |

Save more with
200 gloves per box

Packaging:
200 gloves/box; 10 boxes/case

Product specifications compliant with current standard ASTM D 5250