

**solari**  
MEDICAL

COMFORT & SUPPORT

# PATIENT CARE BEDS & FURNITURE



AVAILABLE IN  
A VARIETY  
OF COLOURS  
& STYLES



## COMFORT AND SUPPORT

Solari Medical specialises in providing quality products to the aged and health care industries of Australia.

Manufactured to Australian Standards our extensive product range aims to provide the highest level of comfort and support.

# *The* Adelphi



AVAILABLE IN  
A VARIETY  
OF COLOURS

# The Adelphi

WARRANTY  
10 YEARS  
STEEL FRAME

WARRANTY  
3 YEARS  
ELECTRONICS



## SPECIFICATIONS

MAXIMUM USER WEIGHT: 200 KILOS  
CENTRAL LOCKING

Overall Length	2150mm
Overall Width	930mm
Lift Range	270 - 660 mm
Mattress Platform	860 x 2000mm

**AVAILABLE IN:**  
**KING SINGLE SIZE**  
AND WITH  
**REMOVABLE BED EXTENSION**

## WALL BUFFER




## CENTRAL LOCKING



## SIDE RAILS (OPTIONAL)



## POSTURE SETTINGS



## BEDSIDE CABINET

FOUR CASTORS WITH BRAKES  
LOCKABLE TOP DRAW

Width	550mm
Depth	450mm
Height	780mm

# The Como



AVAILABLE IN  
A VARIETY  
OF COLOURS

## The Como

WARRANTY  
10 YEARS  
STEEL FRAME

WARRANTY  
3 YEARS  
ELECTRONICS



### SPECIFICATIONS

MAXIMUM USER WEIGHT: 200 KILOS  
CENTRAL LOCKING

Overall Length	2150mm
Overall Width	930mm
Lift Range	270 - 660 mm
Mattress Platform	860 x 2000mm

### AVAILABLE IN:

**KING SINGLE SIZE**

AND WITH

**REMOVABLE BED EXTENSION**

### WALL BUFFER




### CENTRAL LOCKING



### SIDE RAILS (OPTIONAL)



### POSTURE SETTINGS



### BEDSIDE CABINET

FOUR CASTORS WITH BRAKES  
LOCKABLE TOP DRAW

Width	550mm
Depth	450mm
Height	780mm

# *The* Turino

Floor Line Bed



AVAILABLE IN  
A VARIETY  
OF COLOURS

# The Turino

## Floor Line Bed



### SPECIFICATIONS

MAXIMUM USER WEIGHT: 200 KILOS

Overall Length	2245mm
Overall Width	920mm
Lift Range	160 - 557 mm
Mattress Platform	860 x 2000mm

**AVAILABLE IN:**  
**KING SINGLE SIZE**  
 AND WITH  
**REMOVABLE BED EXTENSION**

### WALL BUFFER




### CENTRAL LOCKING



### SIDE RAILS (OPTIONAL)



### POSTURE SETTINGS



### BEDSIDE CABINET

FOUR CASTORS WITH BRAKES  
 LOCKABLE TOP DRAW

Width	550mm
Depth	450mm
Height	780mm

COMFORT AND SUPPORT

