

---

# Oxford<sup>®</sup> Standing Harness Slings







# Standing Harness

The Oxford® Standing Harness is designed specifically to assist with standing and walking activity, supporting the patient around the upper body as well as under the arms. The Standing Harness is not a general-purpose sling and demands a degree of weight bearing ability from the client.


- ✓ Colour coded binding indicates sling size
- ✓ Coloured attachment loops promote correct/safe positioning
- ✓ Padded body and harness straps for comfort
- ✓ Fully adjustable belt and straps to meet patient needs
- ✓ Can be customised to fulfil patient specific needs



Chair to Chair



Bed to Chair/Chair to Bed  
(From/to seated position)



Chair to Commode/Toilet



Rehab/Assisted Standing



Size	Part No.	SWL
XS	SL26600	227 kg
S	SL26601	227 kg
M	SL26602	227 kg
L	SL26603	227 kg
XL	SPECIAL	227 kg





At Enable we source quality, affordable and innovative products globally and have them available for fast delivery. But we do more than just provide products.

We give the knowledge and support to help our customers grow their business and deliver better life care.



**ENABLE LIFECARE PTY LTD**

**T 1300 370 370**

**[enablelifecare.com.au](http://enablelifecare.com.au)**