




SupraCath™

Two-way suprapubic catheters with integrated balloon

The SupraCath range from Endotherapeutics offers several features to enhance patient comfort and allow maximum drainage through a smaller-sized puncture site.

The SupraCath Advantage

- Soft, medical-grade silicone is hypoallergenic and less prone to encrustation.
- Depth measurements ensure accurate placement during catheter exchange.
- Integrated balloon enables insertion through a smaller introducer and facilitates smooth removal by minimising cuffing difficulties.
- The open end reduces blockage and enables over-the-wire exchange of catheters.
- Large side eyes allow maximum drainage.
- Short tip reduces bladder irritation and is ideal for long-term drainage.



Available Sizes

- Sizes 10 Fr through to 24 Fr
- 14 Fr, 16 Fr, 18 Fr and 20 Fr in standard and short lengths



NEW!

All SupraCath sizes now come with a convenient luer lock spigot at no extra cost. The spigot can be easily removed if another valve to manage flow is preferred.



Endotherapeutics

Suite 1, Level 1, 51 Rawson Street
Epping NSW 2121, Australia

T: 1300 00 ENDO (3636) F: + 61 2 9869 7723 E: info@endotherapeutics.com.au
www.endotherapeutics.com.au

A-Cysta-Cath™

Suprapubic Catheter Introducer

SupraCath catheters are compatible with introducers up to **two sizes smaller**:



Compatibility Chart

18 Fr SupraCath	passes through	16 Fr A-Cysta-Cath
14 Fr SupraCath		12 Fr A-Cysta-Cath

Description	Product Code
14 Fr A-Cysta-Cath Plus 12 Fr A-Cysta-Cath, 14 Fr SupraCath & Scalpel	YM-1034-14
18 Fr A-Cysta-Cath Plus 16 Fr A-Cysta-Cath, 18 Fr SupraCath & Scalpel	YM-1034-18



SupraCath™ Specifications

Product Code	Size (Fr)	Balloon Capacity (cc)	Length (mm)
YM 1031-10	10 Fr	3	400
YM 1032-12	12 Fr	5	400
YM 1032-14	14 Fr	5	400
YM 1032-14s		5	280
YM 1032-16	16 Fr	10	400
YM 1032-16s		10	280
YM 1032-18	18 Fr	10	400
YM 1032-18s		10	280
YM 1032-20	20 Fr	10	400
YM 1032-20s		10	280
YM 1032-22	22 Fr	10	400
YM 1032-24	24 Fr	10	400



Endotherapeutics

Suite 1, Level 1, 51 Rawson Street
Epping NSW 2121, Australia

T: 1300 00 ENDO (3636) F: + 61 2 9869 7723 E: info@endotherapeutics.com.au
www.endotherapeutics.com.au